

How to Select Your Heating Tape



Use the following guidelines to determine which heating tape will fit your application.

Guidelines To Select Heating Tape

Moisture or Chemical Resistant

If moisture or chemicals will be present then Silicone Rubber extruded heating tape or moisture resistan insulation is recommended to heat your application.

Temperature Requirements

- Silicone Rubber extruded heating tapes are rated up to a maximum exposure temperature 232°C (450°F).
- Heavy and Standard Insulated heating tapes series are rated up to a maximum exposure temperature of either 482°C (900°F) or 760°C (1400°F).
- Heating cords are rated up to a maximum exposure temperature of either 482°C (900°F) or 760°C (1400°F)



Heavy insulated (BIH and BWH)



Standard insulated (B00 and BW0)



Heating Cords (HTC and HWC)

- Length and Width of Heating Tape
 The minimum length of heating tape needed for an application is equal to the length of the object. This will allow the tape to lay flat across the surface for the length of the system.
- For spiral wrapping, you will need to calculate the circumference of the object (Circumference=diameter x $\Pi(3.14)$) then multiply by the length of your object.

NOTE: The closer the pitch of the heating tape the longer the tape you will need to cover the entire surface of the object when spiral wrapping the heating tape. (minimum pitch distance is 1/8")

- For objects 1" diameter or less it is recommended to use 1/2" width heating tape to ensure direct contact.
- For larger diameters a wider tape may be used to maximize coverage.

Temperature Controlling Device
A temperature controlling device is required for all heating tapes. Some models (HSTAT, BSTS, and BSAT) come equipped with a built-in temperature controllèr. BriskHeat offers a wide range of temerature control options to meet your exact requirements.




Insulation

Insulation can be used to meet temperature requirements. Adding insulation to cover your heating tape significantly reduces the wattage needed to achieve the required change of temperature.

NOTE: It is recommended to use at least 1/4" of insulation for every 50°C change in temperature. This will lower the amount of wattage needed to meet your temperature requirements. The maximum amount of insultion to be used is 2-3 inches.





Example Application: 12' gas line, 2"ø, 2 hour heat-up with a Δt 56°C (100°F)

Heating tapes come in a wide range of sizes and properties. Use the table below to determine which heating tape is correct for your application

	Silicone Rubber						Cloth Insulated						Cord	
Product Series	HSTAT	BSTS	BSAT	BSO	BSO-G	CTL	BIH	BIH-G	BWHD	BWH	B00	BW0	HTC	HWC
Maximum Exposure Temperature	(1)	(2)	232°C	232°C	232°C	232°C	482°C	250°C	760°C	760°C	482°C	760°C	482°C	760°C
Moisture and Chemical Resistant	•	•	•	•	•	•								
Suitable for Electrical Conductive Surfaces	•	•	•	•	•	•	•	•	•	•			•	•
Integral tie Downs							•	•	•	•	•	•	•	•
Grounded		•			•			•						
Built-in-Control	•	•	•											

- (1) Adjustable Thermostat: 10°C 218°C (50°F 425°F)
- (2) Fixed Thermostat: 16°C (60°F)

Need help with your application? Contact BriskHeat and we will do a heat loss calculation for you.