
Cole-Parmer®

FSB-200 Series

Fluidized Sandbath

Instruction Manual
7002757-CPB Version 12.1

Cole-Parmer®
essentials

Contents

| | |
|---|----|
| Important supplementary safety information | 3 |
| General Description | 4 |
| Application | 4 |
| Before Use..... | 5 |
| Safety and installation..... | 6 |
| Sécurité et Consignes d'Installation | 8 |
| Sicherheits – Und Installationsinformationen | 10 |
| Información de Seguridad e Instalación..... | 12 |
| Informazioni sulla Sicurezza e l'Installazione..... | 14 |
| Warning..... | 16 |
| Principle of operation | 17 |
| Description of components and front panel controls..... | 18 |
| Installation | 19 |
| Operation | 20 |
| Start-Up..... | 20 |
| 2408 Controller Overview | 21 |
| Calibration of thermal sensors..... | 22 |
| Dead Bed Calibration | 23 |
| Automatic Air System..... | 23 |
| Shut Down..... | 23 |
| TechneWorks PC Software..... | 24 |
| Operator Maintenance | 27 |
| Emptying the dust collection jar | 27 |
| Cleaning the exhaust filter | 27 |
| Periodic Maintenance | 27 |
| To drain the alumina fluidizing medium | 28 |
| Fault finding | 29 |
| Service and repair | 30 |
| Technical Specification..... | 31 |
| Replacement parts | 31 |
| Accessories..... | 32 |
| Certificate of conformity..... | 34 |

IMPORTANT SUPPLEMENTARY SAFETY INFORMATION

Introduction

Cole-Parmer fluidized baths are safe and effective equipment when installed and operated correctly in accordance with the user manual. However, if used incorrectly they can pose a safety risk. Antylia Scientific have designed all models of fluidized baths to protect operators from hazards, but users should pay attention to the following points.

Caution

1. Please read the Instruction Manual before installation and use.
2. Some Cole-Parmer fluidized baths can heat up to 700°C. High temperatures are dangerous and can cause serious burns to operators and ignite combustible material.
3. Use care and wear protective gloves to protect hands and protective glasses to protect eyes.
4. Do not put hot objects on or near combustible objects.
5. Do not operate the unit close to inflammable liquids or gases.
6. Do not place any liquid directly in the unit.
7. Always ensure a suitable, adequate ventilation system is used when equipment is in use.
8. Always install fireproof metal ducting with sufficient airflow where applicable.

Maintenance

1. When performing maintenance, always disconnect from power supply and cool below 50°C.
2. Antylia Scientific recommend regular cleaning of fluidized baths. Externally, wipe with a damp soapy cloth. No abrasive cleaners should be used. Care should be taken to prevent any water entering the unit.
3. Regular internal and external inspection of extraction ducting is recommended to detect any damage and ensure the internals are clean. Any build-up of particles or debris discovered in the extraction ducting requires the ducting to be cleaned or replaced.
4. In fluidized baths used for polymer burn-off, please regularly inspect fluidizing medium, remove any foreign debris and replace with clean fluidizing medium as required.
5. Never top-up a hot fluidized bath with cold fluidizing medium. Always cool below 50°C first.

Please note:

1. Please ensure an adequate risk assessment is performed before use of a fluidized bath.
2. Please ensure the appropriate temperature is used for the application, always stay safely below the combustion temperature of any material or sample in a fluidized bath.
3. Fluidizing airflow must be switched on before heating a fluidized bath and left operational until the bath cools to below 50°C unless performing dead-bed calibration function.
4. Do not overfill fluidizing media. When fluidized, the aluminium oxide level should be approximately 150mm (6") from the top surface of the bath at your maximum operating temperature.
5. In fluidized baths used for polymer burn-off, always remove excess polymer from sample.
6. In applications where materials being treated produce acidic vapours during thermal decomposition, it is recommended a fume scrubber is utilised to ensure fume emission from the plant conforms to local regulations.
7. If you have any questions, please contact cptechsupport@antylia.com.

GENERAL DESCRIPTION

The Cole-Parmer FSB-200-P-AC employs the principle of fluidization of a mass of finely divided particles (alumina) to provide a safe, essentially isothermal environment with a high rate of heat transfer. The FSB-200-P-AC is specially designed to allow the operating temperature of the fluidized bath to be adjusted from a remote source while the bath is unsupervised. An automatic fluidizing air control system is fitted which adjusts the air flow rate accordingly to suit the set temperature of the bath. This automatic air feature is controlled by a control system which switches five solenoid valves. These valves are opened and closed in various combinations, providing eighteen air flow rates corresponding to eighteen different temperatures throughout the operating range.

To achieve a state of fluidization, a gas stream is forced to pass vertically through the bed of fine powder; chromatographic alumina (a specially refined grade of aluminium oxide) is the fluidizing medium used in the FSB-200-P-AC. At and above a certain critical value of flow rate, the particles become separated from each other and the whole mass behaves like a liquid. In this state, the powder will flow and move as though it were a liquid; low density objects will float on it, while more dense ones will sink; more importantly though, the individual alumina particles will circulate within the bath giving an essentially uniform temperature distribution throughout the working volume of the bath and a high degree of thermal coupling between the heat source and load.

In addition, the mass of fluidized powder changes its basic characteristics very little over a wide temperature range, has no freezing point and no practical boiling point and therefore by suitable choice of fluidizing media, the principle of fluidization can be used to temperatures over 1700°C and below -120°C.

The most commonly used fluidizing gas is compressed air. It is important that it should be clean, dry, free from oil and at a constant pressure. Antylia Scientific can supply a filter/regulator assembly (part code F5915) for compressed air which does not comply with the above requirement. Antylia Scientific are alternatively able to offer a suitable free-standing air compressor complete with the necessary filtration system (part code F120D) for installations where a compressed air supply is not available. Any other inert gas, such as nitrogen, could be used for special applications provided appropriate precautions are taken. The FSB-200-P-AC specification applies only when compressed air is used.

APPLICATION

The characteristics of the FSB-200-P-AC in terms of thermal uniformity and heat transfer make it a useful tool in the calibration and testing of a very wide range of temperature sensors; however, the unit is applicable to many other processes, such as thermal testing of sensitive components e.g. semiconductor devices, wire products, delicate transducers - and can be used as a constant temperature environment for chemical reactions.

Because the fluidized bed is composed of a fine dry powder, it does not display the surface tension effects of liquid baths and does not wet any objects immersed in it. The basic electrical insulating property of the alumina used in the FSB-200-P-AC is not affected by fluidization, making it possible to conduct electrical measurements on immersed objects such as assembled printed circuit boards.

Note: Avoid siting the FSB-200-P-AC in a laboratory environment which contains instruments that are sensitive to dust. Although the fluidized bath has its own dust extraction system, a small amount of fine dust may still be emitted during operation.

BEFORE USE

Thank you for purchasing this Cole-Parmer product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use. If there is any doubt relating to the proper use of this equipment, the staff at Anylia Scientific Ltd. Or your supplier will be happy to assist you.

UNPACKING

When unpacking the unit, check that the following have been removed from the packing:

| Part code | Item description |
|------------------|---|
| 12184-94 | FSB-200-P-AC Fluidized bath |
| 6007764 | Probe plate |
| 6007763 | Probe plate holder |
| 6003948 | Drain extension tube |
| 6007557/1 | Spare dust collection jar Instruction manual Eurotherm controller instructions Warranty card |

Please note that the aluminium oxide fluidizing medium must be ordered separately:

| Part Code | Item description |
|------------------|---------------------------------------|
| White/ALO | Alumina - white aluminium oxide, 25kg |

The FSB-200-P-AC requires 16kg of aluminium oxide

The user is advised to keep the original packaging in case the instrument ever needs to be returned for service or repair. Anylia Scientific Ltd. accepts no responsibility for damage incurred unless the unit is correctly packed and transported in its original packaging.

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by the manufacturer, the protection provided may be impaired.

SAFETY AND INSTALLATION

Please read all the information in this Instruction Manual before using the unit.

WARNING

HIGH TEMPERATURES ARE DANGEROUS: they can cause serious burns to operators and ignite combustible material.

Antylia Scientific have taken great care in the design of these units to protect operators from hazards, but operators should pay attention to the following points:

- USE CARE AND WEAR PROTECTIVE GLOVES TO PROTECT HANDS AND SAFETY GLASSES TO PROTECT EYES
- DO NOT use combustible substances near hot objects.
- DO NOT operate the unit in the vicinity of inflammable liquids or gases.
- DO NOT place any liquid directly into the equipment.
- At all times USE COMMON SENSE.

OPERATOR SAFETY

All operators of Cole-Parmer equipment must have available the relevant literature needed to ensure their safety. It is important that only suitably trained personnel operate this equipment in accordance with the instructions contained in this manual and with general safety standards and procedures. If the equipment is used in a manner not specified by Antylia Scientific, the protection provided by the equipment to the operator may be impaired.

All Cole-Parmer units have been designed to conform to international safety requirements and are fitted with an over-temperature cut-out. On some models, the cut-out is adjustable and should be set to suit the application. On all other models, the cut-out is pre-set to protect the unit.

If a safety problem should be encountered, switch off at the mains power supply.

INSTALLATION

Before connecting the instrument to the mains electricity supply, check the voltage against the rating plate (located on the back of the unit). **Please note that the unit must be earthed to ensure proper electrical safety.**

Connections 220V-240V

| | |
|---------|--------------|
| Live | Brown |
| Neutral | Blue |
| Earth | Green/yellow |

Note: The FSB-200-P-AC is classified as "Permanently Connected Equipment" and should be connected to the electricity supply by a qualified electrician.

A suitable supply for the FSB-200-P-AC is rated at 3kW, 220-240V, 50/60Hz~ single phase.

The equipment is supplied with 2m of flexible, triple core, circular cable to the following specification: 2.5mm², to BS 6500 or equivalent and <HAR> or BASEC approved. Connection to the mains electrical supply should be via a double pole 30mA Residual Current Breaker with Overcurrent protection (RCBO) isolating circuit breaker switch with a continuous current carrying capacity of 16A at 250V and overcurrent of 16A.

Do not switch on until the unit is fully installed.

ENVIRONMENTAL CONDITIONS

This equipment is designed to operate under the following conditions:

- Indoor use
- Use in a well-ventilated area
- Ambient temperature range +5°C to +40°C
- Altitude to 2000 m (6500 ft.)
- Relative humidity 80%
- Mains supply fluctuations not exceeding 10%
- Overvoltage category II IEC60364-4-443
- Pollution degree 2 IEC664

GUARANTEE

Notwithstanding the description and specification(s) of the units contained in the operator's manual, Antylia Scientific hereby reserves the right to make such changes as it sees fit to the units or to any component of the units.

This manual has been prepared solely for the convenience of Antylia Scientific customers and nothing in this instruction book shall be taken as a warranty, condition or representation concerning the description, merchantability, fitness for purpose or otherwise of the units or components.

OPERATOR MAINTENANCE

NOTE: THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL.

REMOVING THE SIDE, FRONT OR REAR PANELS EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your unit which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed and remember to return the unit in its original packing. Antylia Scientific accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

1. Cleaning

Before cleaning your unit ALWAYS disconnect it from the power supply and allow it to cool below 50°C.

Your unit can be cleaned by wiping with a damp soapy cloth. Care should be exercised to prevent water from running inside the unit. Do not use abrasive cleaners.

2. Over-temperature cut-out

In the event of no heater power, check the mains plug and lead. Repeated operation of the cut-out indicates a serious fault: you may need to return the unit to your supplier for repair.

3. Fuses

Your unit is protected by one or two fuses. These should only be changed by suitably qualified personnel. If the fuses blow persistently, a serious fault is indicated and you may need to return the unit to your supplier for repair.

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WARNING

Poor fluidization causes hot spots, heater failure and damage to other parts. Follow this Instruction Manual carefully. For correct fluidization, pay attention to the following:

INSTALLATION

- Ensure the bath is level and air supply is adequate.

OPERATION

- Adjust air valve for even fluidization.
- Do not insert objects larger than recommended.
- Ensure objects do not lie in contact with container wall or porous plate. Sample buckets and probe holders are available.

MAINTENANCE

- Regularly inspect and maintain the air filter to eliminate oil vapour in the air supply.

ALUMINA

- Before handling the alumina, consult the MSDS for safe handling. Contact cptechsupport@antylia.com if you require a copy.
- The principle risk is related to the concentration of dust in the air acting as a nuisance dust.
- Alumina is non-hazardous but dust may irritate the eyes and the respiratory system and cause skin irritation.
- Provide appropriate exhaust ventilation at places where dust is formed.
- For storage: protect from moisture. Avoid dust formation. Keep container tightly closed. Store in a well-ventilated place. Guard against dust accumulation of this material.
- Wear correct PPE whilst handling the Alumina material. Wear safety glasses with side shields, appropriate chemical resistant gloves and suitable protective clothing.
- Should the fluidized bath be stored for a long period under damp or humid conditions, moisture may be absorbed by the alumina which is hydroscopic. To avoid violent fluidization which occurs when damp alumina is heated above 100°C, operate the bath for a period of approximately 8 hours at 90°C prior to operation at elevated temperatures.
- NEVER ADD COLD OR DAMP ALUMINA TO A HOT BATH AS THIS WILL ALSO CAUSE VIOLENT FLUIDIZATION WHICH CAN BE DANGEROUS. Allow the bath to cool then add the fresh alumina. If this fresh alumina is a large portion of the charge, then dry the whole charge as above.

FUME EXTRACTION

- When used for processing items which may emit toxic or inflammable fumes, it is essential that an adequate fume extraction system be installed. The extraction system must be correctly sized to ensure that any toxic fumes are removed from the working environment. To eliminate the risk of spontaneous ignition, the concentration of inflammable fumes above the bath and within the exhaust duct work must be kept below the lower explosive limit.

PRINCIPLE OF OPERATION


Figure 1: Schematic of the FSB-200-P-AC

Control for temperature and a flow meter are mounted on a recessed panel on the front of the unit together with on-off switches for the power supply and the dust extraction system. The inner container is well insulated and the outer case is vented so that it remains safe and cool to the touch even when the bath is operating at its maximum temperature. The inner container is filled from the top with alumina. When fluidized, the alumina is heated by the four immersion heaters close to the inner container outer wall with the control thermocouple giving excellent temperature control.

The heater elements are protected from excessive surface temperatures by a pressure switch operated by the fluidized air and the power supply is removed if the air supply to the unit is disrupted. The unit also has an over-temperature cut-out fitted to protect the heater element in the event of an over temperature fault (limit set to 730°C).

Air supplied to the unit passes through two filters to remove oil and moisture from the system. A pressure regulator is also fitted to set the correct pressure for the cyclone extraction system and for optimum fluidization.

An automatic fluidizing air control system is fitted which adjusts the air flow rate accordingly to suit the set temperature of the bath. The fluidized air passes through a flow meter then to a plenum chamber fitted under the fluidized bath which evenly distributes the fluidized air evenly across the bath inner chamber.

Dust extraction is carried out by ambient air being drawn down (venturi effect) past the probe plate and through a peripheral slot around the top edge of the inner container. The air then passes through a cyclone to separate and remove any alumina particles which are collected in a glass collection jar and then passes through the exhaust filter. Entrained alumina removed by the cyclone and collected in the collection jar, which is accessible from the front of the unit, can be re-used and simply emptied back into the bath.

This model is specially designed to allow the operating temperature of the fluidized bath to be adjusted from a remote source while the bath is unsupervised. Basic control programmes are set by the operator allowing for the control of set temperature, incremental temperature steps, dwell times and control of dead bed state.

Where an ultra-stable temperature condition is required a “dead bed” state can be programmed into the control system. During this “dead bed” condition the air and electrical supply to the fluidized bed are switched off. For a period of up to 8 minutes, the fluidized bed becomes an isothermal mass without heat input and very low heat loss. Under these conditions the stability at the centre of the aluminium oxide is $\pm 0.01^\circ\text{C}$ over the range of the unit.

DESCRIPTION OF COMPONENTS AND FRONT PANEL CONTROLS


Figure 2: FSB-200-P-AC components and controls

| Item | Description | Function |
|------|---------------------------------------|--|
| 1 | Fluidized bath chamber | Heated chamber containing the fluidized alumina. |
| 2 | Exhaust filter unit | Collects remaining alumina from air passing through the cyclone. |
| 3 | Dust collection jar (6007557/1) | Collects entrained alumina removed by the cyclone. |
| 4 | Operating instructions | Quick operating guide. |
| 5 | Flow meter | A float indicates the fluidization air pressure set by the auto air control system. |
| 6 | Eurotherm 2408 temperature controller | When the mains switch is set to the ON position the indicator panel will display the measured and set temperature in digital form. The control parameters have been factory pre-set. |
| 7 | Air extraction switch | Turns the dust extraction system on and off. The extraction system should be used whenever the bath is running. |
| 8 | Air extraction indicator | Lights green when the dust extraction system is working. |
| 9 | Heater indicator | Lights amber when the heaters are on. Heaters are interlocked so that they will only turn on when fluidizing air is also on and exceeds a factory pre-set value. |
| 10 | Main power On/Off switch | Controls the mains input to the unit. |
| 11 | Fuse (6106042) | The unit is fitted with a 16A fast blow fuse |
| 12 | Spirit level | Used to ensure the unit is level. |
| 13 | Adjustable feet (x 4) | Adjustable to ensure the unit is level. |
| 14 | Compressed air inlet | The air supply should be connected via suitable flexible hose. |

INSTALLATION

Avoid siting the FSB-200-P-AC in a laboratory environment which contains instruments that are sensitive to dust. Although the fluidized bath has its own dust extraction system, a small amount of fine dust may still be emitted during operation.

Numbers in the text refer to Figure 2.

POWER

The FSB-200-P-AC requires a power supply of 230V, 50/60Hz at 3kW; please see section ELECTRICAL INSTALLATION.

AIR

The FSB-200-P-AC is supplied with integral filters and regulators suitable for a normal clean dry air line supply between 60 and 125 psi (414 to 862 kPa). If the air supply is excessively dirty Antylia Scientific can supply an additional external filter/regulator assembly which, if used, must be adjusted to its maximum setting of 60 psi (414 kPa). The air supply should be connected via suitable flexible hose to the compressed air inlet (14) at the rear of the left-hand end of the bath.

MECHANICAL

The unit should be mounted on a firm level surface. Normal operation only requires access to the front of the unit but periodic maintenance and servicing requires access to the left-hand side also; this should be borne in mind when positioning the unit initially. A spirit level (12) is built into the frame at the base of the front panel; this should be accurately levelled by adjusting the four corner feet (13). After final adjustment, all four feet should be firmly grounded and the unit should be level. Failure to do this can impair the operation of the unit. The surface should also be heat resistant.

ALUMINA

Turn on the dust extraction system (7) during filling to minimise the amount of dust given off. Fill the bath chamber (1) with sufficient of the alumina so that it is approximately 200mm (8") below the top cover of the unit before it is fluidized. Note that the alumina does exhibit expansion with fluidization and heat so you will need to verify the level at your maximum operating temperature. When fluidized, the alumina level should be approximately 150mm (6") from the top surface of the bath at your maximum operating temperature. As shown in Figure 3, the alumina is at the correct level when you just see the tops of the heaters.

Figure 3: Fluidized alumina at the correct operating level

When correctly fluidized, the surface of the alumina should be flowing with the motion like violently boiling cream.

Full fluidization is far from instantaneous and it takes several minutes to develop from a single vigorous ‘plume’ to a more distributed series of plumes and the level slowly rises even after the multiple plumes have formed.

OPERATION

The FSB-200-P-AC Fluidized Bath is capable of precise performance, but for reliable and consistent results it must be used by personnel conversant with its operation.

For temperature stability, the heat input to the bath must exactly equal heat losses; the heat input from the immersed heaters is varied by the controller which has a characteristic time constant for its response to a change in demand. The heat losses from the bath are due to:

1. heating fluidizing air
2. radiation from the top of the bath
3. conduction through the insulation
4. conduction along loads inserted into the bath
5. heating the mass of cold objects placed in the bath.

Changing any of these will affect the temperature of the bath, but it will recover under the controller's influence; this does mean, however, that after any parameter is changed or the load in the bath is changed, there is a delay before the bath temperature re-stabilises.

START-UP

1. Check that the unit has been set up as described earlier and that required services are available.
2. Check that alumina level is correct and adjust if necessary. Never add new, cold alumina (which may contain moisture) to a hot bath. Should the fluidized bath be stored for long periods of time under damp or humid conditions, moisture may be absorbed by the alumina which is hydroscopic. To avoid violent fluidization which occurs when damp alumina is heated above 100°C, operate the bath for a period of approximately 8 hours at 90°C prior to operation at elevated temperatures.
3. Connect the bath to the electrical and compressed air supplies.
4. Turn on the mains power switch (10). The internal neon should illuminate. The bath will begin to fluidise. The temperature controller will be illuminated and will indicate the actual temperature of the bath.
5. Turn on the dust extraction system switch (7); the green air extraction indicator lamp should illuminate.
6. Set the required operating temperature on the controller (6) by depressing the respective UP/DOWN button until the desired value is obtained. The control parameters have been factory pre-set. For comprehensive detailed information about the 2408 controller, refer to the operator's manual.
7. As the bath heats up, check that the surface of the alumina is just at the top of the “hair pin” bend in the heater assembly. Also check that the alumina is correctly fluidized; the surface of the alumina should be flowing with the motion of violently boiling cream. Any fine alumina particles emitted from the surface of the bath should be entrained by the extraction system when the probe plate is positioned above the bath.

2408 CONTROLLER OVERVIEW

For comprehensive detailed information about the 2408 controller, refer to the operator's user manual. This section gives a brief overview of the front panel display and the key parameters used.


Figure 4: Model 2408 front panel layout

- OP1 (Output 1) - When lit, it indicates that the output installed in module position 1 is on. This is normally the heating output on a temperature controller.
- OP2 (Output 2) – Module not installed.
- SP2 (Setpoint2) - When lit, this indicates that set point 2, (or a set point 3-16) has been selected.
- REM (Remote Set point) - When lit, this indicates that a remote set point input has been selected. 'REM' will also flash when communications are active.
- Auto/Manual button - When pressed, this toggles between automatic and manual mode: If the controller is in automatic mode the AUTO light will be lit. If the controller is in manual mode, the MAN light will be lit. The Auto/Manual button can be disabled in configuration level.
- Run/Hold button - Press once to start a program (RUN light on.) Press again to hold a program (HOLD light on). Press again to cancel hold and continue running (HOLD light off and RUN light ON). Press and hold in for two seconds to reset a program (RUN and HOLD lights off). The RUN light will flash at the end of a program. The HOLD light will flash during holdback or when a PDS retransmission output is open circuit.
- Page button - Press to select a new list of parameters
- Scroll button - Press to select a new parameter in a list
- Down button - Press to decrease a value in the lower readout
- Up button - Press to increase a value in lower readout

Basic Operation

Switch on the power to the controller. It runs through a self-test sequence for about three seconds and then shows the measured temperature, or process value, in the upper readout and the target value, called the set point, in the lower readout. This is called the Home display.

You can adjust the set point by pressing the up/down buttons. Two seconds after releasing either button, the display blinks to show that the controller has accepted the new value.

OP1 will light whenever output 1 is ON. This is normally the heating output when used as a temperature controller.

Operating Modes

The controller has two basic modes of operation:

- Automatic mode in which the output is automatically adjusted to maintain the temperature or process value at the setpoint. This is the normal operation mode.
- Manual mode in which you can adjust the output independently of the setpoint.

You toggle between the modes by pressing the AUTO/MAN button.

Two other modes are also available:

- Remote Setpoint mode, in which the setpoint is generated from an external source. In this mode, the REM light will be on.
- Programmer mode, further information on which can be found in the 2408 user manual.

CALIBRATION OF THERMAL SENSORS

When using the FSB-200-P-AC to calibrate thermal sensors such as thermocouples and platinum resistance thermometers, it is important to ensure that the working volume of the bath is thermally stable. Use an independent temperature sensor, such as a reference thermocouple, to check the temperature uniformity and stability. See section TECHNICAL SPECIFICATION.

The temperature stability of the fluidized bath at any operating temperature may be improved by tuning the control parameters of the temperature controller to match the performance of the system under any set conditions. Adjustment details are contained within the manufacturers instruction manual for the Eurotherm 2408 temperature controller.

For precise calibration, the actual temperature of the fluidized bath should be measured with an independent calibrated reference sensor; the temperature indicated by the temperature controller should only be used as a guide and not as an indication of absolute temperature.

The number of sensors calibrated in the fluidized bath at any one time should be kept to a minimum so as not to disturb the flow pattern within the bath. Typically, up to 10 sensors of approximately Ø10mm may be calibrated without loss of performance.

The probe plate supplied with the FSB-200-P-AC may be machined to allow the sensors being calibrated access to the bath. This probe plate should always be positioned above the surface of the bath in the probe plate carrier to inhibit the loss of fine alumina particles from the surface of the bath. See Figure 5. Accessory probe holders with integral dust extraction systems are also available. Probe holders locate the sensors being calibrated at a known position within the bath and act as equalising blocks so that the basic isothermal properties of the fluidized bath are further enhanced. See section ACCESSORIES.


Figure 5: Probe plate with probe plate carrier

Antylia Scientific can also offer a work basket (F7759) which ensures that probes or other items being processed are held within the working volume of the bath.

DEAD BED CALIBRATION

Temperature stability during calibration may be further improved by using a dead bed calibration technique. This technique involves heating the bath to the required calibration point, then allowing the temperature of the bath to stabilise for a period of approximately 20 minutes, dependent upon the actual calibration temperature and the thermal mass of the sensors being calibrated.

The fluidized bed is collapsed by cutting off the fluidizing air to the FSB-200-P-AC so that the sensors being calibrated are surrounded by a solid bed of alumina which acts as a thermal insulator. After a period of approximately 2 minutes at set temperature, a stability of 0.01°C can be achieved for a further period of approximately 8 minutes, dependent upon heat loss along the stem of the sensors being calibrated.

At the end of the required period of dead bed it is necessary to reconnect the fluidizing air to the FSB-200-P-AC.

AUTOMATIC AIR SYSTEM

An automatic fluidizing air control system is fitted which adjusts the air flow rate accordingly to suit the set temperature of the bath. This allows the operating temperature of the FSB-200-P-AC to be adjusted from a remote source e.g. using the TechneWorks PC software while the bath is unsupervised. This automatic air feature is controlled by a control system which switches five solenoid valves. These valves are opened and closed in various combinations, providing eighteen air flow rates corresponding to eighteen different temperatures throughout the operating range.

| Temp (°C) | Temp (°F) | Flow (LPM) |
|-----------|-----------|------------|
| 25 | 77 | 23.0 |
| 50 | 122 | 22.0 |
| 75 | 167 | 21.5 |
| 100 | 212 | 20.5 |
| 125 | 257 | 20.0 |
| 150 | 302 | 18.5 |
| 175 | 347 | 18.0 |
| 200 | 392 | 17.0 |
| 225 | 437 | 16.5 |
| 250 | 482 | 15.5 |
| 280 | 536 | 15.0 |
| 310 | 590 | 14.0 |
| 340 | 644 | 13.5 |
| 370 | 698 | 11.5 |
| 430 | 806 | 11.0 |
| 490 | 914 | 10.0 |
| 550 | 1022 | 9.5 |
| 610 | 1130 | 8.5 |

Table 1: FSB-200-P-AC automatic air flow rate table (factory settings)

When correctly set, the fluidizing air flow rate will vary according to the above table. It should be noted that the air flow requirements are approximate and accuracy is not critical. The air flow rates achieved will, however, vary according to the air supply pressure and whether the extraction system is switched ON or OFF.

SHUT DOWN

1. Turn off the mains power switch (10).
2. Turn off the dust extraction system (7).
3. Disconnect the compressed air supply.
4. Disconnect the mains power.

TECHNEWORKS PC SOFTWARE


The FSB-200-P-AC is specially designed to allow the operating temperature of the fluidized bath to be adjusted from a remote source while the bath is unsupervised. Fitted to the rear of the fluidized bath is an RS232 port which connects directly to the Eurotherm 2408 Comms1 RS232 module. The TechneWorks software is available to download from the Techne website. <http://www.techne-calibration.com/adminimages/TechneWorksv2.zip>.

Software Features:


- Specify either °C or °F.
- Specify ramp rates.
- Log data from the instrument while connected to the computer and export the data to an Excel® spreadsheet.
- Open, save, view and print logged data.
- Run a program in real-time mode.
- Specifying the logging interval from every 5 seconds to 60 seconds.
- Program daily start and stop times to automatically turn unit off and on for weekly schedules.

To set up remote operation:


1. Install the TechneWorks Software. We highly recommend that TechneWorks be installed on a PC that is typically not restarted at night for updates or Antivirus scans, though the software does have an auto restart feature if the PC is rebooted.
2. Connect the FSB-200-P-AC to an available coms port on your PC noting the port number.
3. Launch Techneworks and then select “Fluidized bath” as shown below.

4. The home screen will be displayed; click on “start device”.

5. From the “Serial Port” window select the Com port you are connecting to and click OK.

- After you click OK the screen should display the set point temperature and other conditions being constantly updated.

| Controller - Generic PID | | | | | | |
|--------------------------|-------|-------|-------|-------|------|------|
| | SP | PV | Alarm | wSP | SPrr | rMPU |
| COM16 | 25.00 | 18.20 | 0.00 | 25.00 | 0.00 | 0.00 |
| Start Device | Set | C | C | C | C | C |
| Interface | | | | | | |

- Clicking on the "Set" button next to SP will allow you to remotely change the bath set temperature and ramp rates.

Change setpoint value

Selected setpoint
Techne FB08

25.00

Ramp Rate 0.00 Max

OK Cancel

- Clicking the "On" and "Off" buttons will turn on or shutdown the bath.
- Bath shut down will occur once the Cool to" temperature has been reached. This temperature can be set by dragging the slider bar.
- Clicking the "Timer" button will allow you to set up daily schedules.

Setting up a bath start/stop sequence is fairly self-explanatory on this screen. Once the “On” time has been reached the bath will ramp up to the value set in the “Temp” window. At the “Off” time the bath will start cooling down to the setting in “Cool to value first”. After the “Cool to value first” is reached the bath will completely shut down and enter a sleep mode. The controller will display an FSH alarm when this occurs which is normal.

Note that the “Temp” value window will turn dark green when the program is running. After setting up your daily schedules be sure to click “Save and Apply Settings” to activate your schedule.

Checking off “Auto-Start Techneworks” will restart the program if the PC is rebooted.

OPERATOR MAINTENANCE

Before any servicing is attempted the unit should be disconnected from the mains power supply and allowed to cool down. In all cases maintenance and repair work should be undertaken only by a skilled technician. Untrained personnel should not attempt to dismantle the instrument.

EMPTYING THE DUST COLLECTION JAR

Empty the dust collection jar as follows (this can be done whilst bath is hot):

1. Turn off the fluidization air.
2. Turn off the dust extraction system.
3. Unscrew the jar from the base of the cyclone unit and remove.
4. Screw on the empty spare jar.
5. Turn on the dust extraction system.
6. Carefully and slowly empty the full jar into the top of the bath. (Except for the first jar from a new charge which should be thrown away).

Note: During normal operation, the collection jar should only become partially full during an 8-hour period. An excessive amount of alumina in the collection jar suggests that the fluidized bath has been overfilled with alumina; check the level in the bath. Add enough from the jar to bring the bath to the correct level and keep the remainder as spare material.

7. Increase fluidization up to the correct working level.
8. Check the filter and clean or replace if necessary, when the extraction system is switched off. (Note the screw and the filter may be hot).

CLEANING THE EXHAUST FILTER

1. Turn off the extraction system.
2. Release the filter housing by removing the retaining screw and vacuum the filter chamber every two weeks with daily bath use.
3. Replace the filter every 6 months with daily bath use.

PERIODIC MAINTENANCE

The fluidized bath requires the following periodical inspection to be carried out to ensure trouble free operation:

1. The condition of the alumina should be checked to make sure that there are no lumps and, if practical, sieved occasionally. A fresh charge of alumina should be used if the medium is contaminated in any way.
2. Should the fluidized bath be stored for long periods of time under damp or humid conditions, moisture may be absorbed by the alumina which is hydroscopic. To avoid violent fluidization which occurs when damp alumina is heated above 100°C, operate the bath for a period of approximately 8 hours at 90°C prior to operation at elevated temperatures. The heaters should be checked for discolouration, oxidation or to determine if any build-up of alumina particles has taken place on the surface.
3. The condition of the porous plate should be checked to see whether there are any oil stains or discolouration: if so, and fluidization is impaired, replace the porous plate. Contact your local Cole-Parmer dealer for further information.

TO DRAIN THE ALUMINA FLUIDIZING MEDIUM

This is best done when the bath is cool as hot alumina can cause severe burns. The drain cap may also be hot.

1. Ensure the fluidized air is turned off.
2. Position a suitable container under the drain situated at the bottom rear of the bath (Figure 6a).
3. Remove the drain cap at the bottom rear of the unit and fit the drain extension tube (do not over-tighten) (Figures 6b and c).
4. Set the Eurotherm 2408 controller to ambient (20°C). The bath will start to drain at a steady flow rate of alumina.
5. Once the chamber is empty, remove the drain extension tube, clean the thread of the drain cap and re-fit.


Figure 6: Draining off the alumina.

FAULT FINDING

The following guide has been prepared should a fault develop with the fluidized bath. If the fault cannot be located by using this guide, contact your local Cole-Parmer dealer.

| Fault | Possible cause |
|--|---|
| The unit fails to operate | Check the fuse, electrical and air supplies |
| Poor fluidization | Unit not level Inadequate fluidizing air flow Inadequate feed air supply Air filters blocked Defective air flow meter Bath over-filled with alumina Fluidization restricted by large mass placed in bath Porous plate blocked or damaged |
| Bed will not fluidise | Inadequate air flow Inadequate feed air supply Check pipe work for kinks, leaks or poor connection Fluidization pressure switch inoperative. |
| Poor temperature control | Poor fluidization Eurotherm controller faulty Transformer faulty Thermocouple faulty |
| Power on but air controller does not illuminate. | Inadequate feed air supply Poor hose connection Bath alumina level too low Air flowmeter faulty Inadequate fluidizing air flow |
| The flowmeter float does not rise. | Check the air inlet supply is 60 psi Check pipe work for kinks, leaks or poor connection Air flowmeter faulty Airflow solenoid valves not working Flowmeter valve faulty |
| Air extraction not working | Check the air inlet supply is 60 psi Check the air extraction air pressure is 50 psi Check pipe work for kinks, leaks or poor connection Venturi requires setting up Exhaust filter clogged |
| Poor heating up of the unit | Heater element faulty Thermocouple faulty Front solid state relay faulty Rear solid state relay faulty Eurotherm controller faulty or parameters incorrect Over temperature circuit not set at correct over current value |
| Dust collection jar is filling too quickly | Bath over-filled with alumina Bath is over-fluidizing; reduce air flow |
| Dust leaking from collection jar | Collection jar damaged or chipped; replace jar Seals damaged or loose |
| Remote temperature control not working | Check the RS232 connection between the unit and PC comms connection Check the remote button is illuminated on the 2408 controller Check TechneWorks software installed correctly 2408 controller comms PCB faulty |

SERVICE AND REPAIR

NOTE THAT THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. **REMOVING THE CASE EXPOSES POTENTIALLY LETHAL MAINS VOLTAGE.** THERE ARE NO OPERATOR MAINTAINABLE PARTS WITHIN THE EQUIPMENT.

In the unlikely event that you experience any problems with your fluidized bath which cannot easily be remedied, you should contact your supplier and return the unit if necessary. Please include any details of the fault observed, quoting the serial number and remember to return the unit in its original packing. Antylia Scientific accept no responsibility for damage to units which are not properly packed for shipping: if in doubt, contact your supplier.

Only spare parts supplied by the manufacturer or its agent should be used. Fitting of non-approved parts may affect the performance of the safety features of the equipment. If in doubt, please contact Antylia Scientific.

If you require further technical or application assistance, please contact Antylia Scientific at:

E-mail: cptechsupport@antylia.com

Phone: +44 (0)1785 810433

Fax: +44 (0)1785 810405

TECHNICAL SPECIFICATION

General characteristics

| | | |
|-----------------------|------------------------------------|-------------------------|
| Temperature range | 50°C to 700°C | |
| Temperature stability | $\pm 0.01^\circ\text{C}$ | |
| | Dead bed stability (for 8 minutes) | |
| | At 50°C: | Short term |
| | | $\pm 0.2^\circ\text{C}$ |
| | | Long term |
| | At 600°C: | Short term |
| | | $\pm 0.5^\circ\text{C}$ |
| | | Long term |
| Nominal heater power | | $\pm 0.3^\circ\text{C}$ |
| Heat up time | From 20°C to maximum | 4 x 750W |
| Cool down time | From maximum to 200°C | 1hr 45min |
| Electrical supply | | 2hrs 45min |
| | | 230V, 50/60Hz |

Air supply

| | |
|--------------|------------------|
| Pressure | 60 psi (414 kPa) |
| Maximum flow | 127 l/min |

Nominal dimensions and weights

| | | |
|---------------------|--------------------------|-------|
| Overall | Width | 770mm |
| | Depth | 515mm |
| | Height | 600mm |
| Internal dimensions | Diameter | 165mm |
| | Depth | 493mm |
| | Working depth | 385mm |
| Weight | Unit (including alumina) | 75kg |
| | Chromatographic alumina | 16kg |
| | Basket | 2kg |

REPLACEMENT PARTS

The following parts may be purchased if replacements are required:

| Part number | Description |
|--------------------|-------------------------------------|
| 6007558 | Exhaust filter |
| 6007557/1 | Dust collection jar |
| 6003948 | Drain extension tube |
| 6007764 | Probe plate |
| 6007763 | Probe plate holder |
| 6106042 | 16A fast blow fuse |
| White/ALO | Alumina, white aluminium oxide 25kg |

ACCESSORIES

The following accessories are available for use with the FSB-200-P-AC.

Sample basket

The steel mesh sample basket safely holds samples away from the heating elements of the fluidized bath and assists retrieval of items from the bath. The basket contains a collar to reduce alumina spilling from the chamber.

| Product Code | Description |
|---------------------|--------------------|
|---------------------|--------------------|

| | |
|-------|---|
| F7759 | Steel mesh sample basket with collar and sensor plate |
|-------|---|

Dust suppression system and probe holder

The Dust Suppression System chimney (FFB08DS1) is combined with a Probe Holder (FFB08PR1 or FFB08PR2) to enable the FSB-200-P-AC user to seal the fluidized bath, making it suitable for calibration labs where airborne dust must be avoided. In addition, Probe Holders allow free flow of alumina around the exterior of the holder preventing dead areas. At the base of each Probe Holder is a large steel equalization block surrounding the sensors, assuring constant uniformity, reduced temperature fluctuation and improved calibration accuracy. Stability and uniformity are >0.010°C in dead-bed mode with a Dust Suppression System.

Note: Dust Suppression System and Probe Holder must be purchased together for correct fitting.

| Product Code | Description |
|---------------------|--------------------|
|---------------------|--------------------|

| | |
|----------|---|
| FFB08DS1 | Dust Suppression System with chimney |
| FFB08PR1 | Probe Holder for 8 probes: 8mm x 281mm deep |
| FFB08PR2 | Probe Holder for 8 probes: 3mm, 4mm(x3), 5mm, 6mm(x2), 8mm x 300mm deep |
| FFB08FL1 | Spare filter for Dust Suppression System |

Air Compressor

Antylia Scientific can supply a free-standing moisture and oil free air compressor complete with the necessary filtration system for installations where a compressed air supply is not available.

| Product Code | Description |
|---------------------|--------------------|
|---------------------|--------------------|

| | |
|-------|----------------------------------|
| F120D | Air Compressor 230V, 50Hz, 1500W |
|-------|----------------------------------|

Air pressure regulator and filter

Antylia Scientific can supply a filter/regulator assembly for compressed air which is not moisture or oil-free.

| Product Code | Description |
|---------------------|--------------------|
|---------------------|--------------------|

| | |
|-------|-----------------------------------|
| F5915 | Air Pressure Regulator and Filter |
|-------|-----------------------------------|

NOTES

This product meets the applicable CE Directives and UKCA Legislation for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards and so we cannot guarantee

that interference will not occur in practise. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advise before use, contact the manufacturer.

Declaration of Conformity is available to view online at www.coleparmer.com

EU Representative address
Antylia Scientific GmbH
Futtererstraße 16
97877 Wertheim
Deutschland Tel: +49 9377 9203-0 Email:
sales@coleparmer.de

UK Representative address Antylia Scientific
9 Orion Court
Ambuscade Road
Colmworth Business Park St.
Neots PE19 8YX
United Kingdom
Tel: +44 (0) 1480 277339
Email: enquiries@antylia.com

Ordering Information

| Order No. | Series | Model | Legacy SKU |
|-----------|---------|--------------|------------|
| 12184-94 | FSB-200 | FSB-200-P-AC | F949J |

Warranty Registration

Cole-Parmer®
essentials

Antylia Scientific Ltd.
Beacon Road,
Stone,
Staffordshire,
ST15 0SA,
United Kingdom

UK

T: +44 (0) 1480 272279
E: uk.sales@antylia.com
W: coleparmer.co.uk

India

T: +9122 61394444
E: info@coleparmer.in
W: coleparmer.in

Germany

T: +49 (0) 9377 92030
E:
de.sales@antylia.com
W: coleparmer.de

China

T: +1 847 549 7600
E: sales@antylia.com
W: coleparmer.com

France

T: +33 (0) 1486 37800
E: fr.sales@antylia.com
W: coleparmer.fr

USA

T: +1 847 549 7600
E: sales@antylia.com
W: coleparmer.com

Italy

T: +39 (0) 284349215
E:
it.sales@antylia.com
W: coleparmer.it

Canada

T: +514 355 6100
E: info@antylia.ca
W: coleparmer.ca

Other

T: +1 847 549 7600