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High precision  
Tubing pumps  
with planetary drive

IP		
4	Channels	ISM940
8	Channels	ISM941
12	Channels	ISM942
16	Channels	ISM943
24	Channels	ISM944

IP-N		
4	Channels	ISM945
8	Channels	ISM946
12	Channels	ISM947
16	Channels	ISM948
24	Channels	ISM949

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## Please note

We recommend you read this operating manual carefully. When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump must not be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes

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## Caution

Rotating pumphead creates a pinch and crush hazard.

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## Safety precautions

- ▶ The pump must be switched OFF when cassettes or tubing are inserted or changed.
- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Pumpheads consist of rotating parts. Therefore, the pump must not be operated before all cassettes are inserted and fully snapped-in.
- ▶ Do not manipulate the pumphead or the cassettes before the pump is switched off and disconnected from the mains supply.

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## Safety precautions

- ▶ Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the rotating roller head.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media it is recommended that the pump be placed in a collecting basin.
- ▶ Only new fuses, according to the specifications stated on Page 13 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.

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## Safety precautions

- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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## Warranty

**2 years from date of delivery.**

**Tubes are regarded as consumables and as such they are not covered by the warranty.**

In case of any queries, please contact your ISMATEC® representative.

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## Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

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## Warranty terms

### Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third-party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications

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## Product

### Contents of the package

- ▶ IP or IP-N pump (type as ordered) incl. Click'n'go cassettes (POM-C)
- ▶ 1 power cord with an IEC 320 connector (female plug) and a country specific mains plug
- ▶ Operating Manual

### Please check the package

and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

- ➔ Complaints can only be accepted within 8 days from receipt of the goods.

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
## Rear panel


1. Analog interface  
Input for:
  - speed control  
0–5 V or 0–10 V, and  
0–20 mA or 4–20mA
  - Run/Stop
  - rotation direction
2. Mains switch (ON/OFF)
3. Fuse-holder with voltage  
selector 115/230 V
4. Mains socket

**⚠** For further details on the analog interface please refer to Pages 26–29.

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#### Window for voltage setting

If necessary, the voltage setting must be changed.

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### ⚠ Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage.

**If necessary, the voltage setting must be changed.**

### ⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be connected to earth ground. (protective conductor contact)

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
## Changing the fuses

### Voltage setting 115V/ 230

➔ Switch the pump OFF, pull out the mains plug.

1. Pull out the fuse-holder by opening it with a small screw-driver (size 0).
2. Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
3. Insert 2 new fuses  
230 V<sub>AC</sub>: 2 x T500mA/250V  
115 V<sub>AC</sub>: 2 x T500mA/250V

⚠ Use always 2 slow-blow fuses complying with the local mains voltage.

4. Insert the fuse-holder. The voltage rating is visible in the window.

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
1. LED-Test »8.8.8.8.«
2. Version der Systemsoftware

## Starting the pump

- ▶ Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage (Page 12).
- ▶ Prepare the cassettes (Page 15)
- ▶ Insert the tubing (Page 16)
- ▶ Insert the cassettes (Page 16)
- ▶ Connect the pump tubing.
- ▶ Connect the pump to the mains and switch it ON with the power supply switch.

## Start-up information

After switching the power supply switch ON, the following values are displayed:

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

## Preparing the cassettes

1. ➔ Switch the pump OFF.
2. Remove the cassette by slightly pressing the release lever on both sides of the pump.
3. The cassettes have slits on both sides. Insert a cassette adaptor into the slit on each side of the cassette.

Press the adaptor firmly against the bottom of the slit.

➔ Complete insertion of the adaptor will allow easy loading of the tubing into the cassettes.

### ⚠ **Pump with planetary drive.**

The geared rollers rotate simultaneously with the shaft of the roller-head.

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
## Inserting the tubing

1. ➔ Switch the pump OFF.
2. Insert the 2-stop tubing with one stop into an adaptor. Stretch the tubing slightly so that it can be inserted with the second stop into the other adaptor.

➔ Ensure that the tube is not twisted.

3. Reinsert all cassettes into the pump. Ensure that they snap in on both sides.

➔ The direction of the arrow engraved on the cassette indicates the flow direction providing the least pulsation (see information on Page 30).

4. Switch the pump ON and start it.

➔ **When the pump is idle it is recommended that pressure be released from the tubing by disengaging all cassettes on one side.**

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
## Operating panel

1. Digital LED display
2. Increase speed
3. – Reduce speed  
– settings\*
4. – MIN key (for standby)  
– default\*
5. – MAX key (ideal for filling or  
emptying the tubing)  
– ok\*
6. – Rotation direction  
– escape\*
7. RUN/STOP
- 8. Operating mode**  
– **PUMP (local operation)**  
– REMOTE (control via analog  
interface (Pages 26–29))

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## Control keys

The control keys are not active under remote control

### a. RUN/STOP

- starts the pump at the preset speed (the rotation direction is briefly indicated with a running dash)
- stops the pump

### b. Rotation direction

- changes the rotation direction
- ➔ minus sign = counter-clockwise direction

### c. ▲▼ Keys

- ▲ = increase speed
- ▼ = reduce speed

### d. MAX

- max. rotation speed (keep key pressed while the pump is running)

### e. MIN

- The pump runs at 4% of the max. speed (e.g. for stand-by over night).

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
## Control keys for the basic settings

- a. **settings**  
Allows access to the menu settings
  - b. **ok**  
Confirms a data entry
  - c. **escape**
    - Cancels an entry
    - Quits the menu settings
  - d. **default**  
Sets a value to its default value
  - e. **▲▼ Keys**
    - Select the basic settings
    - Select value
- ➔ **Basic settings see Pages 19–21.**

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
## Selecting the basic settings

1. Push the settings key while switching the pump ON
2. Select adjustable operating settings with the ▲▼ keys
  - FS (footswitch)
  - Speed (rotation speed)
  - Init (default values)
3. ok  
Confirm your selection
4. Select the value by using the ▲▼ key, or
  - press escape for cancelling wrong input
  - certain values can be reset to the default value by using the default key
5. ok  
Confirm entered value
6. escape  
Quit menu basic settings

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## Basic settings

### Footswitch or hand dispenser

- ▶ Toggle (default value)
  - Step on footswitch = pump starts running
  - Press footswitch again = pump stops
- ▶ Direct: the pump runs as long as the footswitch is pressed

### Speed

- ▶ % (default value)
  - Rotation speed in percent
  - Setting: 1.0 – 100.0 %
  - Resolution: 0.1%
- ▶ RPM
  - Rotation speed in rpm
  - Setting: 0.4 – 45.0 rpm (IP)
  - 0.11 – 11.25 rpm (IP-N)
  - Resolution: 0.1 rpm (IP) / 0.03 rpm (IP-N)

### Initializing (default value)

The ok key resets all basic settings to their default values.

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## Pumping

- 1 Set the required rotation speed (in rpm or %) with the ▲▼ ▶ keys.
- 2 Enter the rotation direction with the ◀▶ key.
- 3 Start the pump with the RUN/STOP key.

- ➔ The flow rate can also be adjusted while the pump is running.
- ➔ For fast filling or emptying the system, maintain pressure on the MAX key while the pump is running.
- ➔ For stand-by, e.g. during the night, push the MIN key. The pump runs at 4% of the maximum speed.

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## Pumping against pressure

The IP pump can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures).

In case of any doubts please contact your ISMATEC® agent.

## When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure. Releasing the cassette on the right side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

### Syphoning effect

When pressure is released from the tubing, fluid can flow back to the reservoir.

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## Overload protector

The IP pump features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL I = Current  
OL S = Speed  
OL t = Temperature

**Whenever this situation occurs, the pump must be switched off immediately.**

Before the pump is re-started, it is most important to check the reason for the overload (e.g. too high differential pressure or pressure lever tightened too much, etc.).

➔ Only when the operating conditions have been changed accordingly may the pump be started again.

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## Please Note

Please also refer to our detailed tubing documentation.

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## Running-in period for tubing

Every new tube requires a run-in period. If constant and repeatable flow rates are required, it is recommended that new tubing be run-in with water or the medium to be pumped for at least 1 to 3 hours before starting the application.

## Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog.

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
## Analog interface

### Pin 1, GND (ground)

Reference potential for all other inputs

### Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

### Pin 3, start

- In normal operation (Pin 2 open) Pin 3 operates the footswitch (see operating settings, Page 21).

- In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND)

### Pin 4, direction

In the open position the pump turns clockwise; when connected to Pin 1 (GND) it turns counter-clockwise.

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
Pin 2, remote / Pin 3, start  
 Pin 4, direction / Pin 13, speed intern

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## Analog interface

### Pin 5, speed IN

External speed control  
 (0–5V, 0–10V, 0–20mA, 4–20mA)  
 Input impedance and input range can  
 be selected via a dip-switch inside the  
 pump (see Page 28)

### Pin 7, +30 V<sub>DC</sub>

About +30 V<sub>DC</sub> are available  
 (max. current 0.2 A).

### Pin 9, speed OUT

0–10 V<sub>DC</sub>, proportionally to the speed  
 0–45.0 rpm (IP) or  
 0–11.0 rpm (IP-N).

### Pin 10, +5 V<sub>DC</sub>

About +5 V<sub>DC</sub> are available.  
 (max. current 0.1 A)

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
## Analog interface

### Pin 13, speed intern

Depending on how the pump is operated, Pin 13 has different functions:

Analog interface not activated (Normal operation, i.e. Pin 2 is open) Pin 13 serves as **auto-start function**. If Pin 13 is connected to Pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be ON).

Analog interface activated (Pin 2 on GND)

- Pin 13 open: The rotation speed is adjusted via Pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
## Settings of DIP-Switch

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected for some time.

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(order 2 adaptors per cassette)

<sup>1)</sup>The cassettes with pressure lever are available on request. This type of cassette provide better results at elevated differential pressure conditions (>1 bar).

Preparing these cassettes see Page 32.

<sup>1)</sup>Materials POM-C / PVDF see Page 31

### ⚠ Click'n'go-Cassettes:

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case, it is recommended that the tubing be primed and the pump be run with the tubing inserted for about 15 to 30 minutes.

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### Material POM-C

Polyoxymethylene-Copolymer

- good chemical resistance to many organic solvents and strong alkaline chemicals
- is affected by strong acids and oxidizing substances
- UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

### Material PVDF

Polyvinylidene fluoride

- very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
- not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
- stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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Please refer to our detailed tubing documentation or our website

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## Preparing the cassettes with pressure lever

- 1 ➔ Switch the pump OFF.
- 2 Push the pressure lever down (it must be relaxed).  
Insert the tubing with 2 colour-coded collars.
- 3 Reinsert the cassette(s) onto the rollers. Ensure that it snaps in on both sides.
- 4 ➔ Switch the pump ON and start it.
- 5 Lift the pressure lever step by step until the medium is drawn in and a continuous flow is achieved (no further step is necessary!)

### Please Note

- Under certain circumstances (e.g. new tubing, elevated differential pressure, etc.) we recommend the user to lift the lever by just one more step.
- Never over-tighten the pressure lever.
- The setting of the pressure lever determines the occlusion.
- The optimum occlusion depends on tubing material, inner diameter, wall thickness, age of tubing, and application conditions.

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## Preparing the cassettes with pressure lever

### ⚠ Please Note

- As a side-effect the occlusion setting also affects the flow rate. This, however, should not encourage the user to use the pressure lever for flow rate adjustments because long-term behaviour is not predictable.
- If, however, the occlusion setting is used to adjust the flow rate, the channel-to-channel flow rate conformity may lead to good results in the beginning of the operation, yet the individual flow rates may increasingly diverge at a later stage.
- For an optimum long-term flow rate stability and channel-to-channel conformity we suggest using the Click'n'go cassettes.

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## Accessories

### Footswitch IS 10039

This footswitch serves as a Start/Stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

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## Flow rates

In the tables listed subsequently the values indicated are only approximate. They were determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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## Please Note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®.

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## Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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## Maintenance

Provided the pump is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

## Repairs

For repairs please send the defective pump with detailed description of the failure to your ISMATEC® representative.

Please use the original ISMATEC® packaging or packaging of equal quality.

Please give information on defect, date of purchase, serial-no., and model.

## Spare parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® representative:

- ▶ replacement parts
- ▶ parts lists
- ▶ wiring diagrams

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Motor type DC motor

### Speed

0.4 – 45.0 rpm (IP)  
0.11– 11.25 rpm (IP-N)  
digitally adjustable in steps of 0.1% or  
0.1 rpm (IP)  
0.03 rpm (IP-N)

### Pump rollers

8 rollers, stainless steel  
actively driven

### Differential pressure

max. 1.0 bar (14.5 psi)  
(For short term use, higher pressures are  
possible provided that small tubing I.D.s  
and/or cassettes with pressure levers  
are used)

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## Please Note

We also recommend you observe our **Warranty Terms as well as our Terms and Conditions of Sale**.

In case of any queries, please contact your local ISMATEC® representative.

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max. 30 W

### Operating conditions

Temperature +41 to 104°F  
Rel. humidity max. 80%  
– not condensing, at normal laboratory conditions

### Dimensions/Weight

Channels	DxWxH (inches)	Weight
4	7.1x5.7x5.1	10.1 lb
8	8.7x5.7x5.1	11.2 lb
12	10.2x5.7x5.1	12.8 lb
16	11.8x5.7x5.1	14.3 lb
24	15.0x5.7x5.1	17.4 lb

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conformément à:  
EN 61326-1, EN 61010-1

# ***ISMATEC***

Learn more  
about the  
Ismatec® product  
line by visiting:  
[www.ismatec.com](http://www.ismatec.com)

**For ordering and technical support, please contact:**

**North America** [sales@ismatec.com](mailto:sales@ismatec.com) | 1-800-323-4340 | 1-847-549-7600

**Europe** [sales.europe@ismatec.com](mailto:sales.europe@ismatec.com) | +49 (0) 9377 9203-0