English

## SDC Benchtop Digital Temperature Controller

### Instruction Manual

You must read and understand this manual before installing, operating, or servicing this product. Failure to understand these instructions could result in an accident causing serious injury or death.

Keep these instructions for future reference.


Language	Page
English.....	1
Spanish (Español).....	9
French (Français).....	17
German (Deutsch).....	25
Italian (Italiano).....	33

### TABLE OF CONTENTS

Introduction.....	2
Safety Alert Symbol.....	3
Important Safety Instructions.....	4
Specifications .....	4
Parameters.....	5
Setting The Set Point .....	6
Parameter Descriptions.....	6-7
Indicators.....	7
Access To All Codes Protected Parameter .....	7
Warranty.....	8

### INTRODUCTION

BriskHeat®'s SDC Benchtop Digital Temperature Controller is designed for general-purpose use in indoor environments to control temperatures on applications requiring automatic control. For successful operation of this controller, read and understand these instructions prior to use.

#### SAVE THESE INSTRUCTIONS!

Additional copies of this manual are available upon request.

## SAFETY ALERT SYMBOL

The symbol above is used to call your attention to instructions concerning your personal safety. It points out important safety precautions. It means "ATTENTION! Become Alert! Your Personal Safety is involved!" Read the message that follows and be alert to the possibility of personal injury or death.

Immediate hazards which **WILL** result in severe personal injury or death.

A person who has not read and understood all operating instructions is NOT qualified to operate this product.

Hazards or unsafe practices that **COULD** result in severe personal injury or death.

Hazards or unsafe practices that **COULD** result in minor personal injury or property damage.

- Do not immerse or spray any component of the temperature controller in liquid.
- Keep volatile or combustible material away from controller when in use.
- Keep sharp metal objects away from controller.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment can be impaired.
- The unit must be connected to Protective Earth Grounding.

Failure to observe these warnings may result in electric shock, risk of fire, and personal injury or death.

- Inspect all components before use.
  - Do not use controller if any component is damaged.
  - Do not repair a damaged or faulty controller.
  - Do not crush or apply severe physical stress on any component of system, including cord assembly.
  - Unplug controller when not in use.
  - Only use power cords provided by BriskHeat®.
  - Unit must be mounted away from vibration, impacts, water and corrosive gases.
  - The box containing the controller should be stored indoors in a dry, protected area. Under no circumstances should condensation, rain, snow or water be allowed to come in contact with the box or the controller while in storage.
- Failure to observe these warnings may result in personal injury or damage to the temperature controller.


### End User Must Comply with the Following:

- Only qualified personnel are allowed to connect electrical wiring.
- All electrical wiring must follow local electrical codes and NEC Article 427.
- The person who performs the final installation / wiring must be qualified for this work.
- The end-user is responsible for providing a suitable disconnect device.
- The end-user is responsible for providing a suitable over-protection device. It is highly recommended that a ground-fault circuit breaker be used.

For best accuracy of temperature measurement, do not install probe cable near power cables.

## SPECIFICATIONS

Input Type	Temperature Range		Accuracy
	°C	°F	
J-type Thermocouple (P5)	-40°C to 370°C	-40°F to 698°F	±1% FS
K-type Thermocouple (P5)	-40°C to 482°C	-40°F to 899°F	±1% FS
<b>Conditions</b>			
Storage Temperature	(-20 to 80°C)	(-4 to 176°F)	
Temperature Limits: Ambient	0 to 80°C	32 to 158°F	
Output	16A SPDT relay @ 250VAC resistive		
Horsepower Rating	1 HP		
Control Type	ON/OFF		
Power Requirements	110VAC, 230VAC, 12VAC/VDC or 24VAC/DC		
Display	3-digit, red, 1/2" (12.7mm) digits, plus sign		
Resolution	1°		
Memory Back-up	Nonvolatile		

### NOTE:

The SDC comes with input power fuses of 100mA fast-acting 250VAC rating for the control circuit.

## PARAMETERS

	Description	Units	Range	Default Value
<b>SP</b>	Set Point	Degrees	r1 to r2	100
<b>r0</b>	Differential or Hysteresis	Degrees	1 to 99°	1
<b>r1</b>	Low Value Set Point	Degrees	0 to 999°	0
<b>r2</b>	Higher Value Set Point	Degrees	0 to 999°	600
<b>d0</b>	Heating or Cooling Control	Option	Ht/Co	Ht
<b>c0</b>	Min. Stop Time for Load	Minutes	0 to 59	1
<b>c2</b>	Load Status During Probe Error	0/1	oFF or on	oFF
<b>P1</b>	Probe Adjustment	Degrees	-30° to +30°	*Varies
<b>P5</b>	Probe Type	Option	J, K, S	**tch/tcJ
<b>P6</b>	Probe Response	Numeric	0 to 3	3
<b>H5</b>	Parameter Access Code	Numeric	0 to 255(factor set 0)	118
<b>A0</b>	Alarm 1 Hysteresis	Degrees	1 to 999°	1
<b>A1</b>	Alarm 1 Threshold	Degrees	0 to 999°	0
<b>A2</b>	Alarm 1 Exclusion Time	Seconds	0 to 999	0
<b>A3</b>	Alarm 1 Configuration	Option	oFF, Lo or Hi	oFF
<b>A4</b>	Alarm 2 Hysteresis	Degrees	1 to 999°	1
<b>A5</b>	Alarm 2 Threshold	Degrees	0 to 999°	0
<b>A6</b>	Alarm 2 Exclusion Time	Seconds	0 to 999	0
<b>A7</b>	Alarm 2 Configuration	Option	oFF, Lo or Hi	oFF

\* Value varies— Factory set

\*\* tch notation used for K

tcJ notation used for J

## SETTING THE SP

- Set point (SP) is the only parameter the user can access without code protection.
- Press SET and the SP text will appear on the display.
- Press SET again and the real value is shown on the display.
- The value can then be modified with the UP and DOWN arrows.
- Press SET to enter any new values.
- Press SET and DOWN at the same time to quit programming or wait one minute and the display will automatically exit programming mode.

## PARAMETER DESCRIPTIONS

**SP** = Set Point-Desired Regulation Temperature.

**r0** = Differential or Hysteresis.

**r1** = Lower Set Point Limit.

**r2** = Higher Set Point Limit.

**d0** = Heating or Cooling Control-regulation cycles only performed, neither defrosting nor continuous cycles exist.

**Heating:** To choose heating control: Set d0=Ht (The output is active when TS1 (temperature of probe) is less than or equal to set point.) TS1<=SP. It then disconnects when TS1>=SP+r0.

**Cooling:** To choose cooling control: Set d0=Co (The output is activated when TS1>=SP+r0.) The display will switch off when TS1<=SP.

**c0** = Minimum Time Between Start to Stop.

**c2** = Load Status during Probe Error. In the event of an open or short-circuited probe, the unit will connect or disconnect the load as defined by this parameter.

**P1** = Probe Calibration. Offset degrees to adjust probe. If the probe is not placed in the exact point that is to be measured, use a standard thermometer and adjust the difference with parameter.

**P5** = Probe Type. Select between J, K, or S type thermocouple.

**P6** = Probe Response Rate (0 = 8 sec, 1 = 4 sec, 2 = 2 sec, 3 = 1 sec.).

**H5** = Access to Probe Parameters.

**A0** = Alarm 1 Hysteresis. The differential associated with A1 parameter.

**A1** = Alarm 1 Threshold. Number of degrees to the working set point that initiates an alarm condition.

**A2** = Alarm 1 Exclusion Time. The amount of time the alarm is disabled from instrument activation.

## PARAMETER DESCRIPTIONS (CONTINUED)

**A3** = Alarm 1 Configuration. Determines the alarm type:

A3= Off; alarm is disabled;

A3= Hi; alarm is activated if the temperature  $\geq$ =SP+A1 and deactivated if  $\leq$ =SP+A1-A0.

A3= Lo; alarm is activated if the temperature  $\leq$ =SP+A1 and deactivated if  $\geq$ =SP+A1-A0.

**A4** = Alarm 2 Hysteresis. The differential associated with A5 parameter.

**A5** = Alarm 2 Threshold. Number of degrees to the working set point that initiates an alarm condition.

**A6** = Alarm 2 Exclusion Time. The amount of time the alarm is disabled from instrument activation

**A7** = Alarm 2 configuration. Determines the alarm type:

A7= Off; alarm is disabled;

A7= Hi; alarm is activated if the temperature  $\geq$ =SP+A5 and deactivated if  $\leq$ =SP+A5-A4.

A7= Lo; alarm is activated if the temperature  $\leq$ =SP+A5 and deactivated if  $\geq$ =SP+A5-A4.

## INDICATORS

**Buzzer**—In the event of alarm or error condition, the internal buzzer is activated. To silence the buzzer, press and hold the SET and DOWN keys.

**LED Indications**—OUT – This indicates the load is connected. The system waits for the programmed minimum stop time of the load.

**Display Messages**— In normal operation, the probe temperature will be shown on the display. In case of alarm or error, the following messages will be shown:

**ER** = Memory Error

**--** = Short-Circuit Probe Error (output determined by c2)

**oo** = Open Probe Error (output determined by c2)

## ACCESS TO ALL CODE PROTECTED PARAMETERS

- Press SET for 8 seconds. The access code value 0 is shown on the display. (Unit comes with code set at 0 from factory).
- With the UP and DOWN arrows, set code to 118. Press SET.
- Move to the desired parameter with the UP and DOWN keys.
- Press SET to view the value on the display.
- The value can be modified with the UP and DOWN arrows.
- Press SET to enter the value and exit to text parameter.
- Repeat until all necessary parameters are modified.
- Press SET and DOWN at the same time to quit programming or wait one minute and the display will automatically exit programming mode.

## WARRANTY INFORMATION

BriskHeat warrants to the original purchaser of this product for the period of eighteen (18) months from date of shipment or twelve (12) months from date of installation, whichever comes first.

BriskHeat's obligation and the exclusive remedy under this warranty shall be limited to the repair or replacement, at BriskHeat's option, of any parts of the product which may prove defective under prescribed use and service following BriskHeat's examination, is determined by BriskHeat to be defective. The complete details of the warranty can be found online at [www.briskheat.com](http://www.briskheat.com) or by contacting us at 1-800-848-7673 (toll free, U.S. / Canada) or 1-614-294-3376 (Worldwide).

Corporate Headquarters:  
4800 Hilton Corporate Dr.  
Columbus, OH 43232, U.S.A.

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