
LEICA ZOOM™ 2000

*Zoom Stereomicroscope
Stereozoom-Mikroskop
Stéréomicroscope avec zoom
Estereomicroscopio con zoom*

*Instruction Manual/Gebrauchsanleitung
Manuel d'emploi/Manual de Instrucciones*

Leica

LEICA ZOOM™2000

Zoom Stereomicroscope

Instruction Manual

Table of Contents

Section	Page
1.0 Introduction	2
2.0 Important Features	2
3.0 Initial Set-Up	4
4.0 Magnification Systems	4
5.0 Illumination Systems	4
6.0 Electrical Considerations/Equipment Ratings	4
6.1 Supply Voltage and Range Supply Frequency	4
7.0 Range of Environmental Conditions	4
8.0 Operating Instructions	6
9.0 Cleaning of Exposed Optical Surfaces	6
10.0 Care of the Stereomicroscope	8
11.0 Lamp Replacement	8
11.1 For Transmitted Illumination	8
11.2 For Reflected Illumination	10
11.3 Fuse Replacement	10

Explanation of Warning Symbols

CAUTION: Risk of electric shock

CAUTION: (Refer to accompanying documents)

 230V Models only: 5 X 20 Typ T 250mA 250V

1.0 INTRODUCTION

The LEICA ZOOM™ 2000 Stereomicroscope from Leica is built to withstand years of heavy classroom use. From biology to botany, entomology to geology, ecology to electronics education, the LEICA ZOOM 2000 has applications wherever three dimensional viewing is needed.

2.0 IMPORTANT FEATURES

All LEICA ZOOM 2000 Stereomicroscopes have zoom magnification systems. This means the stereomicroscope can zoom through an entire magnification range of 7X - 30X or 10.5X - 45X depending on the widefield eyepiece magnification (10X or 15X). Total magnification power is the product of objective magnification multiplied by the eyepiece magnification:

Eyepieces	X	Objectives	=	Zoom Magnification
10X	X	0.7X to 3.0X	=	7X to 30X
15X	X	0.7X to 3.0X	=	10.5X to 45X

Interlocked inclined eyepiece tubes allow adjustment for interpupillary distances from 49mm to 73mm. All eyepieces are built-in to eliminate breakage and loss.

The LEICA ZOOM 2000 comes with a three-way illumination system which may be used to provide transmitted light, reflected light or a combination of both. The built-in, dual on/off switches control the type of illumination desired.

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The rugged, tamper-proof design will meet the challenges of today's varied applications. The glass stage plate is cemented in place to prevent leakage or damage to internal components. The ZOOM 2000's locked-in eyepieces and strong, enclosed telescoping pillars all contribute to the stereomicroscope's image stability and durability.

3.0 INITIAL SET-UP

The ZOOM 2000 comes completely assembled and ready for use. Just follow the step-by-step operating procedure for years of use from your stereomicroscope.

4.0 MAGNIFICATION SYSTEMS

There are two magnification ranges offered with the ZOOM 2000: 7X – 30X and 10.5X – 45X. The zoom control knob allows the user to view specimens through the entire magnification range while remaining in focus.

5.0 ILLUMINATION SYSTEMS

The LEICA ZOOM 2000 is equipped with a built-in, three-way illuminator system. The convenient on/off switches are located on the front of the instrument and control both reflected and transmitted illumination. The illumination system utilizes tungsten and halogen lamps which provide cool, white light for superior crisp image viewing.

6.0 ELECTRICAL CONSIDERATIONS/EQUIPMENT RATINGS

CAUTION: As a safety precaution, the power cord has been grounded to the microscope base. Never use a 2-prong adapter between the power cord and the power source; it will render the grounding feature ineffective.

CAUTION: Removal of the bottom cover of the microscope exposes hazardous voltages, which if contacted, could cause injury or death. Please refer servicing to qualified personnel.

6.1 SUPPLY VOLTAGE AND RANGE SUPPLY FREQUENCY

Supply Voltage	Frequency
120V ~ ±10%	60Hz
230V ~ ±10%	50Hz

7.0 RANGE OF ENVIRONMENTAL CONDITIONS

- Indoor use only
- Temperatures 5°C to 40°C
- Maximum Relative Humidity - 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.
- Installation Category II (Overvoltage Category)
- Pollution Degree 2

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8.0 OPERATING INSTRUCTIONS

Follow these instructions to begin using the LEICA ZOOM 2000 Stereomicroscope.

1. Always use your stereomicroscope on a hard, stable surface.
2. Plug the stereomicroscope power cord into a suitable grounded electrical outlet. A grounded 3-wire cord is provided.
3. Select the appropriate light switch for the specimen. The ZOOM 2000 is equipped with three types of illumination:
 - Transmitted Light - to view translucent specimens;
 - Reflected Light - to view opaque specimens;
 - Combination of both - to develop contrasting conditions on any variety of specimens.
4. Place a specimen on the circular stage area directly below the center of the stereomicroscope body. While looking through the eyepieces with both eyes, slowly adjust the distance between the eyepiece tubes until you see a complete circle of light at the same time without moving your head. Eyepiece tubes are inclined at the proper angle for comfortable viewing and do not require individual focusing. They are factory adjusted for normal or corrected-to-normal vision.
5. Turn the zoom control knob until the lowest number on the knob is opposite the index mark on the body. Now, while looking through the eyepieces, bring the specimen into sharp focus by turning the focusing knobs. The distance is about four inches from the specimen to the bottom of the microscope body when focus is obtained.

In practice, it is always best to focus the zoom model microscopes at the highest power and then turn the zoom control knob to the lower powers. This will keep the specimen in close focus at any magnification. A slight turn of the focus knobs may be needed to sharpen the focus.

9.0 CLEANING OF EXPOSED OPTICAL SURFACES

Do not attempt to disassemble the inclined binocular body or focusing mechanisms. These have been factory aligned and sealed for protection against dirt, dust or other particles. Only exposed surfaces should require cleaning. To clean the stereomicroscope:

- Blow dust particles from the lens surface with an ear syringe.
- Brush the surface with a clean camel's hair brush.
- Clean the lenses and eyepieces with moist lens paper and cotton swabs. This technique also may be used to remove oily smears and fingerprints. Such smears detract from the image quality.

Clean glass surfaces very carefully and only when necessary.

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10.0 CARE OF THE STEREOMICROSCOPE

The ZOOM 2000 stereomicroscope is a precision instrument and should be handled with care. Keep the following tips in mind while using your stereomicroscope.

1. Always carry the stereomicroscope with two hands.
2. Clean the lenses and eyepieces with moist lens paper or a cotton swab.
3. When finished using the stereomicroscope, return the body to the lowest position, turn off the illuminator and place the dust cover over the instrument.
4. The ZOOM 2000 has a durable finish. All non-optical parts of the instrument can be cleaned with methanol, naphtha or soap and water. Avoid the use of all other organic solvents when cleaning the instrument.

11.0 LAMP REPLACEMENT

If lamp replacement is necessary, follow these easy steps:

DANGER: The glass envelope of the lamp may be extremely hot. Do not attempt to change the lamp before it is completely cooled or without wearing adequate skin protection.

11.1 FOR TRANSMITTED ILLUMINATION

CAUTION: Lamp is not user serviceable, please refer service to qualified service person.

1. Unplug the electrical cord.
 2. Carefully lay the instrument on its back.
 3. Unscrew the base assembly.
 4. Unscrew the burned-out lamp, after it has cooled.
 5. Screw in the new lamp and reattach the base.
 6. Never operate the instrument unless the base is securely attached.
- Replacement lamps, Part Number 313125 (115V) or 313126 (230V) can be purchased from a Leica dealer.

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11.2 FOR REFLECTED ILLUMINATION

1. Unplug the electrical cord.
2. Raise the body to the highest position by turning the focusing mechanism.
3. Unscrew the lower portion of the illuminator housing being careful to keep the lower portion intact (Figure 1). Set aside.
4. Carefully disconnect the reflector lamp from the connector socket. Replace the burned out lamp with a new lamp making sure not to bend the pins when inserting it (Figure 2).
5. Reattach the lower portion of the illuminator housing by screwing it into the existing upper portion.
6. Rotate the chrome safety latch into the original position to secure the illuminator housing.
7. Never operate the stereomicroscope illuminator unless the lamp illuminator housing is securely in place.

Replacement lamp Part Number 630 (Phillips P/N 13528) can be purchased from a Leica dealer.

11.3 FUSE REPLACEMENT (230V MODELS ONLY)

CAUTION: FUSE REPLACEMENT - THE FUSE IS NOT USER SERVICEABLE. IT MUST BE REPLACED BY QUALIFIED SERVICE PERSONNEL. THE FUSE CAN BE REPLACED ONLY WITH THE SAME SIZE, TYPE AND RATING OF FUSE.

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Leica Microsystems Inc.