

Table of Contents

1	Introduction	5
2	Safety Notes	6
	2.1	Definition of signal warnings and symbols 6
	2.2	Product Specific Safety Notes 6
3	Design and Function	8
	3.1	Overview 8
	3.2	Connections on the instrument 9
	3.3	User interface 9
	3.3.1	Home screen 9
	3.3.2	Icons and buttons 10
4	Putting into operation	11
	4.1	Scope of Delivery 11
	4.2	Installing power supply 12
	4.3	Installing EasyStir 13
	4.4	Installing sensors 14
	4.5	Installing tubes 15
	4.6	Installing printer & balance 16
5	Setup & Tools	17
	5.1	Settings 17
	5.2	Language 17
	5.3	Diagnostics 17
	5.4	Toolbox 17
6	Functional Apps on homescreen	19
	6.1	Na Measure 19
	6.2	Burette 20
	6.3	Stirrer 20
	6.4	Standard 20
	6.5	Sensor 20
7	Operation of the instrument	22
	7.1	Used standards 22
	7.2	Sample preparation 23
	7.3	Performing sensor conditioning 24
	7.4	Performing an analysis 24
	7.4.1	Single determination 24
	7.4.2	Multiple determination 25
	7.5	Results 26
8	Maintenance	28
	8.1	Maintaining and changing the burette 28
	8.2	Cleaning the burette parts 29
	8.3	Maintaining the valve 29
	8.4	Transporting the instrument 30
	8.5	Cleaning the housing 30
	8.6	Disposal 30
9	Accessories	31

10	Technical Data	34
11	Declaration of Conformity	35

1 Introduction

The Sodium Analyzer has been especially designed for simplicity in the determination of sodium in food product using the standard addition technique. Just add your sample and ionic strength adjustment (ISA) solution then start your measurement. No calibration of the system is needed. Save operator time with reduced efforts in sample preparation and calibration.

Increase productivity with fast measurements of your samples. The Sodium Analyzer is modern with a small footprint for use in a wide variety of sodium determination applications. The instrument has been developed primarily for the use in Food and Beverage Labs. The Sodium Analyzer perfectly combines simple, easy-to-understand operation with high precision and reliability. Thanks to the App-oriented user interface and the built-in know-how, operation could not be any simpler or more intuitive. The Sodium Analyzer can be controlled by touch screen and all measured data can be stored in the PC Software EasyDirect.

All main functions can be started directly from the home screen on the touch screen via Longclick™ on the corresponding App, which makes routine use extremely simple. The integrated, proven algorithm specifically designed for Na⁺ determination delivers highly accurate and repeatable results and is dedicated to determining directly the harmful sodium in Food and Beverage products.

Find out more and visit us on:

www.mt.com/SodiumAnalyzer

If you have any additional questions, METTLER TOLEDO is always available to assist you.

2 Safety Notes

2.1 Definition of signal warnings and symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

Signal words

WARNING	for a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	for a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
Attention	(no symbol) for important information about the product.
Note	(no symbol) for useful information about the product.

Warning symbols



General hazard



Electrical shock



Toxic substance



Inflammable or explosive substance



Acid / Corrosion

2.2 Product Specific Safety Notes

Your instrument represents state-of-the-art technology and complies with all recognized safety rules, however, certain hazards may arise in extraneous circumstances. Do not open the housing of the instrument; it does not contain any parts that can be maintained, repaired or replaced by the user. If you ever have problems with your instrument, contact your authorized METTLER TOLEDO dealer or service representative.

Intended use



This instrument is designed to be used in analytical laboratories and is suitable for the processing of reagents and solvents.

The use therefore requires knowledge and experience in working with toxic and caustic substances as well as knowledge and experience working with application-specific reagents, which may be toxic or hazardous.

The manufacturer shall not be held liable for any damage resulting from incorrect usage divergent to the operating instructions. Furthermore, the manufacturer's technical specifications and limits must be adhered to at all times and in no way exceeded.

Location



The instrument has been developed for indoor operation and may not be used in explosive environments.

Place the instrument in a location which is suitable for the operation, protected from direct sunlight and corrosive gases. Avoid powerful vibrations, excessive temperature fluctuations and temperatures below 5 °C and above 40 °C.

Protective Clothing

It is advisable to wear protective clothing in the laboratory when working with hazardous or toxic substances.



A lab coat should be worn.



Suitable eye protection such as goggles should be worn.



Use appropriate gloves when handling chemicals or hazardous substances, checking their integrity before use.

Safety notes



WARNING

Risk of electric shock

Use only 3-pin grounded electrical outlet and extension cables to connect the instrument.

- a) Only 3-pin grounded electrical outlet and extension cables for connecting your instrument must be used.
 - b) Intentional disconnection of the equipment grounding conductor is prohibited.
-



WARNING

Risk of corrosion

Leaks in tubing connections and loose titration vessels are a safety risk.

- a) Tighten all connections well by hand, avoid applying excessive force to tubing connections.
 - b) Always test the titration vessel for firm seating in the titration head.
-



WARNING

Flammable solvents

All relevant safety measures must be observed when working with flammable solvents and chemicals.

- a) Keep all sources of flame away from the workplace.
 - b) When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules.
-



WARNING


Chemicals

All relevant safety measures are to be observed when working with chemicals.

- a) Set up the instrument in a well-ventilated location.
 - b) Any spills should be wiped off immediately.
 - c) When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules.
-

3 Design and Function


3.1 Overview



1	On/Off button	2	Status LED
3	Touch screen	4	Valve
5	Burette	6	Release button
7	Analyzer head Na	8	DX222-Na (Na-ISE sensor)
9	DX205-SC (Reference sensor)	10	Beaker
11	Stirrer (EasyStir)		

3.2 Connections on the instrument

Rear view




1	Socket for stirrer (Mini-DIN)	2	Socket for temperature sensor (RCA - Cinch)
3	Socket for sensors	4	USB B interface for PC connection (EasyDirect PC Software)
5	USB A interface for printer and balance	6	DC socket for power supply

3.3 User interface

3.3.1 Home screen

The home screen is the main screen and appears after startup of the instrument. Tap and hold any of the home screen apps (except Standard) to start the last analysis without further navigation (LongClick™).



1	Functional Apps: Apps for various functions.	2	Results: Select this menu item to display the results of the previous analysis.
3	Setup & Tools: Select this menu item to make changes to the system settings and to perform diagnostics.		

3.3.2 Icons and buttons



Tap this menu item to return to the home screen.



Tap this menu item to start an action.



Tap the arrows to page through parameter sets.



Tap this icon to accept and close any entry screen.



Tap this icon to print parameters or results.



Tap this icon on a parameter screen to change the parameter set. Parameter set with number two is active in this example.



Tap this menu item to go back to the previous screen.



Tap this menu item to stop a running action.



Tap the help icon to change to help mode. Then tap any of the parameters to get a specific help description.



Tap this icon to reject and close any entry screen.


















Indicates which analysis will be started when you tap and hold an app (LongClick™). Parameter set with number two is active in this example.

4 Putting into operation

Before getting started, please make sure you have received all parts listed in the section scope of delivery.
Read the following sections carefully for a correct and faultless installation.

4.1 Scope of Delivery

Check the completeness of the delivery. The following accessories are part of the standard equipment of your new instrument:

	Sodium Analyzer		EasyStir
	Burette 10 mL		Magnetic stirrer bar
	Analyzer head Na		Bottle head incl. flat seal
	Beaker PP (100 mL, 10 pcs.)		DX222-Na (Na-ISE sensor)
	Cable reference sensor Na		DX205-SC (Reference sensor)
	Tubing set instrument		AC/DC power adapter
	Country specific power cable		Application brochure
	EasyDirect Titration Trial Version		


4.2 Installing power supply

The instrument is supplied with a universal AC adapter. The AC adapter is suitable for all line voltages in the range of 100 to 240 V, 50/60 Hz.

Attention

- Before installing, check cables for damage.
- Only 3-pin grounded electrical outlet and extension cables for connecting your instrument must be used.
- Ensure the cables are arranged so that they cannot be damaged or interfere with the operation.
- Take care that the AC adapter does not come into contact with liquids.
- The power plug must be accessible at all times.

Procedure




- 1 Connect the plug of the AC adapter with the DC socket of the instrument.
- 2 Secure the plug by firmly tightening the knurled nut.
- 3 Connect the 3-pin grounded power cable to the AC adapter and then connect the power cable to the electrical outlet.
- 4 Push the **On/Off** button to switch on the instrument. The LED next to the button flashes as the system starts up and then remains permanently lit.


4.3 Installing EasyStir

This stirrer is powered by the instrument and will automatically be switched on/off according to the settings.

- 1 Connect EasyStir to the instrument stirrer socket, observing the arrow on the connector (arrow must be top-side).




- 2 Push the blue release button to install the analyzer head on stirrer rod.



4.4 Installing sensors

The Easy Na works with two sensors, a reference sensor (DX205-SC) and a Na-ISE sensor (DX222-Na). Both sensors need to be installed and connected to the instrument as follows:



1	BNC plug	2	Socket for adapter cable for sensors
3	Reference sensor	4	Plug for Na-ISE sensor (dark blue-colored)
5	Plug for Reference sensor (cyan-colored)	6	Analyzer head Na
7	DX222-Na (Na-ISE sensor)	8	DX205-SC (Reference sensor)


- 1 Insert Reference sensor into the cyan-colored rubber ring of the Analyzer head Na.
- 2 Insert Na-ISE sensor into the dark blue-colored rubber ring of the Analyzer head Na.
- 3 Connect the cyan-colored plug to the Reference sensor.
- 4 Secure the plug by firmly tightening the knurled nut.
- 5 Connect the dark blue-colored plug to the Na-ISE sensor.
- 6 Secure the plug by firmly tightening the knurled nut.
- 7 Connect the BNC plug to the socket for adapter cable for sensors, observing the nuts on the BNC plug.
- 8 Secure the BNC plug by turning the knurled nut a 1/4 turn.

4.5 Installing tubes



Note

The tube connectors of the tubes are marked in different colors (green, red, cyan).
An incorrect connection may lead to a malfunction!



1	Standard bottle	2	Burette
3	Drying tube	4	Tube burette-valve
5	Suction tube for standard	6	Tube titration
7	Cap nut burette	8	Cap nut valve

- 1 Screw the cap nut valve on tightly.
- 2 Screw the cap nut burette on tightly.
- 3 Connect Suction tube for standard to standard bottle and the cyan-colored valve inlet.
- 4 Connect burette-valve tube to burette and the red-colored valve inlet.

- 5 Connect titration tube to the green-colored valve outlet and the green-colored rubber ring on the Analyzer head Na.
- 6 Screw on all tube connectors tightly.
- 7 Ensure that all tubes are firmly connected.

4.6 Installing printer & balance


Printers and balances can be connected to the USB interface found on the back of the instrument. Printers and balances are recognized automatically when connected. They can be used immediately by the instrument without any special settings.

Supported devices


Manufacturer	Type	Model
METTLER TOLEDO	Balance	MS/ML (NewClassic)
METTLER TOLEDO	Printer	USB-P25

A standard USB-hub can be used if more than one device is to be connected to the USB port of the instrument.

Connecting printer or balance:



Connecting printer and balance, using a hub:



5 Setup & Tools



Setup & Tools

Select this menu item on the home screen to make changes to the system settings and to perform diagnostics.

5.1 Settings



Navigation: **Setup & Tools > Settings**

On tab **Settings** you can define general settings such as date-/time format, temperature unit or blank value.

On tab **System information** you will find general information about the instrument's hardware and software.

5.2 Language



Navigation: **Setup & Tools > Language**

Defines the user interface language.

5.3 Diagnostics



Navigation: **Setup & Tools > Diagnostics**

Performs diagnostics for the peripheral devices and the printer.

Available functions



Peripheral check

A connection test of the available peripheral devices is performed.



Printer check

A test printout is generated on the connected printer.

5.4 Toolbox



Navigation: **Setup & Tools > Toolbox**

This screen contains different tools for maintaining the firmware, adjusting the screen or performing a factory reset.

Available functions



Factory reset

This function will reset all data and settings of the instrument.



FW Update

This function will update the firmware of the instrument. The following firmware can be updated:

- **FW instrument:** This is the application firmware.
- **FW mainboard:** This is the firmware of the micro-controller.



Adjust screen

This function will start the touch screen adjustment.

Adjustment of the touch screen is necessary when your touch screen responds inaccurately e.g. after a firmware update has been performed.

6 Functional Apps on homescreen

6.1 Na Measure



Na Measure

Select this menu item to define analysis parameters and to start a measurement

Parameters

Only the most important parameters are listed here. To obtain additional information on all other parameters, you may refer to the online help by clicking [?] on the instrument.

Parameters	Description	Values
Potential diff. dE	Defines the targeted potential difference for each standard addition. Standard will be added until the potential difference is reached. Please make sure that (dE) x (No. of additions) is between 30 and 60 mV. Please refer the application notes to obtain a suggestion on the optimum values for your sample.	1.0 - 60.0
Analysis type	Selects the type of the analysis. For most applications Direct can be used. If your sample preparation contaminates your solution with Na ⁺ (e.g. acid hydrolysis), determine the blank value first and use Blank compensated afterwards to obtain more accurate results.	Direct Blank determination Blank compensated
Control	Definition of internally fixed parameters for the value acquisition.	Normal Fast Cautious User defined
Sampling	Select this option to choose the way the sample is prepared for the analysis. Direct Preparation of the sample directly in the beaker. Note: For solid samples, the sample volume is not taken into account ($V_{\text{solvent}} = V_{\text{water}} + V_{\text{ISA}}$). This sample preparation type is therefore recommended for solid samples which do not dissolve in the solvent. Liquid samples are assumed to be soluble in the solvent and, therefore, the sample volume is taken into account ($V_{\text{solvent}} = V_{\text{water}} + V_{\text{ISA}} + V_{\text{sample}}$). Aliquot incl. ISA Select "Aliquot incl. ISA" if you create a dilution solution with your sample, water and ISA. Take a known amount (aliquot) out of your dilution solution and put it in the analysis beaker. Aliquot excl. ISA Select "Aliquot excl. ISA" if you create a dilution solution with your sample and water (without ISA!). Take a known amount (aliquot) out of your dilution solution and put it in the analysis beaker. Then add the ISA directly to the beaker.	Direct Aliquot incl. ISA Aliquot excl. ISA
Multiple determination	This option will calculate the mean value and standard deviation (relative and absolute) in the samples. Maximum number of samples: 5.	Yes No
Report	None No Report is printed at the end of the analysis. Short A summary is printed at the end of the analysis. Long A full report is printed at the end of the analysis.	None Short Long

6.2 Burette



Burette

Select this menu item to rinse the burette or to dispense a defined quantity of standard. Change the burette size according the actual installed burette.

Properties

Parameters	Description	Values
Burette size	Select the size of burette you are using. The standard size of the burette is 10 mL. The adaption of a 20 ml burette is also possible.	10 mL 20 mL

Actions

Parameters	Description	Values
Rinse	Define the rinse cycles of the burette. A maximum of 3 rinse cycles is normally sufficient.	1 - 20
Dispense	Define the dispensing volume of the burette.	0.1 - 100.0

6.3 Stirrer



Stirrer

Select this menu item to switch the stirrer on or off at a definable stirring speed.

Parameters	Description	Values
Stir speed	Select the required stirrer speed.	Low Medium High Maximum

6.4 Standard



Standard

Select this menu item to define the standard properties. It is possible to define 3 standards (**S1, S2, S3**).

Parameters

Parameters	Description	Values
Standard type	The normally used standard is Na ⁺ (Sodium). Defining a user-specific standard name is also possible.	Sodium User defined
Actual concentration	Enter the Actual concentration of Na ⁺ of the used standard.	0.01 - 400'000.00

6.5 Sensor



Sensor

Select this menu item to view the sensor properties and to start a sensor conditioning. In the two tabs **Properties|Conditioning** you will find the sensor properties and the parameters to execute a sensor conditioning. If no temperature sensor is defined for the conditioning, the system settings are taken into account for temperature compensation (Global temperature).

Properties

Parameters	Description	Values
Sensor	Name of the measuring half-cell. The Mettler Toledo DX222-Na is part of delivery.	-
Reference sensor	Name of the reference electrode. The Mettler Toledo DX205-SC is part of delivery.	-

Temperature sensor	Define whether you are working with a temperature sensor or not. Mettler Toledo TEMPERATURE PROBE NTC 30K is optional available. Please note that temperature sensor must be activated if the temperature sensor is connected.	On Off
Cond. Date/Time	Shows the date and time of the last sensor conditioning.	-
Cond. result	Shows the status of the last sensor conditioning.	-

Sensor conditioning

We recommend that daily sensor conditioning is performed before starting measurements.

See also

- Performing sensor conditioning (page 24)

7 Operation of the instrument



Note

Risk of spilling harmful chemicals!

While working on the instrument wear safety goggles, a lab coat and suitable gloves at any time!

Attention:

Before starting any measurements with the instrument please ensure the startup procedure has been completed. Read the following sections carefully for a correct and faultless analysis.

7.1 Used standards

Recommended ISA solution 1 M

- ▶ Ideally, prepare the solutions in a hood or a well-ventilated room.
 - ▶ ISA (Ionic strength adjuster) solution is added to the preconditioning solution, to the standard sodium solution as well as to the sample solution.
 - ▶ The final ISA-concentration in the analysis beaker should be 0.5 M ISA.
- 1 Prepare a 1'000 mL volumetric flask filled with approximately 140 mL deion. H₂O (dH₂O).
 - 2 Add 724 mL 0.5 M HCl and agitate to mix.
 - 3 Add 140.15 mL concentrated diisopropylamine (DIPA, 99.5 %).
 - 4 Fill dH₂O up to the 1'000 mL volumetric mark. The DIPA will dissolve slowly under intensive agitation.

Note

- The ratio of standard or sample to ISA should be 1:1.

Standard sodium solution with 0.5 M ISA

The standard solution 2'000 ppm Na Std. with 0.5 mol/L DIPA is prepared as follows:

- 1 Use 5.084000 g NaCl for 1 L of solution.
- 2 Dissolve the solution in a volumetric flask using 400 mL dH₂O + 500 mL 1 M ISA (see above).
- 3 Fill with dH₂O up to the volumetric mark 1000 mL.

The standard solution 5'000 ppm Na Std within 0.5 mol/L DIPA is prepared as follows:

- 1 Use 12.71000 g NaCl for 1 L solution.
- 2 Dissolve the solution in a volumetric flask using 400 mL deion H₂O + 500 mL 1 M ISA (see above).
- 3 Fill with dH₂O up to the volumetric mark 1000 mL.

Note

- The ratio of standard or sample to ISA should be 1:1.

Electrolyte for the reference sensor


- Use 0.5 M ISA solution.

7.2 Sample preparation

You can prepare samples in different ways. The used values in the descriptions are only for example.


Analysis type Direct:

- 1 Add 5 g (solid sample) or 10 mL (liquid sample) into the analysis beaker.
- 2 Add 20 mL of water.
- 3 Add 20 mL of 1 M ISA.




Analysis type Aliquot incl. ISA:

- 1 Add 5 g (solid sample) or 10 mL (liquid sample) into a volumetric flask.
- 2 Fill up to the volumetric mark e.g. 1000 mL with 0.5 M ISA.
- 3 Take out 40 mL aliquot and fill in the analysis beaker.



Analysis type Aliquot excl. ISA:

- 1 Add 5 g (solid sample) or 10 mL (liquid sample) into a 1000 mL volumetric vessel.
- 2 Fill up with water exactly to the volumetric mark of 1000 mL.
- 3 Take out 20 mL aliquot and fill in the analysis beaker.
- 4 Add 20 mL of 1 M ISA.



Note

- Detailed information about specific applications can be found in our application brochure which is part of the delivery.

7.3 Performing sensor conditioning



Note

Damage to the Sensors!

The Na-ISE sensor and the reference sensor will be damaged if they are allowed to dry out. Make sure that both sensors are kept in liquid at all times. For detailed information see the sensor manuals.

Before starting your daily measurements, the sensors have to be conditioned. The sensor conditioning is a pre-conditioning with an internal test to check that the sensor values are stable. The instrument measures the voltage and verifies that the value is stable over a minimum time of 60 s (The maximum deviation of the measured values is 0.1 mV). Sensor conditioning will stop automatically when the sensors are ready to use. If the sensor values are not stable, sensor conditioning will stop automatically after 30 min with the result = failed.

- 1 Remove the Reference sensor from the Analyzer head Na.
 - 2 Check by hand that the diaphragm moves smoothly. Rinse diaphragm with distilled water.
 - 3 Open the rubber cap of the filling aperture and fill up the Reference sensor with electrolyte solution (0.5 M DIPA-HCL-ISA).
 - ⇒ The rubber cap must be open during sensor conditioning and measurements!
 - 4 Put the Reference sensor back into the Analyzer head Na.
 - 5 Put the Reference sensor into a beaker with Na⁺ solution (including ISA).
 - ⇒ Make sure the diaphragm is fully immersed into the conditioning solution.
 - 6 Make sure the Na-ISE sensor is immersed into the conditioning solution.
 - 7 Tap [🟢] to start sensor conditioning.
 - ⇒ Sensor conditioning will normally take about max. 5 min and stops automatically.
 - ⇒ Tap [🔴] at any time to stop sensor conditioning.
 - ⇒ A timeout will be set automatically after 30 min (result = failed).
- ⇒ **Result = Passed**, the sensors are ready for use (Threshold = 0.1 mV).
- ⇒ **Result = Failed**, please check the sleeve diaphragm and the electrolyte level of the reference sensor. Make sure no air bubbles are inside the electrolyte of the reference sensor. Check the cable connections and repeat sensor conditioning.

Note

- We recommend that a daily sensor conditioning is performed before starting measurements.

7.4 Performing an analysis

Attention

- Before performing an analysis, ensure that all tube connections are firmly tightened.
- Make sure that the analysis tube is directed into an analysis beaker that is a multiple of the volume of the burette.
- Rinse the burette and tubes with the function **Burette**. Perform this function each day before the first analysis is run.

Note

- Rinse the burette until no more air bubbles are present in the tube connections. Air bubbles may effect the measurement results!

7.4.1 Single determination

- ▶ The sample is prepared and added to the beaker.
- ▶ The magnetic stirrer bar is placed in the beaker.

- ▶ All parameters for the method are checked and defined. **Multiple determination = No.**
- 1 Place the beaker on the center of the stirrer unit.
 - ⇒ Verify that the magnetic stirrer bar is turning correctly.
- 2 On the home screen, tap and hold [🟩] to start the analysis.
 - ⇒ You will be prompted if any parameters are required by the analysis.
- 3 Follow the instructions on the screen. With the progress of the measurement, the displayed curve is automatically rescaled so that the entire measurement is visible.
 - ⇒ The result screen is displayed.

Note

- Tap [🟢] to stop the measurement for early data transfer. The actual results will be shown in the menu **Results**
- Tap [🔴] to stop the measurement without storing data.
- If an error occurs, measurement will stop immediately.

7.4.2 Multiple determination

- ▶ The sample is prepared and added to the beaker.
- ▶ The magnetic stirrer bar is placed in the beakers.
- ▶ All parameters for the method are checked and defined. **Multiple determination = Yes.**
- 1 Place the beaker on the center of the stirrer unit.
 - ⇒ Verify that the magnetic stirrer bar is turning correctly.
- 2 On the home screen, tap and hold [🟩] to start the analysis.
 - ⇒ You will be prompted if any parameters are required by the analysis.
- 3 Follow the instructions on the screen. With the progress of the measurement, the displayed curve is automatically rescaled so that the entire measurement is visible.
 - ⇒ When the first analysis is finished, the result screen is displayed with the two options **Last sampleOverview**.
- 4 Prepare the next sample and place the beaker on the stirrer.
- 5 Tap [🟢] to start the next analysis and follow the instructions on the screen.
 - ⇒ Repeat these steps for a maximum of 5 samples.
- 6 To end the series before the maximum of 5 samples is reached, tap [🔵]
- ⇒ The result screen is displayed with the options **Last sampleOverview**. The results of the last sample and the overview including the statistical data for the sample series are displayed.


Note

- Tap [🟢] to stop the measurement for early data transfer. The current results will be shown in the menu **Results**
- Tap [🔴] to stop the measurement without storing data. Accumulated data will automatically be excluded from the statistics.
- If an error occurs, measurement will stop immediately and accumulated data will also automatically be excluded from the statistics.

7.5 Results

Results of a single determination

Content (R1)	126.08 ppm	Na Measure M1
Content (R2)	0.126 g/L	
Sample Size	23 mL	
Volume	3.258 mL	
Number of additions	3 (fixed)	
Slope	55.7 mV/log(c)	
Coeff. of determ. (R ²)	0.99997717	
Date/Time	12/08/2013 13:47	

← 1/2 → 

All available results for the sample are shown.

Content(R1)	Entered value in Calculation 1
Content(R2)	Entered value in Calculation 2
Sample size	Entered value in Sample size (mL) Only for Sampling = Aliquot incl. ISA/Aliquot excl. ISA : The sample size is the amount of sample in the dilution volume.
Sample size in aliquot	Only for Sampling = Aliquot incl. ISA/Aliquot excl. ISA : The sample size aliquot is the amount of the sample which is currently in the measuring beaker.
Volume	Total volume of the standard added to the sample.
Number of additions	Number of additions performed. Fixed Terminated according to additions in method Manual Manually terminated additions [✔]
Slope	When measuring sodium ions, the slope factor at 298 K (25°C) has a value of 59.16 mV. This is termed the Ideal Slope Factor and means that for each tenfold change in Sodium concentration, an ideal measuring system will sense a mV change of 59.16. The measurement of slope factor gives an indication of the performance of the sensor.
Coeff. of determ. (R²)	R ² denotes the strength of the relationship between the concentration (logarithmic) and the mV signal and indicates how well data points fit the linear regression. R ² represents the percent of the data that is the closest to the line of best fit. For example, if R ² = 0.950, this means that 95% of the total variation in y can be explained by the linear relationship between x and y.
Date/Time	Date and time of the displayed analysis-results.

Results of a multiple determination

No.	R1 [%]	R2 [ppm]		Na Measure M1
1	9.98	1.829	✓	
2	9.99	1.832	✓	
3	10.02	1.842	✗	
4	10.01	1.839	✓	
RSD1 [%]		AVG1	AVG2	
9.99		0.511	0.2	

1 —

The results of the individual samples are displayed as well as the statistics (Average, absolute and relative standard deviation). By tapping on the individual sample, all available results for the respective sample are shown. If two calculations are defined in the method type, tap on the bottom bar (1) to get all statistics listed.



Tap this icon to exclude a result from the statistics.



Tap this icon to include a result, which was previously excluded.

Results of stopped or erroneous analyses are excluded automatically.

8 Maintenance



Note

Risk of spilling harmful chemicals!

While cleaning the instrument, wear safety goggles, a lab coat and suitable gloves at any time!



Note

Before maintaining any parts of the instrument, the following safety precautions must be observed:

Rinse and empty the burette, valve and tubes using dH₂O.

Switch off the instrument and unplug the AC adapter from the instrument.


8.1 Maintaining and changing the burette



Note

Damage to the O-rings!

Do not pull the piston out of the glass cylinder, reinserting the piston will damage the O-rings!



Disassembling the burette cylinder

- ▶ The burette, valve and tubes are rinsed and emptied.
 - ▶ The AC adapter is unplugged from the instrument.
- 1 Unscrew the connection tube on top of the burette cap and clean the connections with a tissue.
 - 2 Unscrew the blue cap nut.
 - 3 Lift off the burette cover including cap nut.
 - 4 Carefully lift the glass cylinder, until you are able to slide off the burette from the piston rod.

⇒ The burette can now be exchanged or reused.

Assembling the burette cylinder

– Reassemble the burette in reverse order.

8.2 Cleaning the burette parts

Depending on the frequency of usage, you should clean the burette cylinder, piston, valve and tubing relatively often. It is important to use high quality ethanol for the cleaning procedure.

- 1 Depending on the contamination caused by the standard, rinse cylinder, valve and tubes with deionized H₂O then with ethanol.
- 2 Dry the parts with oil-free compressed air.

Note

- Never place O-rings in organic solvents.
- Never attempt to remove any crystals in the cylinder by scratching with a hard object. Pipe cleaners or Q tips™ are more suitable.
- Never put the parts in a drying oven with a temperature higher than 40 °C.
- Replace the burette if the piston leaks or is badly scored at the edge. Pay special attention to crystal formation between the seals of the piston if you work with certain solutions.

8.3 Maintaining the valve

Cycle counter


The instrument permanently counts the cycles of the valve. When the life time of 5000 cycles is reached, a message box will open. At this time it is recommended that the valve is exchanged. It is possible to reset the cycle counter manually by using the function **Reset cycles valve** in the menu **Setup & Tools**. Use this function only if the valve has been replaced before reaching 5000 cycles.

- Tap [**Reset cycles**] to reset the valve cycles to 0 (when you have exchanged the valve).
– or –
Tap [**Continue**] to continue working with the valve.

The counter can be viewed at any time in:

Setup & Tools > Settings > No. cycles valve.

Disassembling the valve



- ▶ The burette, valve and tubes are rinsed and emptied.
 - ▶ The AC adapter is unplugged from the instrument.
 - 1 Unscrew the tube connections on top of the valve and clean the tube connections with a tissue.
Note Ensure that all tube connections are unscrewed before losing the cap nut!
 - 2 Unscrew the cap nut, holding the valve.
 - 3 Remove the valve.
- ⇒ The valve can now be exchanged or reused.

Assembling the valve

- 1 Insert the valve, observing the correct positions of the driver bar and the guide.
- 2 Fasten the cap nut.
- 3 Clean the tube connections with a tissue.
- 4 Screw the tube connections to the valve.

8.4 Transporting the instrument

Note the following instructions when transporting the instrument to a new location.

- ▶ The burette, valve and tubes are rinsed and emptied.
 - ▶ The AC adapter is unplugged from the instrument.
- 1 Remove all tube connections from burette and valve.
 - 2 Remove all cable connections from the instrument.
 - 3 Remove the burette as described above.
- ⇒ The instrument is ready to be transported.

Procedure if the burette and valve cannot be emptied



Note

Risk of spilling harmful chemicals!

Do not remove any tubes attached to the burette or valve. Harmful chemicals may leak during transport.

- 1 Remove all cable connections from the instrument.
 - 2 Carefully remove the burette as described above, without removing the tube.
 - 3 Carefully remove the valve as described above, without removing any tubes.
 - 4 Place burette and valve in a suitable container for storage or further handling.
- ⇒ The instrument is ready to be transported.

8.5 Cleaning the housing



Note

Damage to the instrument!

Ensure that no liquid enters the interior of the instrument.

Wipe off any spills immediately.

The housing is made of Polypropylene (PP GF30). This material is sensitive to certain acids and organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK).

- Clean the housing of the instrument using a cloth dampened with water and a mild detergent.

8.6 Disposal










In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.


Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.




9 Accessories

	Description	Order number
	Tubing set instrument	30065464
	EasyStir GT	30065467
	Analyzer Head Na	30079612
	Insert set Analyzer head Na	30079618
	Valve	30042860
	Glass titration vessels, 100 mL (set of 20 pcs.)	101446
	Titration vessels polypropylene, 80 mL (set of 120 pcs.)	51109388
	Burette 10 mL Burette 20 mL	30043900 30043901
	Siphon tips (5 pcs.)	23240




Burette cover set

30057633




DX222-Na

30079616




DX205-SC

30066675




Temperature probe NTC 30K

51300164



Bottle Head incl. flat seal


30060023



Bottle adapter
Merck, Germany
Fisher, USA

23774

23787




Bottle adaptor China (D34) incl. flat seal
Bottle adaptor China (D28.5) incl. flat seal
Bottle adaptor India incl. flat seal

30079456


30060026

30060024




USB-P25 compact printer (including USB-cable)

USB-P25




AC/DC power adapter

51105795



Cable Reference sensor Na

30079614



EasyDirect Titration Trial Version

30096988

10 Technical Data

Sodium Analyzer

Interfaces	USB A	USB full / low speed
	USB B	USB full / low speed
Stirrer output	Voltage	0 - 9 V DC
	Socket	4-pin Mini-DIN
Sensor mV Input	Measuring range	± 2000 mV
	Socket	BNC
Temperature Input	Sensor type	NTC 30 KΩ at 25 °C
	Measuring range	0 - 100 °C
	Socket	Cinch (RCA)
Touch screen	Display	4.3 inch color, 480 x 272 pixel
	Input Technology	Full-coverage touch screen
Ambient conditions	Ambient temperature	5 °C - 40 °C
	Relative humidity	Max 80% (non condensing) at 31 °C, linearly descending to 50% at 40 °C
	Overvoltage category	Class II
	Pollution degree	2
	Range of application	For indoor use only
	Maximum operating altitude	Up to 2000 m
Dimensions	Width	170 mm
	Depth	220 mm
	Height	350 mm
Weight	Instrument	1850 g
	Stirrer	800 g
Power rating instrument	Input voltage	24 V $\overline{=}$
	Input current	1.25 A
Power rating AC adapter	Line voltage	100 - 240 V $\sim \pm 10 \%$
	Input frequency	50/60 Hz
	Input current	0.8 A
	Output voltage	24 V $\overline{=}$
	Output current	1.25 A
Materials	Housing	PP GF30
	Metal parts	Stainless steel
	Touch screen cover	Polyester

EasyStir

Stirrer motor	Motor type	DC
	Voltage	0 - 9 V DC
	Cable connection	4-pin Mini-DIN
Materials	Housing	PP GF30
	Metal parts	Stainless steel

11 Declaration of Conformity

EC - DECLARATION OF CONFORMITY

EG-Konformitätserklärung

KD-Nr.: 30065475

Doku-Nr.: 20120038

The undersigned, representing the following manufacturer

Die Unterzeichnenden vertreten das folgende Unternehmen

Mettler-Toledo AG (MTANA)

Sonnenbergstrasse 74

CH-8603 Schwerzenbach, Switzerland



herewith declares that the product

hiermit deklarieren wir, dass das Produkt

Titrator/Analyzer

EasyPlus Easy Pro (EasyPlus xxxxx - Series)

For additional types, see page type code

For optional equipment, see page accessories

certified model:

--

Modell für Eichprüfung

is in conformity with the provisions of the following EC directives (incl. all applicable amendments)

mit den folgenden EG-Richtlinien (Inkl. Änderungen) übereinstimmt

2006/95/EC

Low voltage (LVD)

2004/108/EC

Electromagnetic compatibility (EMC)

and that the standards have been applied.

und die Normen zur Anwendung gelangen.

Last two digits of the year in which the CE marking was affixed: **13**

Die letzten zwei Zahlen des Jahres der Erst-CE-Kennzeichnung des Produkts mit dem CE Zeichen.

CH-8603 Schwerzenbach

11.01.2013

.....
Chris Radloff

General Manager GBA Analytical

.....
Christian Walter

Manager Strategic Business Unit Anachem

References of standards for this declaration of conformity, or parts thereof:

Harmonized standards of Europe and Switzerland:

Safety standards:

IEC/EN61010-1:2010

EMC standards (* Emission; ** Immunity):

IEC61326-1:2005 / EN61326-1:2006 (class B *)

IEC61326-1:2005 / EN61326-1:2006 (Industrial requirements **)

Metrological standards:

--

Environmental standards:

--, --, --

Standards for Canada, USA and Australia:

CAN/CSA-C22.2 No. 61010-1-12

UL Std. No. 61010-1 (3rd Edition)

FCC, Part 15, class A (Declaration) *

AS/NZS CISPR 11 *, AS/NZS 61000.4.3 **

EC - DECLARATION OF CONFORMITY

EG-Konformitätserklärung
KD-Nr.: 30065475

Doku-Nr.: 20120038

Type code

Typenschlüssel

other types of same construction:

andere Typen/Modelle mit der gleichen Konstruktion:



EasyPlus Easy Pro

Tested type

EasyPlus Easy pH, EasyPlus ET18

Titration General pH

EasyPlus Easy Cl, EasyPlus ET28

Titration General Cl

EasyPlus Easy Ox, EasyPlus ET38

Titration General Ox

EasyPlus Easy Pro, EasyPlus ET58

Titration General Pro

EasyPlus Easy KfV, EasyPlus ET08

Titration Karl Fischer

EasyPlus Easy Na

Analyser Sodium

EasyPlus ETy8 – Series is similar to EasyPlus Easy
xxx Series, except software/initial language.

Remarks

Bemerkungen:

EC - DECLARATION OF CONFORMITY

EG-Konformitätserklärung
KD-Nr.: 30065475

Doku-Nr.: 20120038

Accessories
Zubehör und Optionen



EasyPump
Sensors
Temperature sensor
External stirring unit

released: 5.06.2012 / CR

Seite 3 von 3

printed: Jan. 2013

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

www.mt.com/sodiumanalyzer

For more information

Mettler-Toledo AG, Analytical

CH-8603 Schwerzenbach, Switzerland

Tel. +41 (0)44 806 77 11

Fax +41 (0)44 806 73 50

www.mt.com

Subject to technical changes.

© Mettler-Toledo AG 10/2014

30079624B

