orbe^{co}
HELLIGE

SC450

Single Parameter Colorimeter Iron

For use with SC450-70
SC450-71
SC450-75

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Product name / Nom du fabricant / Nombre del productor:

SC450, Single Parameter Colorimeter

Declaration of EC-Conformity according to DIRECTIVE **2004/108/EG** OF THE
EUROPEAN PARLIAMENT AND OF THE COUNCIL of 2004, December the 15th.
The manufacturer declares that this product meets the requirements of the
following product family standard:

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
EN 61326-1:2006

Basic immunity test requirements (Table1) / Emission according to the
requirements for class B equipment

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Bradley K. Martell, President

Contents

• General notes	4
Guidelines for photometric measurements	4
Method notes	4
Replacement of batteries	5
• Functional description	6
Operation	6
Display backlight	7
Recall of stored data	7
Countdown	7
• Methods	8
Iron, with tablet	9
Iron with powder packs	10
Iron, (TPTZ) with powder packs	11
• Menu options	12
Menu selections	12
Recall of stored data	12
Transmitting stored data	12
Setting date and time	13
• Calibration Mode	13
User calibration	13
Factory calibration reset	15
• Technical data	16
Operating messages	17
Error codes	17

 **CAUTION** 

The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326.

Wireless devices, e.g. wireless phones, must not be used near the instrument.

General notes

Guidelines for photometric measurements


1. Vials, caps and stirring rods should be cleaned thoroughly **after each analysis** to prevent interference. Even minor reagent residues can cause errors in the test result.
2. The outside of the vial must be clean and dry before starting the analysis. Clean the outside of the vials with a towel to remove fingerprints or other marks.
3. Zero calibration and test must be carried out with the same vial as there may be slight differences in optical performance between vials.
4. The vials must be positioned in the sample chamber for zeroing and test with the Δ mark on the vial aligned with the ∇ mark on the instrument.
5. Always perform zeroing and test with the vial cap tightly closed. Only use the cap with a sealing ring.
6. Bubbles on the inside wall of the vial lead to incorrect measurements. To prevent this, remove the bubbles by swirling the vial before performing the test.
7. Avoid spillage of water into the sample chamber because this can lead to incorrect test results.
8. Contamination of the transparent cell chamber can result in incorrect readings. Check at regular intervals and – if necessary – clean the transparent cell chamber using a moist cloth or cotton swab.
9. Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the cell chamber or on the vial.
10. To avoid errors caused by stray light do not use the instrument in bright sunlight.
11. Always add the reagent tablets to the water sample straight from the foil without touching them with fingers.
12. The reagents must be added in the correct sequence.

Method notes


- Prior to measurement ensure that the sample is suitable for analysis (no major interferences) and does not require any preparation i.e. pH adjustment, filtration etc.
- Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children.
- Ensure proper disposal of reagent solutions.
- Material Safety Data Sheets are available on request
(Internet: www.orbeco.com)

General notes


Correct position of the vial (\varnothing 24 mm):

Correct filling of the vial:

Replacement of batteries:

CAUTION:

To ensure that the instrument is water proof:

- seal ring (E) must be in position
- battery compartment cover (B) must be fixed with the four screws

If the batteries are removed for more than one minute the date and time menu starts automatically when the photometer is switched on the next time.

Functional description

Operation

METHOD

METHOD

METHOD

0.0.0

METHOD

RESULT

Turn the unit on using the [ON/OFF] key.

The display shows the following:

Select the required test using the [MODE] key.

Scroll Memory (SM)

To avoid unnecessary scrolling for the required test method, the instrument memorizes the last method used before being turned off. When the instrument is turned on again, the scroll list comes up with the last used test method first.

The display shows the following:

Fill a clean vial with the water sample up to the 10 ml mark, screw the cap on and place the vial in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

The "Method" symbol flashes for approx. 8 seconds.

The display shows the following:

After zero calibration is completed, remove the vial from the sample chamber. The characteristic coloration appears after the addition of the reagents.

Replace the cap on the vial and place in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

(For Countdown/reaction period see page 7)

The "Method" symbol flashes for approx. 3 seconds.

The result appears in the display.

The result is saved automatically.

Repeating the test:

Press the [ZERO/TEST] key again.

Repeating the zero:

Press the [ZERO/TEST] key for 2 seconds.

Functional description

Display backlight

Press the [!] key to turn the display backlight on or off. The backlight is automatically turned off during the measurement.

Recall of stored data

If the instrument is turned on, press the [!] key for more than 4 seconds to access the recall menu.

Countdown / reaction period

If a reaction period is included in a method a countdown function can be used:

Zero
Test

Press the [!] key and hold.

Press the [ZERO/TEST] key.

Release the [!] key; the countdown starts.


After the countdown is finished the measurement starts automatically.

It is possible to interrupt the countdown by pressing the [ZERO/TEST] key. Measurement starts immediately.

Caution:

An incomplete reaction period can lead to incorrect test results.

Methods

Digestion procedure for the determination of total dissolved and undissolved iron:

1. Add 1 ml of concentrated sulfuric acid to 100 ml water sample. Heat and boil for 10 minutes or until all particles are dissolved. After cooling down the sample is set to a pH-value of 3 to 6 by using ammonia solution. Refill with deionised water to the previous volume of 100 ml and mix well. 10 ml of this pre-treated solution is used for the following analysis. Perform as described by the selected test method.
2. Water which has been treated with organic compounds like corrosion inhibitors must be oxidised where necessary to break down the iron. Therefore add 1 ml concentrated sulfuric acid and 1 ml concentrated nitric acid to 100 ml water sample and boil to approx. half volume. After cooling down proceed as described above.

Methods

Iron LR with Tablet 0.02 – 1.0 mg/l Fe

Fill a clean vial (24 mm Ø) with **10 ml of the water sample** and perform zero calibration (see "Operation").

Add **one IRON LR tablet** straight from the foil to the 10 ml water sample and crush the tablet using a clean stirring rod.

Close the vial tightly with the cap and swirl gently several times until the tablet is dissolved.

Place the vial in the sample chamber making sure that the \times marks are aligned.

Wait for a reaction period of 5 minutes.
(Countdown can be activated, see page 7)

The method symbol flashes for approx. 3 seconds.

The result is shown in the display as mg/l iron.

Tolerance: ± 0.05 mg/l Fe

Notes:

1. This method determines the total dissolved Iron as Fe^{2+} and Fe^{3+} .
2. The IRON (II) LR tablet is used for differentiation – as described above – instead of the IRON LR tablet.
3. For the determination of total dissolved and undissolved iron digestion is required.

Reagent re-order information

Iron Reagent Set	RS5370
Contains 1 each of the following:	
Iron LR Tablets, 100 ct.	RT537-OBT
Iron (II) LR Tablets, 100 ct.	RT542-0

For additional accessories, see page 17

Methods

FE 1
Tablet ▼
Powder

Iron with Powder Pack 0.02 – 3.0 mg/l Fe

0.0.0
Tablet ▼
Powder

Fill a clean vial (24 mm Ø) with **10 ml of the water sample** and perform zero calibration (see "Operation").

Add the contents of **one Ferro F10 Powder Pack** straight from the foil into the water sample.

Close the vial tightly with the cap and invert several times to mix the contents (Note 5).

Place the vial in the sample chamber making sure that the X marks are aligned.

! Zero Test
FE 1
Tablet ▼
RESULT Powder

Wait for a reaction period of 3 minutes (Note 2).

(Countdown can be activated, see page 7)

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Iron.

Tolerance: ± 0.05 mg/l Fe

Notes:

1. This method determines all dissolved and most undissolved forms of iron in the water sample.
2. Water samples containing visible rust should be allowed to react for at least five minutes.
3. For determination of total iron a digestion is required (see page 8).
4. Very strong alkaline or acidic water samples must be adjusted to a pH value between 3 and 5 before analysis.
5. Accuracy is not affected by undissolved powder.

Reagent re-order information

Iron Ferro Reagent Powder Pack, 100 ct.

RP056-0

For additional accessories, see page 17

Methods

Iron, (TPTZ) with Powder Pack 0.02 – 1.8 mg/l Fe

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

Fill a clean 24 mm vial with **10 ml deionised water** (this is the blank).

Fill a second clean 24 mm vial with **10 ml of the water sample** (this is the sample).


Add the contents of **one IRON TPTZ F10 Powder Pack** straight from the foil into each vial.

Close the vials tightly with the caps and swirl several times to mix the contents (30 sec.).

Wait for a reaction period of 3 minutes.

(Countdown can be activated, see page 7)

After the reaction period is finished proceed as follows:

Place the vial (the blank) in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 8 seconds.

Remove the vial from the sample chamber.

Place the vial (the sample) in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Iron.

Tolerance: ± 0.05 mg/l Fe

Notes:

1. Iron in the sample, including precipitated or suspended iron (e.g. rust) is reduced to Fe²⁺. The TPTZ reagent reacts with Fe²⁺ to a colored complex.
2. For determination of total iron a digestion is required (see page 8).
3. Rinse all glassware with 1:1 Hydrochloric acid solution first and then rinse with deionised water to remove iron deposits that can cause slightly high results.
4. Strong alkaline or acidic water samples must be adjusted between pH 3 and pH 8 before the reagent is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).

Reagent re-order information

Iron Ferro Reagent Powder Pack, 100 ct.

RP055-0

For additional accessories, see page 17

Menu options

Menu selections

Press the [MODE] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Allow the 3 decimal points to be displayed before releasing the [MODE] key.

The [!] key allows for selection of the following menu points:

- ▲ diS recall stored data
- ▲ Prt printing stored data
- ▲ setting the date and time
- Cal user calibration

▲ diS – Recall of stored data

After confirming the selection with the [MODE] key the photometer shows the last 16 data sets in the following format (automatically proceeds every 3 seconds until result is displayed):

Number n xx (xx: 16...1)
Year YYYY (e.g. 2010)
Date mm.dd (monthmonth:dayday)
Time hh:mm (hourhour:minuteminute)
Test Method
Result x,xx

The [ZERO/TEST] key repeats the current data set.

The [MODE] key scrolls through all stored data sets.


Quit the menu by pressing [!] key.

▲ Prt – Transmitting stored data (to Printer or PC)

Note: To print data, or to transmit to a PC, the optional IRIM (Infrared Interface Module) is required.

The IRIM Module and the connected printer/PC must be ready. Press the [MODE] key to start the transmitting, the instrument displays "PrtG" (Printing) for approx. 1 second followed by the number of the first data set and its transmission. All data sets will be transmitted one after the other. After finishing the instrument switches to test mode.

The print job can be cancelled by pressing the [On/Off] key. The instrument will turn off.

Menu options – Calibration Mode

E 132

If the instrument is not able to communicate with the IRIM, a timeout occurs after approx. 2 minutes. The error E 132 is displayed for approx. 4 seconds. Subsequently, the instrument returns to test mode (see also IRIM manual).

2 Setting date and time (24-hour-format)

After confirming the selection with the [MODE] key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. The same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the [!] key to continue setting the minutes in steps of 1.

Increase the value by pressing the [MODE] key.

Decrease the value by pressing [ZERO/TEST] key.

Proceed to the next value to be edited by pressing [!] key.

After setting the minutes and pressing the [!] key the display will show "IS SET" and the instrument returns to the measurement mode.

Cal User calibration

Note:

user calibration (Display in calibration mode)

factory calibration (Display in calibration mode)

After confirming the selection with the [MODE] key the instrument will show CAL/"Method".

Scroll through methods using the [MODE] key.

Fill a clean vial with the standard up to the 10 ml mark, screw the cap on and place the vial in the sample chamber making sure that the X marks are aligned.

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 8 seconds.

The display shows the following in alternating mode:

Perform calibration with a standard of known concentration (see "Operation").

Press the [ZERO/TEST] key.

The method symbol flashes for approx. 3 seconds.

Calibration Mode

RESULT

The result is shown in the display, alternating with CAL.

CAL

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the [ON/OFF] key.

 **Mode**

Changing the displayed value:

 **Zero**
 **Test**

Pressing the [MODE] key once increases the displayed value by 1 digit.

CAL

Pressing the [ZERO/TEST] key once decreases the displayed value by 1 digit.

RESULT + x

 **On**
 **Off**

Press the corresponding key until the reading equals the value of the calibration standard.

: **:**


By pressing the [ON/OFF] key, the new correction factor is calculated and stored in the user calibration software.

Confirmation of calibration (3 seconds).

Calibration Mode

Factory calibration reset

Resetting the user calibration to the original factory calibration will reset all methods and ranges.


A user calibrated method is indicated by a "Cal" symbol while the test result is displayed.

To reset the calibration press both the [MODE] and [ZERO/TEST] key and **hold**.


Switch the unit on using the [ON/OFF] key.
Release the [MODE] and [ZERO/TEST] keys after approx. 1 second.

The following messages will appear in turn on the display:


The factory setting is active.
(SEL stands for Select)

or:

Calibration has been set by the user.
(If the user calibration is to be retained, switch the unit off using the [ON/OFF] key).

Calibration is reset to the factory setting by pressing the [MODE] key.

The following messages will appear in turn on the display:

Turn the unit off using the [ON/OFF] key.

Technical Data

Technical Data

Instrument	triple wavelength, automatic wavelength selection, direct reading colorimeter
Light source:	LEDs, interference filters (IF) and photosensor in transparent cell chamber. Wavelength specifications of the IF: 530 nm $\Delta \lambda = 5$ nm 560 nm $\Delta \lambda = 5$ nm 580 nm $\Delta \lambda = 5$ nm
Wavelength accuracy	± 1 nm
Photometric accuracy*	3% FS (T = 20° C – 25° C)
Photometric resolution	0.01 A
Power supply	4 batteries (AAA/LR 03) lifetime: approx. 17 hours or 5000 tests
Auto-OFF	automatic switch off 10 minutes after last keypress
Display	backlit LCD (on keypress)
Storage	internal ring memory for 16 data sets
Interface	IR interface for data transfer
Time	real time clock und date
Calibration	user and factory calibration resetting to factory calibration possible
Dimensions	155 x 75 x 35 mm (LxWxH)
Weight	approx. 260 g (incl. batteries)
Ambient conditions	temperature: 5–40°C rel. humidity: 30–90 % (non-condensing)
Waterproof	floating; IP 68 (1 hour at 0.1 meter)
CE	Certificate for Declaration of CE-Conformity at www.orbeco.com

*measured with standard solutions

To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

Operating messages – Error codes

Operating messages

Hi	Measuring range exceeded or excessive turbidity.
Lo	Result below the lowest limit of the measuring range.
	Replace batteries, no further tests possible.
btLo	Battery capacity is too low for the display backlight; measurement is still possible.
Store Print Cal	Date Time
RESULT	

Error codes

E27 / E28 / E29	Light absorption too great. Reasons: e.g. dirty optics.
E 10 / E 11	Calibration factor "out of range"
E 20 / E 21	Too much light reaching the detector.
E23 / E24 / E25	Too much light reaching the detector.
E 22	Battery capacity was too low during measurement. Change battery.
E 70	FE (Tablet): Factory calibration incorrect / erased
E 71	FE (Tablet): User calibration incorrect / erased
E 72	FE 1 (Powder Pack): Factory calibration incorrect / erased
E 73	FE 1 (Powder Pack): User calibration incorrect / erased
E 74	FE 2 (Powder Pack): Factory calibration incorrect / erased
E 75	FE 2 (Powder Pack): User calibration incorrect / erased

Instrument accessory information

IRiM	MC500-60
24 mm Sample Cells, 12 pack	L197620
Light shielding rings	L197626
Cleaning cloth for sample cells	L197635
Sample cell rack	9999

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Orbeco-Hellige, Inc.