07519-05 Pump Head with 07519-70 Large Cartridge
Tête de pompe 07519-05 à cartouche épaisse 07519-70
Pumpenkopf 07519-05 mit großer Schlauchkassette 07519-70
Cabeza de bomba 07519-05 con cartucho grande 07519-70
Testa pompante 07519-05 con cartuccia grande 07519-70

OPERATING MANUAL

CARTRIDGE PUMP HEAD SYSTEM

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07519-70

**PUMP FOR LIQUIDS
ORIGINAL INSTRUCTIONS**

CONTENTS

Title	Page
SAFETY PRECAUTIONS	2
INTRODUCTION	3
(1) Application Data.	3
(2) General Description	3
INSTALLATION	4
SETUP	5
(1) Select Cartridges	5
(2) Load Cartridges	5
(3) Cartridge Removal	7
(4) Bi-Directional Pumping	7
(5) Partial Bank Pumping	7
OPERATION	8
(1) Select Tubing and Pump Heads	8
(2) Select Pump Speed	9
(3) Adjust Occlusion Settings on Cartridges	9
MAINTENANCE	10
Service Parts	10
SPECIFICATIONS	10
WARRANTY	11
PRODUCT RETURN	11
TECHNICAL ASSISTANCE	11
APPENDIX A MASTERFLEX® L/S® Pump Drive Types	12
APPENDIX B Cartridge Loading Capacity	12

C-FLEX, PHARMED, NORPRENE, TYGON—Reg TM Saint-Gobain Performance Plastics Corp.
SANTOPRENE — Reg TM AES.
VITON — Reg TM E.I. duPont DeNemours & Co.
Trademarks bearing the © symbol in this publication are registered in the U.S. and in other countries.

SAFETY PRECAUTIONS

WARNINGS: *Tubing breakage may result in fluid being sprayed from pump. Use appropriate measures to protect operator and equipment.*

Turn off drive before removing or installing Cartridges. Safety guards are provided to minimize risk of fingers getting caught between the roller mechanism and the base of the module. However, be safe—Keep your fingers away from these areas.

CAUTIONS: *Excessive occlusion can cause high pump temperatures and early tubing and roller failure. Do not set occlusion tighter than what is recommended in Table 4.*

Use only tubing/cartridge combinations defined by Tables 1, 2, and 3. Use of other combinations could cause malfunction or damage to the pump.

Not using the BOOT will damage the pump.

Explanation of Symbols

CAUTION: *Risk of Danger. Consult Operator's manual for nature of hazard and corrective actions.*

CAUTION: *Risk of crushing. Keep fingers away from rotor while pump is in operation. Stop pump before loading or unloading tubing.*

WARNING: PRODUCT USE LIMITATION

These products are not designed for, nor intended for use in, patient-connected applications, including, but not limited to, medical and dental use and, accordingly, have not been submitted for FDA approval.

INTRODUCTION

The instructions in this manual are task-oriented for easy reference. You can go directly to a particular section and quickly find the answers. Appendix A and B list MASTERFLEX® L/S® pump drive types which can be used with this system.

Pump Heads

MODEL	TYPE
Model 07519-05	3-Roller, 4 or 8 Channel
Model 07519-06	4-Roller, 4 or 8 Channel

Cartridges

MODEL	TYPE
Model 07519-70	Large
Model 07519-80	Small

These Cartridge Pumps are designed to provide up to 8 simultaneously driven pump channels and incorporate the following features.

- Simultaneous dispensing into 8 separate containers with one set-up and pump cycle.
- The four-roller design provides reduced pulsation.

(1) Application Data

The Cartridge Pump System is designed for multi-channel use where close occlusion control of a variety of different flows is required. Tubing found in Table 3 may be used to deliver from .074 to 1700 mL/min.

(2) General Description

Each Pump Head System is composed of a heavy duty Pump Head and two sizes of snap-in tubing Cartridges, each with individually adjustable occlusion settings. The occlusion “wedge” design provides accurate and more balanced flow control.

The Pump Heads contain either three or four rollers. Both models can operate with a variety of tubing and with either the small or the large Cartridges.

The Pump Head quickly mounts to most existing MASTERFLEX® L/S® Pump Drive types (see Appendix A and B).

INSTALLATION

Tool Required: Allen Key (provided).

Mount the Pump Head to a MASTERFLEX® L/S® Pump Drive as follows:

1. Install the plastic Boot over the tang at the end of pump shaft, (Figure 1). This Boot prevents metal-to-metal contact, and thus reduces noise. This Boot should be inspected and replaced when worn or at 500 hour intervals.

CAUTION: *Not using the boot will damage the pump.*

NOTE: The boot with closed ends (A-3480) will fit most drives. Use of the boot with open ends (A-3514) may be required on some drives. Both boots are supplied with the pump.

2. Connect the Pump to Drive by aligning the tang on Pump Head, (including installed boot) with the slot in the motor Drive shaft.
3. Attach Pump Head to drive unit with four screws (provided) using the Allen Key (also provided).

WARNING: *Tubing breakage may result in fluid being sprayed from pump. Use appropriate measures to protect operator and equipment.*

4. If the Pump Head requires support, install the appropriate length support screw with plastic foot and locknut (three provided). Level the Pump Head with the drive unit by adjusting the foot. Tighten the locknut against the bottom of the pump body.




FIGURE 1 PUMP MOUNTING

A-3480

A-3514

FIGURE 2 PLASTIC BOOTS

SETUP

Two Cartridge sizes (0.72 in and 0.36 in thick) (see Figure 3) accommodate a broad range of tubing sizes as noted in Tables 1 and 2.


FIGURE 3 SMALL AND LARGE CARTRIDGES

(1) Select Cartridges

NOTE: These Pump Heads are designed to work with the Model 07519-70 and 07519-80 Cartridges. Do not use Model 07519-75 or 07519-85 Cartridges, which are designed for different Pump Heads.

Depending on which Pump Head is used, up to four large or eight small tubing Cartridges can be teamed up with the Pump Head. But there is a limit depending on available drive power and pump load. For full tubing flow rates at 0 psi, the table in Appendix B shows cartridge capacities according to tubing size and drive Type (see Appendix A for description of Drive Types). Cartridge sizes can be intermixed as long as the total load does not exceed the drive capability.

The small Cartridge is designed to operate with MASTERFLEX® precision tubing sizes L/S® 13 and L/S® 14, and will also accommodate specially designed Microbore Tubing Sets (.9 mm wall — see Tables 1 and 3). The small Cartridge can be used with either the three- or the four-roller Pump Head.

The large Cartridge is designed to operate with MASTERFLEX® L/S® precision tubing sizes L/S® 14, L/S® 16, L/S® 25 and L/S® 17, with either the three- or four-roller Pump Head (see Tables 2 and 3).

Use only MASTERFLEX Precision Tubing and Microbore Tubing Sets with MASTERFLEX Pumps to ensure optimum performance. Use of other tubing may void applicable warranties.

(2) Load Cartridges

Easy, snap-in mounting and quick-release removal come with a Cartridge system—no tools required. Since the Cartridges function uniformly in either direction, Cartridge orientation is not important. However, you may wish to alternate Cartridges so that the Adjustment Knob on each is easier to turn when several small Cartridges are mounted close together.

WARNING: Turn off drive before removing or installing Cartridges. Safety guards are provided to minimize risk of fingers getting caught between the roller mechanism and the base of the module. However, be safe—Keep your fingers away from these areas.

CAUTION: Use only tubing/cartridge combinations defined by Tables 1, 2, and 3. Use of other combinations could cause malfunction or damage to the pump.

With the drive stopped, follow this procedure to install Cartridges.

1. Turn the Occlusion Adjustment Knob (see Figure 4) counterclockwise to maximum open. This step is not necessary when the occlusion is set at the factory setting or if the occlusion has already been set appropriately during previous pumping.
2. Set the Tubing Retainer according to tubing size (see Figure 6A) when using the large Cartridge. No retainer adjustment is required when using the small Cartridge. (**Note:** The retainer button located on the Occlusion Adjustment Knob-side of the large Cartridge is a loose component and may fall out. To reassemble, push it back into place while pressing in the button.)


FIGURE 4 CARTRIDGE ASSEMBLY




FIGURE 5 MICROBORE TUBING SET


3. Select tubing and place in Cartridge (see Figure 6B). If using a Microbore Tubing Set (see Figure 5), place the Collar up against the Tubing Retainer (see Figure 6C) on both sides of the Cartridge.
4. Attach the Hook Leg of the Cartridge onto one of the Support Rods while holding the tubing between the Cartridge legs and against the Tubing Retainer (see Figure 6D).
5. Swing the Cartridge with the tubing over the rotor and push down on the Tab (see Figure 4) until the Latch snaps closed. (See Figure 6E.) **Note:** Depress the Cartridge as required in the location indicated in Figure 6E to facilitate latching and unlatching.
6. Lightly pull the tubing at the outlet and inlet to remove the slack when using continuous tubing. No slack removal is required when using a Microbore Tubing Set.
7. Adjust the occlusion setting (see OPERATION section). For a nominal setting with MASTERFLEX tubing, turn adjustment screw until inside edge of the wedge is aligned midway between #3 and #4 on the label.


6A


6B

6C

6D

6E

FIGURE 6 TUBING AND CARTRIDGE LOADING

(3) Cartridge Removal

Follow this procedure to unload the Cartridge.

1. Turn the Drive OFF.
2. Push down on latch side of Cartridge (see Figure 6E) and pull Latch Leg outward to unhook from rod. The Latch Release (see Figure 4) can be depressed to aid in this process as desired.
3. Pivot Cartridge upwards and unhook Hook Leg from rod and remove Cartridge.

(4) Bi-Directional Pumping

Fluids can be pumped into and back out of containers by reversing the drive. On drives that cannot be reversed, use longer tubing lengths and simply reverse the loaded Cartridge.

WARNING: *Tubing breakage may result in fluid being sprayed from pump. Use appropriate measures to protect operator and equipment.*

(5) Partial Bank Pumping

The Pump can be operated with either a partial or a full bank of Cartridges. There is an elastic, VITON® fluoroelastomer O-ring on each of the two metal rods of the pump (see Figure 7) that holds the Cartridges in place on the pump. These O-rings should be adjusted as follows. If the pump is partially loaded with Cartridges, slide the O-rings along the rods, up against the front of the outermost Cartridge, to keep it in place (note position A, for example, in Figure 7). If the pump is fully loaded with Cartridges, slide the O-rings to the front of the pump as far as possible (note position B in Figure 7).


FIGURE 7 O-RING POSITION

OPERATION

This section describes the procedures for obtaining the desired performance.

(1) Select Tubing and Pump Heads

Use Tables 1, 2, and 3 to select the tubing, Cartridge, and Pump Head to provide the required single channel flow rate.

CAUTION: Use only tubing/cartridge combinations defined by Tables 1, 2, and 3. Use of other combinations could cause malfunction or damage to the pump.

TABLE 1 SMALL CARTRIDGE FLOW RATES—07519-80

Pump Head	Drive rpm	Units	MICROBORE TUBING (mm ID)				MASTERFLEX® L/S® TUBING	
			0.89	1.42	2.06	2.79	L/S® 13	L/S® 14
			FLOW RATES					
3-Roller	6 600	mL/min	.44	1.0	2.2	3.8	.30	1.3
		mL/min	44	100	220	380	30	130
3-Roller	1 100	mL/min	.074	.17	.37	.63	.050	.21
		mL/min	7.4	17	37	63	5.0	21
4-Roller	6 600	mL/min	.44	1.0	2.2	3.8	.28	1.2
		mL/min	44	100	220	380	28	120
4-Roller	1 100	mL/min	.074	.17	.37	.63	.047	.20
		mL/min	7.4	17	37	63	4.7	20

TABLE 2 LARGE CARTRIDGE FLOW RATES—07519-70

Pump Head	Drive rpm	Units	MASTERFLEX® L/S® TUBING			
			L/S® 14	L/S® 16	L/S® 25	L/S® 17
			FLOW RATES			
3-Roller	6 600	mL/min	1.3	4.6	10	17
		mL/min	130	460	1000	1700
3-Roller	1 100	mL/min	.21	.76	1.7	2.8
		mL/min	21	76	170	280
4-Roller	6 600	mL/min	1.2	3.8	8.3	14
		mL/min	120	380	830	1400
4-Roller	1 100	mL/min	.20	.63	1.4	2.3
		mL/min	20	63	140	230

NOTE: The Pumps operate at the same speed as the drive.

NOTE: Improper tubing stretch could result in reduced flow rate (see LOAD CARTRIDGE section).

TABLE 3 COMPATIBLE TUBING

Tubing Material	*MICROBORE TUBING SET (mm ID)	MASTERFLEX® L/S® TUBING	
	.89, 1.42, 2.06, 2.79	L/S® 13, 14, 16, 25	L/S® 17
Soft			
Silicone (platinum)	X	X	X
Silicone (peroxide)		X	X
TYGON® Silicone		X	X
C-FLEX® (50 A)		X	X
Stiff			
PHARMED® BPT (65)		X	
TYGON LFL		X	
NORPRENE® food (A 60 F)		X	
TYGON lab (R-3603)		X	
NORPRENE (A 60 G)		X	
TYGON fuel (F-4040-A)		X	
PVC	X		
VITON	X		
SANTOPRENE®	X		

*Use only specially designed MICROBORE TUBING SETS with Cartridges (Figure 5).

(2) Select Pump Speed

Select fixed speed drive or adjust speed of variable speed drive to provide desired nominal flow rate within the rpm speed range shown for the pump model selected.

(3) Adjust Occlusion Settings on Cartridges

CAUTION: *Excessive occlusion can cause high pump temperatures and early tubing and roller failure. Do not set occlusion tighter than what is recommended in Table 4.*

NOTE: With large and small Cartridges, the scale provides nominal occlusion for MASTERFLEX® L/S® precision tubing at the #3–#4 setting.

With a multi-channel Cartridge system, flow rate, pressure sensitivity and tubing life can be fine-tuned—even while the pump is running. No lost time, and with operating speeds up to 600 rpm, you get higher maximum flows with the same size tubing.

Figure 8 (large and small Cartridge) shows the wedges adjusted at #3 on the occlusion setpoint scale for the MASTERFLEX tubing.

The Occlusion Adjustment Knob (see Figure 8) controls a screw mechanism to move the occlusion wedges. Once a Cartridge is adjusted to a desired setting for a particular application, there is no need to reset. Or, you can record settings and quickly reproduce them at a later time.

(a) Normal Occlusion Settings

The following table shows the appropriate occlusion settings to provide satisfactory performance for various tubing sizes. (For optimum occlusion settings, see following section.)

TABLE 4 RECOMMENDED OCCLUSION SETTINGS

Tubing	Large Cartridge Scale	Small Cartridge Scale
Microbore Tubing Set	N/A	#1 – #2
MASTERFLEX® L/S® Precision	#3 – #4	#3 – #4

(b) Occlusion Setting Procedure:

1. Select the recommended occlusion value from the table.
2. Turn the Occlusion Adjustment Knob to align the inside edge of the white Wedge with the scale number. (Clockwise rotation increases the occlusion.)

(c) Optimized Occlusion Settings

Some applications require additional fine-tuning of the occlusion setting to vary the flow rate for a particular tubing, to reduce flow variations caused by changes in system pressure, or to increase tubing life.

1. Adjust the Occlusion Wedges as described in the preceding table.
2. Refine this setting, depending on your objectives:

NOTE: To Maximize Tubing Life

While running the pump, simply reduce the adjusted occlusion setting by turning the Occlusion Adjustment Knob counterclockwise to move the wedges to a higher setpoint scale number (toward 5). Continue turning the knob counterclockwise until the flow drop-off or the pressure sensitivity becomes unacceptable, then turn knob slightly clockwise.

NOTE: To Reduce or Eliminate Flow Drop-Off Rate with Pressure Fluctuations

While running the pump alternate the backpressure between the expected high and low values and increase the adjusted occlusion setting. Turn the Occlusion Adjustment Knob clockwise to move the wedges to a lower setpoint scale number (toward 1) until the flow drop-off is minimized.

NOTE: To Fine-Adjust the Flow Rate

Increase or decrease the occlusion value to vary the flow.


FIGURE 8 OCCLUSION INDICATIONS

MAINTENANCE

No maintenance required. Wipe Pump with a clean cloth and a mild detergent. Never immerse or use excessive fluid.

Service Parts:

Boot (closed ends, see Figure 2)	A-3480
Boot (open ends, see Figure 2)	A-3514
Allen Key	A-4376
Screws 8-32 x 1 Socket Head	B-1079-0526
Level Feet Assembly	A-4128-0002, (2 in)
Level Feet Assembly	A-4128-0003, (1-1/4 in)
Level Feet Assembly	A-4128-0004, (1-5/8 in)

SPECIFICATIONS

PUMP HEAD

Speed Range:	0 to 600 rpm
Torque Load:	5 oz-in (max.) with no Cartridges
Operating Temperature Range:	0° to 40°C (32°F to 104°F)
Storage Temperature Range:	-40°C to 110°C (-40°F to 230°F)
Noise level:	<70 dBA @ 1 meter
Construction and Chemical Resistance:	Fiberglass-reinforced polysulfone housing; stainless steel, anodized aluminum, sealed stainless steel ball bearings and Buna N safety guard.
Compliance:	EN809 (EU Machine Directive)

CARTRIDGES

Construction:	Durable polycarbonate frame with filled nylon and acetal components. Anodized aluminum adjustment knob on stainless steel screw, VITON O-ring.
Dimensions:	
Large	4 in Tall x 5 in W x 0.72 in Thick (102 mm x 127 mm x 18.3 mm)
Small	4 in Tall x 5 in W x 0.36 in Thick (102 mm x 127 mm x 9.1 mm)

WARRANTY

Use only MASTERFLEX Precision Tubing and Microbore Tubing Sets with MASTERFLEX Pumps to ensure optimum performance. Use of other tubing may void applicable warranties.

The Manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the warranty period, the problem will be corrected at no charge if it is not due to misuse or abuse on your part as determined by the Manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced to you.

The warranty period for this product is two (2) years from the date of purchase.

PRODUCT RETURN

To limit charges and delays, contact the seller or Manufacturer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. When returning the product, please state the reason for the return. For your protection, pack the product carefully and insure it against possible damage loss. Any damages resulting from improper packaging are your responsibility.

TECHNICAL ASSISTANCE

If you have any questions about the use of this product, contact the Manufacturer or authorized seller.

US & Canada only

Toll Free 1-800-MASTERFLEX | 1-800-637-3739

Outside US & Canada

1-847-549-7600 | 1-847-381-7050

***EN809 manufactured by:**

Cole-Parmer Instrument Company

28W092 Commercial Avenue, Barrington, IL 60010

techinfo@masterflex.com | www.masterflex.com

APPENDIX A
MASTERFLEX® L/S® PUMP DRIVE TYPES
FOR USE WITH CARTRIDGE PUMP SYSTEMS

Max rpm	Power (hp)	Torque oz-in	Drive Type
600	1/10	180	I
120	1/20	180	
600	1/20	90	II
100	1/10	360	III

Note: Check Max rpm and Power of the Drive or use the Torque rating to determine the Drive Type.

APPENDIX B
CARTRIDGE LOADING CAPACITY

Cartridge Size	Tubing Size	Tubing Material*	Torque** (Avg. oz-in)	Max. No. of Cartridges		
				Drive Type I	Drive Type II	Drive Type III
Small	Microbore Tubing Set .89 mm dia.	Soft	‡	8	8	8
		Stiff	‡	7	3	8
	Microbore Tubing Set 1.42 mm dia.	Soft	‡	8	8	8
		Stiff	‡	—	—	—
	Microbore Tubing Set 2.06 mm dia.	Soft	‡	8	8	8
		Stiff	‡	7	3	8
	Microbore Tubing Set 2.79 mm dia.	Soft	‡	8	8	8
		Stiff	‡	7	3	8
	L/S® 13	Soft	10	8	8	8
		Stiff	26	7	3	8
	L/S® 14	Soft	10	8	8	8
		Stiff	31	6	3	8
Large	L/S® 14	Soft	10	4	4	4
		Stiff	31	4	3	4
	L/S® 16	Soft	14	4	4	4
		Stiff	31	4	3	4
	L/S® 25	Soft	43	4	2	4
		Stiff	80	2	1	4
	L/S® 17	Soft	43	4	2	4

*See Table 3 for a listing of soft and stiff tubings.

**Torque requirements are the same for both 3- and 4-roller rotor pump heads.

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