

Outstanding Weighing Performance Even under Harsh Conditions



Reducing the impact of air drafts

This white paper describes the influences of air flow during the weighing process. Recent innovations in the design of the weighing pan for precision balances and scales help to minimize any errors due to air drafts. The innovative new weighing pan makes it possible to deliver results up to twice as fast, even under harsh conditions. What's more, the results show an outstanding 86% improvement in repeatability (precision) for a 1 mg readability balance under harsh conditions (inside a safety cabinet).

Content

1	Introduction	2
2	Considerations when Weighing with Precision Balances	2
3	Experimental Procedure	5
4	Summary of Results	10
5	Conclusion	16
6	Tips and Recommendations for Reliable Weighing	16
7	Further Information on External Influences	19
8	About the Author	19

1. Introduction

Weighing is one of the most common activities carried out in the laboratory. An electronic balance does not directly measure a mass, rather the force of its weight. The uncertainty of this force depends not only on the technical specification of the measuring instrument but also on many environmental effects such as temperature variations and air drafts.

The environmental conditions have an important impact on precision weighing. Air drafts from open doors or windows, air conditioning units or fans, as well as air escape in a fume hood or safety cabinet when working with malodorous, toxic or aggressive substances, and even movement of people around the balance, can influence the weighing performance. Drafts can cause slower settling times and drift of measurements. In these situations, the operators of a scale or balance cannot rely on accurate and reliable measurements. For more than ten years, state-of-the-art analytical balances have been equipped with glass draft shields using motorized opening for precise weighing, even in unstable environments. The user has access to the weighing pan via a partially opening glass draft shield. The balance stabilizes faster because there is less disturbance from air currents than in an unshielded environment. Additionally, there is the option to use cleverly designed hanging weighing pans, SmartGrid™, with a grid structure. These pans present less surface area for air turbulence in the weighing chamber to act upon, than conventional weighing pans. The result is shorter settling times and faster availability of the results. But what measures can be taken to optimize weighing on precision balances?

2. Considerations When Weighing with Precision Balances

Precision balances are mainly used in open environments of laboratories and production sites where the user has only limited possibilities to shield his equipment against air drafts.

In general, modern precision balances provide a high weighing performance. This means fast and reliable weighing results, coupled with a level of security that meets the demands of regulatory bodies. Air flow from different sources like people moving, opening of doors, heating radiators, and fans from computer and laboratory equipment influences this performance. The most critical metrics to determine the effect of these external influences on a balance are the "average settling time" and the "average repeatability". The inherent risk in this situation is that not enough care will be taken to prevent disturbing influences.

The air currents can be minimized using weighing pans which are specifically designed. An open grid structure may seem, on the face of it, to address the issue of air drafts on precision balances. However, there are inherent drawbacks to this approach. One issue, in particular, is the difficulty to clean and maintain a weighing pan with an open grid structure. The latest innovation in design of weighing pan is a sturdy quadruped platform. This simple, yet proven geometry, was developed based on consideration of four primary design criteria:

- 1) Performance (speed and precision)
- 2) Placement of tare container
- 3) Weighing pan material and construction
- 4) Ease of cleaning

2.1 Performance

Air currents act as a force on the weighing pan and will influence the weighing results. Minimizing the effect of this perturbing force with a clever geometry of the pan improves the above mentioned performance (Fig. 1). The pan is primarily designed to minimize the negative effects of air currents without the need for any draft shield. The drip-tray beneath also includes a draft ring which provides an additional mechanical protection against air flow influence on the weighing pan. Experiments under different environmental conditions demonstrate the positive effect

on the weighing performance (see Chapter 3). Speed and precision are key performance metrics of interest to the balance user. The parameters used to measure these elements are settling time of the balance (speed) and repeatability of the balance measurement (precision).




Figure 1: Forces caused by air currents on standard pan versus SmartPan™. The SmartPan™ displays less air resistance

2.2 Placement of Tare Container

The SmartPan™ is compatible with any kind of tare container without any danger of spillage (Fig. 2). As precision balances are used in a wide range of applications in different industry segments and workplaces, various sizes of tare container, such as vials, flasks, beakers, petri capsules, cylinders, boxes or buckets could be used every day to contain and weigh samples.




Figure 2: The SmartPan™ is compatible with many different types of tare container

Whether small, large or bulky tare containers need to be used, the SmartPan™ geometry offers a considerable flat support surface where any kind of tare container or sample can be easily placed in a stable position without risk (Fig. 2). To prevent or minimize eccentricity errors, the cross-shaped design of this weighing pan provides the user with an easy and intuitive guide to correctly place the tare container or sample right in the center of the weighing pan.

2.3 Weighing Pan Material and Construction

Weighing pan material and construction were carefully considered in the design of the innovative SmartPan™.

A zinc alloy was selected as the ideal material for the new weighing pan design, because of its better machinability and superior tensile and impact strength, compared to the conventional aluminium weighing pan. This leads to a more rugged construction.

Finite Element (FE) Analysis is a computer simulation technique that is often used in mechanical design processes to predict and optimize the behaviour of complex objects. It allows a highly detailed analysis of any product or equipment design, in order to carry out many physical tests, such as stress, vibration or heat transfer analysis, for example.

In this case FE was employed to determine the force distribution, and therefore strength and robustness, of the new weighing pan design.

Figures 3 and 4 show a finite element simulation with the static force distribution of the conventional pan compared to the new weighing pan design. The finite element method calculates component displacements, strains, and stresses under internal and external loads. Thus it is ensured that the geometry remains in the linear elastic range.

The simulations were executed with a 5 kg weight put in the middle of the balance XPE6003SD5 (6100 g maximum capacity and 5 mg readability).




Figure 3: Finite element simulation of the force distribution of the conventional aluminium pan geometry




Figure 4: Finite element simulation of the force distribution of the new zinc weighing pan geometry

The observed maximum stress on the SmartPan™ is approximately 7.45 Mpa (Megapascal) which is a non-critical and negligible value. This value is far below the critical proof stress value of 300 Mpa which means irreversible deflection of the material. One can easily observe, from the colour maps in Figs. 3 & 4, that the quadruped pan has the same force distribution as the standard pan; therefore the SmartPan™ is very robust and the factor of safety is retained.

2.4 Ease of Cleaning

This innovative new weighing pan with quadruped design (SmartPan™) has many advantageous features for improving the weighing performance and the ease of cleaning. In order to allow an excellent cleaning performance the pan is made of die cast zinc. The choice of material offers a high resistance to acid and aggressive cleaning agents. The design incorporates a drip-tray underneath, similar to that on the analytical balances, to collect any substances which are accidentally spilled. The rounded edges and the gap-free design of all components enable convenient cleaning operations.

3 Experimental Procedure

The experiments described in Table 1 were carried out to demonstrate the performance of the new geometry of weighing pan compared with the conventional weighing pan design (standard pan). The key performance metrics measured were settling time and repeatability of the balance.

Each test was carried out was repeated with (a) the balance fitted with a standard weighing pan and (b) the same balance fitted with a SmartPan™. For every test carried out, three replicates were each weighed 20 times and the reported results are the average values.

Tests were made in different environmental situations (standard – on open lab bench; harsh – under air conditioning unit or in safety cabinet), in order to represent specific environmental situations that may be realistically encountered in a laboratory environment, as well as with different balance models (5 mg readability and 1 mg readability).

The balance settings used during these experiments were:

- Weighing mode: universal
- Environment: standard
- Value release: reliable and fast

Experiment	1	2	3	4
Conditions	Standard	Harsh	Standard	Harsh
Description	On open lab bench	Under air conditioning unit	On open lab bench	In safety cabinet
Weighing pan	(a) Standard pan (b) SmartPan™	(a) Standard pan (b) SmartPan™	(a) Standard pan * (b) SmartPan™	(a) Standard pan * (b) SmartPan™ *
Balance model	XPE6003SD5	XPE6003SD5	XPE1203S	XPE1203S
Capacity of balance	6100 g	6100 g	1210 g	1210 g
Readability of balance	5 mg	5 mg	1 mg	1 mg
Test weights	400 g 2 kg 5 kg	400 g 2 kg 5 kg	200 g 500 g	200 g 500 g
Number of weighings	60	60	60	60

* equipped with draft shield

Table 1: Description of Experiments comparing Standard weighing pan with new SmartPan™ design

3.1 Experiment 1: Weighing on a 5 mg readability balance under standard conditions

Experiments 1(a) and 1(b) were carried out using a 5 mg readability precision balance (XPE6003SD5) on an open lab bench. The experimental set-up is illustrated in Fig. 5. All necessary factors were taken into account to ensure that the balance was working in the best possible environment, such as positioning the balance: on a stable lab bench; in a laboratory temperature of 20 – 23 °C; in a relative humidity of 35 – 45%; in vibration-free surroundings (i.e. not located near any vibration sources such as fans, vacuum pumps or generators); and in a suitable location (i.e. not too close to doors and with minimal passing foot-traffic).

This weighing environment was defined as being under "Standard conditions".

For recommendations on the ideal environment for the best weighing results, see Chapter 6.

In Experiment 1(a) the balance was equipped with a standard weighing pan. Test weights of 5 kg, 2 g, and 400 g were placed repeatedly on the weighing pan of the balance under test. Three replicates were each weighed 20 times, with a total of 60 measurements made. In Experiment 1(b) the balance was equipped with the new SmartPan™ design and the same tests were carried out. The settling time and the first stable digit were automatically recorded and average values are reported, as shown in Fig. 6.




Figure 5: Experiment 1(b): 5 mg readability precision balance with SmartPan™ on open lab bench

Weighing on a 5 mg readability balance with a SmartPan™ compared to standard weighing pan in a standard environment (on a open lab bench)




Figure 6: Results of Experiment 1: Effect of SmartPan™ on settling time (Speed) and repeatability (Precision) of a 5 mg readability precision balance on open lab bench (Standard Conditions)

The results of Experiment 1, shown in Fig. 6, demonstrate a decrease in the global average settling time and an improvement in the repeatability achieved using the new weighing pan design. This clearly shows that the SmartPan™ offers faster and more precise weighing performance than a standard weighing pan under standard laboratory conditions for weighing on a 5 mg readability precision balance.

3.2 Experiment 2: Weighing on a 5 mg readability balance under harsh conditions

Experiment 2 was carried out using a 5 mg readability precision balance (XPE6003SD5) placed directly beneath an air conditioning unit. The experimental set-up is illustrated in Fig. 7. Although this type of weighing environment is seen quite frequently, it is a less than ideal for achieving stable weighing results, due to disruptive air currents and is termed "Harsh conditions" for the purpose of this white paper. In this situation, the new weighing pan design has the potential to have a significant positive impact on the speed and precision of weighing.

In Experiment 2(a) the balance was equipped with a standard weighing pan. Test weights of 5 kg, 2 g, and 400 g were placed repeatedly on the weighing pan of the balance under test. Three replicates were each weighed 20 times, with a total of 60 measurements made.

In Experiment 2(b) the balance was equipped with the new SmartPan™ design and the same tests were carried out. The settling time and the first stable digit were automatically recorded and average values are reported, as shown in Fig. 8




Figure 7: Experiment 2(b): 5 mg readability precision balance with SmartPan™ exposed to air currents caused by air conditioning unit

Weighing on a 5 mg readability balance with a SmartPan™ compared to standard weighing pan in a HARSH environment (under air conditioning unit)




Figure 8: Results of Experiment 2: Effect of SmartPan™ on settling time (Speed) and repeatability (Precision) of a 5 mg readability precision balance under an air-conditioning unit (Harsh Conditions)

The results of Experiment 2 (Fig. 8) demonstrate a decrease in the global average settling time and a marked improvement in the repeatability achieved using the new weighing pan design. This proves that the SmartPan™ offers faster and significantly more precise weighing performance than a standard weighing pan under harsh laboratory conditions for weighing on a 5 mg readability precision balance.

3.3 Experiment 3: Weighing on a 1 mg readability balance under standard conditions

As the results of SmartPan™ tests for the 5 mg balance were so positive, the effect of the new weighing pan design was investigated further. More demanding weighing processes (in terms of balance type and conditions) were explored.

So, Experiment 3 was carried out using a 1 mg readability precision balance (XPE1203S) under the same standard environmental conditions as described in Experiment 1. The experimental set-up is illustrated in Fig. 9. In Experiment 3(a) the balance was equipped with a standard weighing pan AND a draft shield. Test weights of 500 g and 200 g were placed repeatedly on the weighing pan of the balance under test. Three replicates were each weighed 20 times, with a total of 60 measurements made.

In Experiment 3(b) the balance was equipped with the new SmartPan™ design, but had NO draft shield. The same tests were carried out. The settling time and the first stable digit were automatically recorded and average values are reported, as shown in Fig. 10.




Figure 9: Experiment 3(b): 1 mg readability precision balance with SmartPan™ on open lab bench

Weighing on a 1 mg readability balance with a SmartPan™ compared to standard weighing pan in STANDARD environment (on a open lab bench)




Figure 10: Results of Experiment 3: Effect of SmartPan™ on settling time (Speed) and repeatability (Precision) of a 1 mg readability precision balance on an open lab bench (Standard Conditions)

The results of Experiment 3 (Fig. 10) demonstrate a decrease in the global average settling time and an improvement in the repeatability achieved using the new weighing pan design. This proves that the SmartPan™ offers faster and more precise weighing performance than a standard weighing pan under standard laboratory conditions for weighing on a 1 mg readability precision balance.

3.4 Experiment 4: Weighing on a 1 mg readability balance under harsh conditions

Experiment 4 was carried out using a 1 mg readability precision balance (XPE1203S) placed inside a safety cabinet under a 0.35 m/s laminar air flow. Under these harsh conditions, weighing using a balance with a readability of 1 mg is especially challenging, and requires a draft shield to be fitted both for the conventional weighing pan and the SmartPan™. The experimental set-up is illustrated in Fig. 11.

In Experiment 4(a) the balance was equipped with a standard weighing pan and equipped with a draft shield. Test weights of 500 g and 200 g were placed repeatedly on the weighing pan of the balance under test. Three replicates were each weighed 20 times, with a total of 60 measurements made.

In Experiment 4(b) the balance was equipped with the new SmartPan™ design and a draft shield and the same tests were carried out. The settling time and the first stable digit were automatically recorded and average values are reported, as shown in Fig. 17.

Note: The balance was equipped with a draft shield in tests 4(a) and 4(b), because it is unrealistic to use a 3-place balance without a draft shield in a safety cabinet. However, the draft shield had the right-hand door open throughout these tests (as illustrated in Fig. 11).




Figure 11: Experiment 4(b): 1 mg readability precision balance with SmartPan™ in a safety cabinet

Weighing on a 1 mg readability balance with a SmartPan™ compared to standard weighing pan in HARSH environment (in a safety cabinet)




Figure 12: Results of Experiment 4: Effect of SmartPan™ on settling time (Speed) and repeatability (Precision) of a 1 mg readability precision balance in a safety cabinet (Harsh Conditions)

The results of Experiment 4 (Fig. 12) demonstrate a remarkable effect on both the global average settling time and the repeatability achieved using the new weighing pan design. This proves that the SmartPan™ offers much faster and significantly more precise weighing performance than a standard weighing pan under harsh laboratory conditions for weighing on a 1 mg readability precision balance.

4 Summary of Results

4.1 Results from 5 mg readability balance testing (Experiments 1 and 2)

The results of Experiments 1 and 2 carried out on a 5 mg readability precision balance (XPE6003SD5) are shown in Table 2 for Settling time (Speed) and Table 3 for Repeatability (Precision).



Figure 13: Precision balance with **standard weighing pan** (5 mg readability, XPE6003SD5) – used in Tests 1(a) and 2(a)



Figure 14: Precision balance with **SmartPan™** (5 mg readability, XPE6003SD5) – used in Tests 1(b) and 2(b)

4.1.1 Effect on Speed (Settling Time)

Experimental Conditions	Settling Time (Speed)		Effect of SmartPan™
	(a) Standard Pan	(b) SmartPan™	
Experiment 1 – Standard conditions (on open lab bench)	2.3 sec	1.7 sec	26% faster than standard pan
Experiment 2 – Harsh conditions (under air conditioning unit)	3.6 sec	1.7 sec	53% faster than standard pan

Table 2: Summary of Settling Time results when weighing on a SmartPan™ compared to standard pan under standard and harsh weighing conditions (average values)

This proves that the SmartPan™ has a very positive impact on the speed of the weighing process. A stable weighing result can be obtained 26% faster on an open lab bench using the SmartPan™ compared to the standard weighing pan. However, in harsh conditions where the air currents from an air conditioning unit are disrupting the weighing process, a stable weighing result can be obtained 53% faster using the SmartPan™ compared to the standard weighing pan. These results from Table 2 are displayed graphically in Fig. 15.

Effect of SmartPan™ on Settling Time (Speed) for 5 mg readability balance




Figure 15: Performance comparison of settling time (average values) for 1 mg balance with SmartPan™ compared to standard pan under standard and harsh weighing conditions

Fig. 15 shows that a stable weighing result can be obtained 1.4 times faster using the SmartPan™ under standard conditions. However, in harsh conditions a stable weighing result can be obtained 2.1 times faster using the SmartPan™ compared to the standard weighing pan.

4.1.2 Effect on Precision (Repeatability)

Experimental Conditions	Repeatability (Precision)		Effect of SmartPan™
	(a) Standard Pan	(b) SmartPan™	
Experiment 1 – Standard conditions (on open lab bench)	4.5 mg	2.9 mg	36% more precise than standard pan
Experiment 2 – Harsh conditions (under air conditioning unit)	14.7 mg	6.2 mg	57% more precise than standard pan

Table 3: Summary of Repeatability results when weighing on a SmartPan™ compared to standard pan under standard and harsh weighing conditions (average values)

This proves that the SmartPan™ has a very positive impact on the precision of the weighing process. Repeatability is 36% better on an open lab bench using the SmartPan™ compared to the standard weighing pan. However, in harsh conditions where the air currents from an air conditioning unit are disrupting the weighing process, repeatability is 57% better using the SmartPan™ compared to the standard weighing pan. The results from Table 3 are displayed in Fig. 16.

Effect of SmartPan™ on Repeatability (Precision) for 5 mg readability balance




Figure 16: Performance comparison of repeatability (average values) for 5 mg balance with SmartPan™ compared to standard pan under standard and harsh weighing conditions

Fig. 16 shows that a weighing precision (repeatability) is 1.5 times better using the SmartPan™ under standard conditions. However, in harsh conditions, the repeatability is 2.4 times better using the SmartPan™ compared to the standard weighing pan.

With the 5 mg readability balance, the new weighing pan design has a superior performance compared to the standard weighing pan, in terms of speed (settling time) and precision (repeatability) of the results, under the experimental conditions investigated. The difference is even more significant when weighing under harsh environmental conditions (such as under an air conditioning unit).

In this case, the weighing result can be obtained at least twice as fast (Fig. 15) and more than twice as precisely (Fig. 16) by using the new SmartPan™ instead of a standard weighing pan.

4.2 Results from 1 mg readability balance testing (Experiments 3 and 4)

The results of Experiments 3 and 4, carried out on a 1 mg readability precision balance (XPE1203S) are shown in Table 4 for Settling time (Speed) and Table 5 for Repeatability (Precision). The balance configurations used for these tests are shown in Figs. 17 – 19.



Figure 17: Precision balance with standard weighing pan and draft shield (1 mg readability, XPE1203S) – used in Tests 3(a) and 4(a)



Figure 18: Precision balance with SmartPan™ (1 mg readability, XPE1203S) – used in Test 3(b)



Figure 19: Precision balance with SmartPan™ and draft shield (1 mg readability, XPE1203S) – used in Test 4(b)

4.2.1 Effect on Speed (Settling Time)

Experimental Conditions	Settling Time (Speed)		Effect of SmartPan™
	(a) Standard Pan	(b) SmartPan™	
Experiment 3 – Standard conditions (on open lab bench)	2.7 sec	1.9 sec	30% faster than standard pan
Experiment 3 – Harsh conditions (in safety cabinet)	5.2 sec	1.9 sec	63% faster than standard pan

Table 4: Summary of Settling Time results when weighing on a SmartPan™ compared to standard pan under standard and harsh weighing conditions (average values)

This proves that the SmartPan™ has a very positive impact on the speed of the weighing process. A stable weighing result can be obtained 30% faster on an open lab bench using the SmartPan™ compared to the standard weighing pan. However, in harsh conditions where the air currents from working inside a safety cabinet are disrupting the weighing process, a stable weighing result can be obtained 63% faster using the SmartPan™ compared to the standard weighing pan. The results from Table 4 are displayed in Fig. 20.

Effect of SmartPan™ on Settling Time (Speed) for 1 mg readability balance




Figure 20: Performance comparison of settling time (average values) for 1 mg balance with SmartPan™ compared to standard pan under standard and harsh weighing conditions

Fig. 20 shows that a stable weighing result can be obtained 1.4 times faster using the SmartPan™ under standard conditions. However, in harsh conditions a stable weighing result can be obtained 2.7 times faster using the SmartPan™ compared to the standard weighing pan.

4.2.2 Effect on Precision (Repeatability)

Experimental Conditions	Repeatability (Precision)		SmartPan™ Improvement
	(a) Standard Pan	(b) SmartPan™	
Experiment 4 – Standard conditions (on open lab bench)	1.9 mg	0.8 mg	58% more precise than standard pan
Experiment 4 – Harsh conditions (in safety cabinet)	5.6 mg	0.8 mg	86% more precise than standard pan

Table 5: Summary of Repeatability results when weighing on a SmartPan™ compared to standard pan under standard and harsh weighing conditions (average values)

This proves that the SmartPan™ has a very positive impact on the precision of the weighing process. Repeatability is 58% better on an open lab bench using the SmartPan™ compared to the standard weighing pan. However, in harsh conditions where the air currents from an air conditioning unit are disrupting the weighing process, repeatability is 86% better using the SmartPan™ compared to the standard weighing pan. The results from Table 5 are displayed in Fig. 21.

Effect of SmartPan™ on Repeatability (Precision) for 1 mg readability balance




Figure 21: Performance comparison of settling time (average values) for 1 mg balance with SmartPan™ compared to standard pan under standard and harsh weighing conditions

Fig. 21 shows that a weighing precision (repeatability) is 2.7 times better using the SmartPan™ under standard conditions. However, in harsh conditions, the repeatability is 7 times better using the SmartPan™ compared to the standard weighing pan.

With the 1 mg readability balance, the SmartPan™ has a far superior performance compared to the standard weighing pan, in terms of speed (settling time) and precision (repeatability) of the results, under all experimental conditions investigated. The difference is even more significant when weighing under harsh environmental conditions (such as in a safety cabinet). In this case, the weighing result can be obtained more than twice as fast (Fig. 20) and seven times more precisely (Fig. 21) by using the new weighing pan design instead of a standard weighing pan.

5 Conclusion

Using a new design of weighing pan, with a quadruped structure, has proved to have a dramatic effect on the optimization of balance performance, in terms of speed and precision, compared to using a standard weighing pan.

Weighing on a 5 mg readability balance:

Experiments presented in this white paper show that in harsh weighing environments (such as under an air conditioning unit), using a 5 mg readability balance equipped with a SmartPan™ gives a significant improvement in speed and precision of weighing results.

- Weighing is more than twice as fast using a SmartPan™.
- Weighing results are more than twice as precise using a SmartPan™.

Weighing on a 1 mg readability balance:

In tests designed to push the limits even further, using a 1 mg readability (3-place) balance equipped with a SmartPan™ in harsh weighing environments (such as inside a safety cabinet) gives an even greater improvement in speed and precision of weighing results.

- Weighing is more than twice as fast using a SmartPan™.
- Weighing results are seven times more precise using a SmartPan™.

Weighing on an open lab bench without a draft shield

It is even possible to successfully use a 3-place balance on an open lab bench without a draft shield! In weighing comparisons on an open lab bench (standard conditions), the 3-place balance equipped with a SmartPan™ but no draft shield is still faster and more precise than the same balance equipped with a standard weighing pan WITH a draft shield.

- Weighing is almost 1.5 times faster using a SmartPan™, even without a draft shield.
- Weighing results are more than twice as precise using a SmartPan™, even without a draft shield.

Weighing pan design and ease of cleaning

The SmartPan™ has been specifically optimized to minimize the disruptive effects of air drafts and currents and has minimal flow resistance. It has static force distribution similar to a standard weighing pan structure for the precision balance. Furthermore, the new weighing pan design is much easier to clean than a grid structure.

Therefore, the new design of weighing pan is proven to be less sensitive to air drafts and presents an innovation with clear benefits to any users of precision balances.


6 Tips and Recommendations for Reliable Weighing

This white paper explains how balance users now have a simple, elegant solution for overcoming frequent problems due to air drafts in harsh environments. Nevertheless, there are useful standard tips and recommendations that the operator can observe to help obtain reliable weighing results.

Recommendations for location of the balance

The precision of weighing results is closely associated with the location of the balance. To ensure that your balance can work under the best conditions, please observe the following guidelines:


Weighing bench



- Stable (lab bench, lab table, stone bench). Your weighing bench should not sag when work is carried out on it and should transfer as few vibrations as possible.
- Antimagnetic (no steel plate).
- Protected against electrostatic charges (no plastic or glass).
- Wall or floor installation: The weighing bench should be fixed either to the floor or on the wall. Mounting the bench on both places at once transfers vibrations from wall and floor.
- Reserved for the balance.

The place of installation and the weighing bench must be stable enough that the balance display does not change when someone leans on the table or steps up to the weighing station. Do not use soft pads underneath, such as writing mats. It is better to position the balance directly over the legs of the bench, since the area is subject to the fewest vibrations.


Work room



- Vibration-free
- Free from drafts

Place the weighing bench in a corner of a room. These are the most vibration-free areas of a building. Ideally, the room should be accessed through a sliding door to reduce the influence of door movements.


Temperature



- Keep the temperature of the room as constant as possible. Weighing results are influenced by temperature! (Typical drift: 1 – 2 ppm/°C)
- Do not weigh near radiators or windows.


METTLER TOLEDO balances with "FACT" (fully automatic motorized self-calibration) can compensate virtually all the remaining temperature drift. For this reason, "FACT" should always be switched on.

Atmospheric humidity




- Ideally, the relative humidity (% RH) should be between 45% and 60%. Balances should never be operated above or below the measuring range of 20% to 80% RH.
- Constant monitoring is advisable with micro balances. Changes should be corrected whenever possible.

Light



- If possible, place the balance on a window-free wall. Direct sunlight (heat) will influence the weighing result.
- Place the balance a significant distance from lighting fixtures to avoid heat radiation. This especially applies to light bulbs. Use fluorescent tubes.

Air



- When possible do not place the balance in the air flow of air conditioners or devices with ventilators, such as computers or large laboratory devices.
- Place the balance a sufficient distance from radiators. In addition to the potential temperature drift, strong currents of air could interfere with operation of the balance. Do not place the balance next to a door.
- Avoid places with high foot traffic. Passers-by will usually create a draft at the weighing location.

7. Further Information on External Influences

7.1 eLearning course "Lab Balances: External Influences and Cleaning"

This course focuses on the six most common external influences which can seriously affect weighing results, and also offers solutions and tips to eliminate or at least reduce their effects. It also offers guidelines for the proper location of a balance in the lab, as well as proper methods for cleaning. During this course you will learn about:

- The effects of typical external influences acting on the balance
- Choose the proper location for a balance in your lab
- Taking simple corrective measures to eliminate or at least reduce these problems
- Best practices for cleaning your balance

Following successful completion of the final test, a qualification certificate will be automatically generated. This certificate can be printed and stored for internal training records, especially for GMP qualification, and is valid for three years.

For more information on this topic:

► www.mt.com/lab-elearning-influences

7.2 "Weighing the Right Way – Proper Weighing with Laboratory Balances"

This 36-page pdf guide is available to download free at:

► www.mt.com/lab-properweighing

8. About the Author



Dr. Peter Ryser is Professor of Microengineering at the Swiss Federal Institute of Technology in Lausanne. He has over three decades of research and teaching experience from various corporate and academic institutions. He was previously a Director at Siemens Building Technologies where he was responsible for R&D, product innovation and patents. Dr. Ryser has a Ph.D. in applied Physics from the University of Geneva, a Master's degree in Experimental Physics and a MBA.

Disclaimer

METTLER TOLEDO provides this White Paper as a service to its customers. In reading or making any use of this document, you acknowledge and agree to the following:

This document may contain inaccuracies and errors of both a substantive and/or typographical nature. METTLER TOLEDO does not guarantee the accuracy or completeness of the information or the reliability of any advice, opinion or statement in this document. If you rely on the information or any advice, opinion or statement, you are doing so at your sole risk. METTLER TOLEDO does not guarantee that this document or its contents are accurate, complete, reliable, truthful, current or error-free.

METTLER TOLEDO will not be liable for any decision made or action taken by you or others in reliance on the information in this document. METTLER TOLEDO and its affiliates are not liable for any Damages based on claims arising out of or in connection with your use of this document.

METTLER TOLEDO DOES NOT ASSUME ANY RESPONSIBILITY OR RISK FOR YOUR USE OF THE INFORMATION PROVIDED IN THIS DOCUMENT. THIS INFORMATION IS PROVIDED WITHOUT ANY REPRESENTATIONS, ENDORSEMENTS, OR WARRANTIES OF ANY KIND WHATSOEVER, EITHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF TITLE OR ACCURACY AND ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT, WITH THE SOLE EXCEPTION BEING WARRANTIES (IF ANY) WHICH CANNOT BE EXPRESSLY EXCLUDED UNDER APPLICABLE LAW. In no event will METTLER TOLEDO or its affiliates be liable for any Damages, even if METTLER TOLEDO is aware of the possibility of such Damages, arising in connection with the information provided herein. "Damages" includes but is not limited to all losses and all direct, indirect, incidental, special, consequential and punitive damages arising under a contract, tort or other theory of liability (including reasonable legal and accounting fees and expenses).

No part of this publication may be reproduced or distributed for any purpose without written permission from METTLER TOLEDO.

©2015 METTLER TOLEDO. All rights reserved.

www.mt.com/precision-balances

For more information

Mettler-Toledo GmbH

CH-8606 Greifensee, Switzerland
Tel. +41-44-944 22 11

Subject to technical changes
© 12/2015 Mettler-Toledo GmbH
Printed in Switzerland 30216952A