

COMING SOON!



BioPaddles™ Colony Identification App

# MAC/EMB

Code 5544

MacConkey Agar (**MAC**)

Eosin Methylene Blue Agar (**EMB**)

## USE:

Isolation and differentiation of Gram (-) enteric bacilli (**MAC**) / Coliform Testing / Recovery of Stressed Coliforms (**EMB**)

Side 1: MacConkey Agar (**MAC**) (pink)

(Side 1 is marked with an indented laser line)


Side 2: Eosin Methylene Blue Agar (**EMB**) (dark red)

## APPLICATION

Enumeration of Gram (-) rods, in particular, coliforms, and specifically, *E. coli*.

## PADDLE AGARS

**MacConkey Agar (MAC)** – Both selective and differential; used to differentiate between Gram negative bacteria while inhibiting the growth of most Gram positive bacteria. The medium also differentiates between lactose-fermenting coliforms Lac (+) and lactose non-fermenters Lac (-), which include potential pathogens. Addition to the nutrient agar base of bile salts and crystal violet will inhibit the growth of most Gram positive bacteria, making MacConkey agar selective. Lactose, a fermentable carbohydrate, and neutral red, a pH indicator, are added to differentiate the lactose positive coliforms from the potentially pathogenic lactose non-fermenters. When lactose is fermented, acid products lower the pH below 6.8, with the resulting colonial growth turning pinkish-red. If an organism is unable to ferment lactose, the colonies will be colorless. Bile salts mixture and crystal violet are the selective agents, inhibiting Gram positive cocci and allowing Gram negative organisms to grow. Sodium Chloride maintains the osmotic environment. Agar and a proprietary polymer are the solidifying agents.

**Eosin Methylene Blue Agar-Levine (EMB)** - Differential and slightly selective medium for the isolation, cultivation and differentiation of Gram negative enteric microorganisms (bacilli) isolated from both clinical and nonclinical specimens<sup>1</sup>. Gram (+) bacteria are inhibited. It is widely used for the examination of materials for the presence of coliforms. Eosin Methylene Blue Agar, abbreviated EMB, was developed by Holt-Harris and Teague.<sup>2</sup> This formula contains lactose and sucrose with two indicator dyes, Eosin Y and Methylene Blue. Levine modified the formula by removing sucrose and doubling the concentration of lactose.<sup>3 4</sup>

<sup>1</sup> Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (eds.). Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

<sup>2</sup> Holt-Harris, J. E., and O. Teague. 1916. A new culture medium for the isolation of *Bacillus typhosa* from stools. J. Infect. Dis. 18:596.

<sup>3</sup> Levin, M. 1918. Differentiation of *E. coli* and *A. aerogenes* on a simplified eosin-methylene blue agar. J. Infect. Dis. 23:43-47.

<sup>4</sup> Levin, M. 1921. Bacteria fermenting lactose, the significance in water analysis. Bull. 62. Iowa State College Eng. Exp. Sta., Ames, Iowa.

NOTE: EMB Agar is moderately inhibitory. Some staphylococci, streptococci and yeast may grow. They will appear as small, pinpoint colonies. Not all strains of *E. coli* produce a green metallic sheen. The presence of the green metallic sheen is not diagnostic for *E. coli*.<sup>5</sup>

## CULTURE CONTROLS

10-300 inoculum (CFU)


	MAC Agar	EMB Agar
<i>Enterococcus faecalis</i>	INHIBITED	PARTIAL INHIBITION (colorless colonies)
<i>Escherichia coli</i>	GROWTH	GROWTH (blue-black w/ green metallic sheen)
<i>Proteus mirabilis</i>	GROWTH (partial inhibition of swarming)	GROWTH
<i>Salmonella typhimurium</i>	GROWTH	GROWTH

## STORAGE / EXPIRATION

Store tightly sealed BioPaddles™ in a cool, dry location (less than 68°F/20°C). Avoid temperature changes. BioPaddles™ may be refrigerated, but it is not necessary. Do not freeze. If freezing occurs, thaw (3-6 hours) under refrigeration temperatures (40°F; 4.4°C). Freezing can promote excess water loss and variation in media surface due to crystal formation. The average shelf-life is one year. Refer to expiration date (See: BBE stamped on vial). Discard if paddle agar appears oxidized (darkens from expected color). The expiration date applies to the medium in an intact container when stored as directed.

## SAMPLING

**Liquids:** Twist to remove paddle from vial. Fill vial to 40 mL fill line with the liquid to be sampled. The 40 mL volume can be used to calculate Total Viable Count (TVC) and/or Total Colony Count (TCC). Replace paddle. Allow a contact time of 15 seconds. Remove the paddle. Empty the vial. Replace the paddle in the vial.



**Surfaces:** Twist to remove paddle from vial. Allow the paddle surface (10 cm<sup>2</sup>) to come into physical contact with the test surface. Recovery rate is about 50%. To insure an accurate recovery, gently sweep (or touch) the paddle to cover a 20 cm<sup>2</sup> area. Replace paddle in vial


## INCUBATION

Incubate at 35°C ±2°C for 18-24 hours.




<sup>5</sup> MacFaddin, J. F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1. Williams & Wilkins, Baltimore, MD.

## COLONY MEASURING

Each BioPaddles™ paddle has molded media attachment points that are 4mm in length (point-to-point). This feature provides a useful guidepost to estimating nearby colony size.







## IDENTIFICATION

ORGANISM		MAC			EMB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Bacillus spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole (-)</li> <li>• Oxidase (-)</li> <li>• Catalase (+)</li> <li>• Gram (+) Rod</li> </ul>	INHIBITED	---	---	INHIBITED	---	---
<i>Candida albicans</i>	<ul style="list-style-type: none"> <li>• Catalase (+)</li> <li>• Actinomycetes</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Cream</li> <li>• CVEG</li> <li>• 1-2mm</li> </ul>		INHIBITED	---	---
<i>E. coli</i>	<ul style="list-style-type: none"> <li>• Lactose (+)</li> <li>• Indole (+) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Pink / Red</li> <li>• CVEG</li> <li>• 2-4mm</li> </ul>		+++	<ul style="list-style-type: none"> <li>• Blue-black bulls eye w/ green metallic sheen</li> <li>• CVEG</li> <li>• 2-4mm</li> </ul>	

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles. All rights reserved.




LaMotte\_BioPaddles\_MAC\_EMB

ORGANISM		MAC			EMB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Enterobacter aerogenes</i>	<ul style="list-style-type: none"> <li>• Lactose (+)</li> <li>• Indole (+) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Colorless</li> <li>• Thick, round, raised to low-convex; Spreading</li> <li>• 2-3mm</li> </ul>		+++	<ul style="list-style-type: none"> <li>• Pink to purple colonies often with a central dark purple dot (bulls eye); mucoid</li> <li>• CVEG</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING
<i>Enterococcus spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (-) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (+) Sphere</li> </ul>	PARTIAL - COMPLETE INHIBITION	---	---	PARTIAL - COMPLETE INHIBITION	---	---
<i>Klebsiella spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (+)</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (+) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Colorless / Light Pink</li> <li>• Spreading</li> <li>• 4-5mm</li> </ul>		+++	<ul style="list-style-type: none"> <li>• Pink to purple colonies often with a central dark brown dot (bulls eye); mucoid</li> <li>• CVEG</li> <li>• 4-5mm</li> </ul>	IMAGE PENDING
<i>Proteus spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• SEE: INDOLE ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (+) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Colorless; Spreading / swarming</li> <li>• 2-3mm</li> </ul>		+++	<ul style="list-style-type: none"> <li>• Colorless</li> <li>• CVEG</li> <li>• 0.5-1mm</li> </ul>	
<i>Pseudomonas aeruginosa</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole (-) ◆</li> <li>• Oxidase (+) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> <li>• Fluoresces blue under long-wave UV light (400-nm)</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Cream / off-white</li> <li>• Irregular / spreading</li> <li>• Diffusible green-blue pigment</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING	+++	<ul style="list-style-type: none"> <li>• Colorless</li> <li>• Irregular / spreading / glossy</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles. All rights reserved.

LaMotte\_BioPaddles\_MAC\_EMB

ORGANISM		MAC			EMB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Pseudomonas fluorescens</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole (-) ◆</li> <li>• Oxidase (+) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> <li>• Fluoresces blue-green under long-wave UV light (400-nm)</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Colorless to pink</li> <li>• Center target</li> <li>• Glossy</li> <li>• Umbonate</li> <li>• 2-4mm</li> </ul>		++	<ul style="list-style-type: none"> <li>• Colorless to pink</li> <li>• Center target</li> <li>• Glossy</li> <li>• Umbonate</li> <li>• 2-4mm</li> </ul>	
<i>Salmonella typhimurium</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Clear, transparent or amber</li> <li>• Circular, Smooth, Flat</li> <li>• 2-3mm</li> </ul>	IMAGE PENDING	++	<ul style="list-style-type: none"> <li>• Colorless to gray, circular</li> <li>• FED</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING
<i>Salmonella epidermidis</i>	<ul style="list-style-type: none"> <li>• Lactose (-) ◆</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> </ul>	PARTIAL - COMPLETE INHIBITION	---	---	++	<ul style="list-style-type: none"> <li>• Colorless to gray,</li> <li>• Circular</li> <li>• Umbonate</li> <li>• 2-4mm</li> </ul>	
<i>Serratia spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (-) ◆</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (+) ◆</li> <li>• Gram (-) Rod</li> </ul>	++	<ul style="list-style-type: none"> <li>• Pink / Red</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING	++	<ul style="list-style-type: none"> <li>• Colorless</li> <li>• 1-3mm</li> </ul>	MAGE PENDING
<i>Shigella spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (-)</li> <li>• Indole - mixed ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Urease (-) ◆</li> <li>• Gram (-) Rod</li> </ul>	+++	<ul style="list-style-type: none"> <li>• Colorless or transparent / faintly pink</li> <li>• Raised</li> <li>• 1-3mm</li> </ul>	IMAGE PENDING	++	<ul style="list-style-type: none"> <li>• Colorless to pink; circular</li> <li>• FED</li> <li>• 2-4mm</li> </ul>	IMAGE PENDING

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles. All rights reserved.

LaMotte\_BioPaddles\_MAC\_EMB

ORGANISM		MAC			EMB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Staphylococcus aureus</i>	<ul style="list-style-type: none"> <li>• Lactose (-) ◆</li> <li>• Indole (-) ◆</li> <li>• Oxidase (-) ◆</li> <li>• Catalase (+) ◆</li> <li>• Gram (+) Sphere</li> </ul>	PARTIAL - COMPLETE INHIBITION	---	---	PARTIAL - COMPLETE INHIBITION	---	---
<i>Streptococcus spp.</i>	<ul style="list-style-type: none"> <li>• Lactose (+)</li> <li>• Indole (+)</li> <li>• Oxidase (-)</li> <li>• Catalase (-)</li> <li>• Gram (+) Sphere</li> </ul>	PARTIAL - COMPLETE INHIBITION	---	---	PARTIAL - COMPLETE INHIBITION	---	---
Gram (+) Bacteria		PARTIAL - COMPLETE INHIBITION	---	---	PARTIAL - COMPLETE INHIBITION	---	---

NOTE: EMB Agar is moderately inhibitory. Some staphylococci, streptococci and yeast may grow. They will appear as small, pinpoint colonies. Not all strains of *E. coli* produce a green metallic sheen. The presence of the green metallic sheen is not diagnostic for *E. coli*.<sup>6</sup>

<sup>6</sup> MacFaddin, J. F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1. Williams & Wilkins, Baltimore, MD.




## ENUMERATION

### MAC




Very Light    Light    Moderate

### EMB



Very Light    Light    Moderate

### TVC / TCC



100    1000    10000    100000    1000000

Approximate Colony Count per 100 mL

### **TVC/TCC** (Total Viable Count/ Total Colony Counts)

Colony Counts < 1000

Count colonies

TVC/TCC Count = Count x 2.5

Colony Counts > 1000

Use chart

TVC/TCC Count = Count x 2.5

(Based on a 40 mL sample)



Example:  
Inoculated NUT/TTC  
paddle showing  
approximately 1000  
CFU/100 mL.

## DISPOSAL

Twist to remove paddle from vial. Fill vial to 40 mL fill line with 1:9 dilution of household bleach (5.25% sodium hypochlorite). Replace paddle in vial. Allow 15 minute contact time. Remove paddle. Discard bleach solution. Replace paddle in vial and dispose. Alternatively, loosen cap and microwave for 30 seconds, autoclave, or incinerate.

## GLOSSARY:

<b>CVEG</b>	Convex, Entire, Glossy
<b>FED</b>	Full, Entire, Dull
<b>Catalase</b>	Cat (+) contains catalase enzymes that degrade cellular H <sub>2</sub> O <sub>2</sub> .
<b>Lactose</b>	Lac (+) bacteria can ferment available lactose in the medium producing an acid which lowers the pH. Lac (-) are non-fermenting.
<b>Indole</b>	Biochemical test to determine the ability of an organism to split indole from the amino acid tryptophan. <i>P. vulgaris</i> is indole (+) while <i>P. mirabilis</i> is indole.
<b>Oxidase</b>	OX (+) contains cytochrome c oxidase. In contact with an indicator turns dark blue if OX (+); colorless if OX(-).
<b>Urease</b>	UR (+) presence of enzyme urease which hydrolyzes urea into CO <sub>2</sub> & NH <sub>4</sub> .
<b>β-D-Glucoronidase</b>	Glu (+) Found in 97% of <i>E. coli</i> strains. The presence of <i>E. coli</i> is determined when both β-D-Glucoronidase and Indole are (+), and the organism is Gram (-).
<b>Gram</b>	Gram reaction