



Cole-Parmer®

SAFESKIN® Nitrile Gloves

SAFESKIN NITRILE GLOVES FROM COLE-PARMER ARE A SUPERIOR ALTERNATIVE TO LATEX



These High Quality gloves from Kimberly-Clark are puncture - resistant, and provide a comfortable feel to the user - especially when worn over longer periods, without compromising on dexterity.

SAFESKIN® gloves are purple so you can easily identify them—reduce your chances of putting on latex gloves. These 100% nitrile gloves are soft, latex-free, powder-free, and non-sticky. Low chemical content reduces risk of allergic reactions. Textured fingertips enhance wet or dry grip. The gloves are Ambidextrous & have a Beaded cuff.

Specifications :

- Material : Nitrile
- Thickness : 5 mil
- Compliance : ASTM D 6319, ISO 9002

Catalog Number	Length	Size	Qty/bx	Box/cs	Price/cs
YK-33466-04	9.5"	Small	100	10	Rs. 9,000
YK-33466-05		Medium	100	10	Rs. 9,000
YK-33466-06		Large	100	10	Rs. 9,500
YK-33466-07	12"	Small	50	10	Rs. 7,000
YK-33466-08		Medium	50	10	Rs. 7,000
YK-33466-09		Large	50	10	Rs. 7,500

Prices are Ex-Warehouse, Taxes Extra, Valid upto Dec 2008

* Trademark of Kimberly-Clark Worldwide, Inc.