

Animal Origin Statement

Date: 04/03/2012

Catalog #'(s): 06402-35

Cole-Parmer certifies that, according to the manufacturer, the product described meets the specifications set forth in this Certificate:

- **The Manufacturer does not add any animal derived ingredients into the Norprene[®] tubing product formulation. However, during the manufacture of some of the product components, our suppliers have informed us that some contain tallow-based additives derived from bovine material. They certify that the materials originate from the United States, Canada and Mexico.**
- **Additionally, upon review of their certification statements regarding the processing conditions,**

we can state that they meet the rigorous processing requirements published in the *Official Journal of the European Communities*, (2004/C 24/03) and the *EMEA/410/01, Rev. 3*, dated 23 Jun 04. Please see page 9, section 6.4 “Tallow Derivatives.”

This certification applies only to the condition of the above-described product in its unopened package. Cole-Parmer cannot and does not assume any responsibility for the failure of the product to meet this document that results from any handling or use subsequent to opening the package.

Cole-Parmer relies upon the information provided by the manufacturer of this product(s), and Cole-Parmer does not assume any responsibility for inaccuracies in the manufacturer's statements.

Cole-Parmer expressly warrants that the product listed above conforms only to the specifications herein, subject to the disclaimers listed above. Except as otherwise set forth in writing by Cole-Parmer, there are no other warranties, express or implied, that apply to this certificate or the subject product(s).

COLE - PARMER