

Cole-Parmer® Blades and Scalpels



Animal Physiology

Dissecting Blades

Ideal for dissection these blades are available in carbon or stainless steel and individually sealed in foil packs. Use with handles at right.



Photo	Blade type	Qty/ bx	Carbon steel blades		Stainless steel blades	
			Catalog number	Price/ bx	Catalog number	Price/ bx
	#10	100	GH-10822-02		GH-10822-03	
	#11	100	GH-10822-04		GH-10822-05	
	#12	100	GH-10822-06		GH-10822-07	
	#15	100	GH-10822-08		GH-10822-09	
	#20	100	GH-10822-10		GH-10822-11	
	#22	100	GH-10822-12		GH-10822-13	

Dissecting Blade Handles

Use with blades at left.



Description	Use with blade type	Handle length	Catalog number	Price/ea
Siegel, stainless steel	#10 through #15	6.5" (16.5 cm)	GH-10822-20	
No. 3, stainless steel	#10 through #15	5" (12.7 cm)	GH-10822-22	
No. 3, economy, nickel-plated	#10 through #15	5" (12.7 cm)	GH-10822-24	
No. 4, stainless steel	#20 through #25	5.375" (13.6 cm)	GH-10822-26	
No. 4, economy, nickel-plated	#20 through #25	5.375" (13.6 cm)	GH-10822-28	
No. 7, stainless steel	#10 through #15	6.5" (16.5 cm)	GH-10822-30	
No. 6, plastic with stainless steel tip	#20 through #25	5.25" (13.3 cm)	GH-10822-32	

Disposable Dissecting Scalpels

Get superior dissecting with stainless steel razor-sharp blades. Single-use scalpels are individually packaged in sealed pouches. Each includes a blade safety sleeve to protect you and the blade.




Photo	Blade type	Catalog number	Qty/bx	Price/bx
	#10	GH-10818-31	10	
	#11	GH-10818-33	10	
	#12	GH-10818-35	10	
	#15	GH-10818-37	10	
	#20	GH-10818-39	10	
	#21	GH-10818-41	10	
	#22	GH-10818-43	10	
	#23	GH-10818-45	10	

Disposable Safety Dissecting Scalpels

These single-use scalpels come with fully retractable blades. Blades lock securely into place—eliminating accidents.



Photo	Blade type	Catalog number	Qty/bx	Price/bx
	#10	GH-10818-90	10	
	#11	GH-10818-92	10	
	#12	GH-10818-93	10	
	#15	GH-10818-94	10	

Warning—These products are not approved or intended for, and should not be used for, medical, clinical, surgical or other patient oriented applications.

Find even more online! Go to ColeParmer.com to search more than 100,000 products.