

Ansell Neoprene Apron

Won't tear

- Resists acids and oils—neoprene design is made for frequent use
- Provides additional chest protection with wider bib



Thickness	Size (L x W)	Color	Catalog number	Price per apron	
				1 ea	12 ea
Lightweight	45" x 35" (114.3 x 89.0 cm)	Yellow	GH-86479-03		

PVC/Nylon Apron

Protects against acids and caustics

- No break throughs—puncture resistant 20-mil thick design
- Feature cotton ties
- Made in USA



Thickness	Size (L x W)	Qty/cs	Catalog number	Price/cs
20 mil	50" x 35" (127.0 x 89.0 cm)	12	GH-81501-11	

Cryo-Apron® Aprons

Withstand temperatures to -256°F (-160°C)

- Waterproof—keeps liquids and cold equipment from contacting the torso
- Stands up to use with LN₂ vapor phase and liquified natural gas (LNG)
- Breathable nylon on sides, resists most chemicals
- Machine washable



Adjustable plastic fasteners at waist and neck snap easily into place, even when you're wearing gloves.

Size (L x W)	Catalog number	Price/ea
36" x 24" (91.4 x 60.9 cm)	GH-09113-36	
42" x 24" (106.6 x 60.9 cm)	GH-09113-38	
48" x 24" (121.9 x 60.9 cm)	GH-09113-40	

Ansell Vinyl Aprons

Won't trap contaminants

- Seamless design eliminates areas for contaminants to gather
- No purchasing extra protection—sleeves protect shoulders and arms too
- Heavy weight construction—40-mil PVC grommets secure for strength and durability



Description	Thickness	Size	Cat. no.	Price/ea*
Coat apron with open sleeves	6 mil	X-large	GH-81620-52	
Coat apron with elastic sleeves	8 mil	Medium	GH-81620-54	
		Large	GH-81620-56	
		X-large	GH-81620-58	
Bib apron	6 mil	45" L x 35" W (114.3 x 89.0 cm)	GH-81620-60	

*Quantity discounts available—call for details.

DuPont™ Frocks

Superior particulate and light barrier protection

IsoClean® Tyvek® Frocks have five snap closures, elastic wrists, and a three-snap adjustable mandarin collar. Serged seams.



Size	Qty/cs	Catalog number	Price/cs
Medium	30	GH-33683-02	
Large		GH-33683-04	
X-large		GH-33683-06	



Pro/Clean® Frocks have a zipper front and elastic wrists. Serged seams.

Size	Qty/cs	Catalog number	Price/cs
Small	30	GH-33683-40	
Medium		GH-33683-42	
Large		GH-33683-44	
X-large	30	GH-33683-46	
XX-large		GH-33683-48	

Wait!

There's More at
ColeParmer.com

Find all the
product specs
you need!