

Polyethylene-Coated Polypropylene Lab Coats

Get splash protection without the bulk

These polyethylene (PE)-coated polypropylene (PP) lab coats are a low-cost disposable alternative to Tyvek® coats. The durable, lightweight, and breathable fabric is ideal for general-purpose laboratory applications. The PE-coated spunbonded PP provides excellent splash protection from water-based liquids.



Size	Length	Qty/cs	Catalog number	Price per case	
				1 cs	10 cs
Medium	41 1/2" (105.4 cm)	30	GH-86227-12		
Large			GH-86227-14		
X-large			GH-86227-16		

Wait! There's More...

Protect your hands with our full selection of gloves. [see pages 356-367.](#)

Isolation Barrier Gowns

Full protection at a low cost

Disposable gowns provide full frontal protection from particulates or light splashes

Lightweight, breathable gowns feature ties in back of neck and at waist for a secure, comfortable fit. Elastic cuffs at wrists keep the sleeves firmly in place and gather softly under gloves. Select 100% spunbound polypropylene blue gowns for durable, particulate protection. Select polyethylene-coated polypropylene yellow gowns for added protection against water-based liquids and certain chemicals.



Material	Length	Qty/cs	Catalog number	Price per case	
				1 cs	10 cs
Polypropylene	45 1/2" (115.5 cm)	50	GH-86566-00		
Polyethylene-coated polypropylene			GH-86566-10		

Wait!

There's More at ColeParmer.com

Find all the product specs you need!

DuPont™ Tyvek® Coveralls

Tear-free, whole body protection

- Rugged fabric resists tears, snags, punctures, and abrasions
- Easy-zip front for quicker donning

Non-woven fabric provides protection against particulates. Tyvek is chemically and biologically inert, low linting, and static-free until laundered.

Compliance: sizing requirements of ANSI/ISEA 101-1996



86225-01

Size	Catalog number	Price per unit	
		1 ea	25 ea†
Coveralls			
Medium	GH-86225-01		
Large	GH-86225-02		
X-large	GH-86225-03		
XX-large	GH-86225-04		
Coveralls with elastic wrists			
Medium	GH-86225-11		
Large	GH-86225-12		
X-large	GH-86225-13		
XX-large	GH-86225-14		
Coveralls with elastic wrists and hood			
Medium	GH-86225-21		
Large	GH-86225-22		
X-large	GH-86225-23		
XX-large	GH-86225-24		
Coveralls with elastic wrists, hood, and boots			
Medium	GH-86225-35		
Large	GH-86225-36		
X-large	GH-86225-37		
XX-large	GH-86225-38		

†Discount: Save 10% on 25 or more coveralls of any combination of sizes and styles

Static Control Lab Coats

Get control of static

- Made from non-linting 94% polyester and 6% carbon suffused
- Material maintains static resistance even after multiple washings

Coats close easily with snap closures. Adjustable sleeves feature a grounding snap at each wrist. Each coat has two lower pockets and one left breast pocket.



Compliance: ESD STM2.1-1997

Size	Length	Catalog number	Price
Small	33" (83.8 cm)	GH-86297-40	
Medium		GH-86297-41	
Large		GH-86297-42	
X-large		GH-86297-43	
Small	41 1/2" (105.4 cm)	GH-86297-44	
Medium		GH-86297-45	
Large		GH-86297-46	
X-large		GH-86297-47	