

Zip-Lip Bags, LDPE

Won't let in contaminants

- No need to worry about dust or moisture getting to your dry samples
- Single track seal not for use with liquids
- Easily identify samples by using bags with write-on labels
- FDA-compliant



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
Zip-lip bags without labeling area					
2 x 3	5.1 x 7.6	2	GH-60101-00	1000	
3 x 5	7.6 x 12.7		GH-60101-10		
4 x 4	10.2 x 10.2		GH-60101-20		
4 x 6	10.2 x 15.2	2	GH-60101-30	1000	
5 x 8	12.7 x 20.3		GH-60101-40		
6 x 6	15.2 x 15.2		GH-60101-50		
6 x 9	15.2 x 22.9	4	GH-60101-60	1000	
8 x 10	20.3 x 25.4	2	GH-60101-70		
9 x 12	22.9 x 30.5	3	GH-60101-80		
12 x 15	30.5 x 38.1	2	GH-60101-90		
Zip-lip bags with labeling area					
2 x 3	5.1 x 7.6	2	GH-60102-00	1000	
3 x 5	7.6 x 12.7		GH-60102-10		
4 x 6	10.2 x 15.2		GH-60102-30		
5 x 8	12.7 x 20.3	2	GH-60102-40	1000	
6 x 9	15.2 x 22.9		GH-60102-60		
Assortment pack, zip-lip bags w/o label area; 100 ea. of 60101-10, -40, -50, -60, & -80					
—	—	—	GH-06503-90	500	

Heavy-Duty Bag, PE

No more weak bags

- Conforms to ordinances for disposal of dead animals
- Ideal for animal research labs, and humane societies



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases


Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
36 x 50	91.4 x 127	4	GH-06504-10	25	

Autoclavable Bags, PP

No worrying about contaminants

- Stands up to your demands—sturdy 1.5-mil polypropylene construction.
- Autoclavable at 250°F (121°C)

NEW



Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
4 x 6	10.2 x 15.2	1.5	GH-06505-01	100	
4 x 10	10.2 x 25.4	1.5	GH-06505-03	100	
6 x 10	15.2 x 25.4	1.5	GH-06505-05	100	
6 x 12	15.2 x 30.5	1.5	GH-06505-07	100	
10 x 16	25.4 x 40.6	1.5	GH-06505-09	100	
12 x 18	30.5 x 45.7	1.5	GH-06505-11	100	
12 x 24	30.5 x 70	1.5	GH-06505-13	100	
16 x 18	40.6 x 45.7	1.5	GH-06505-15	100	
18 x 23	45.7 x 58.2	1.5	GH-06505-17	100	
22 x 30	55.9 x 76.2	1.5	GH-06505-19	100	

[GH-06831-15](#) Plastic twist ties, 4"L (10.2 cm). Pack of 1000

Bags, LDPE

Stands up to your needs with strong bottom seal

- Seal with portable heat sealer or bag ties



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
3 x 5	7.6 x 12.7	2	GH-06505-02	1000	
3 x 8	7.6 x 20.3		GH-06505-04		
4 x 5	10.2 x 12.7		GH-06505-06		
5 x 7	12.7 x 17.8		GH-06505-08		
6 x 9	15.2 x 22.9		GH-06505-10		
6 x 18	15.2 x 45.7	3	GH-06505-12	500	
8 x 10	20.3 x 25.4		GH-06505-14		
10 x 12	25.4 x 30.5		GH-06505-16		
12 x 18	30.5 x 45.7	2	GH-06505-18	500	
18 x 24	45.7 x 61		GH-06505-20		
18 x 30	45.7 x 76.2		GH-06506-22		

Portable heat sealer 03014-00



[GH-03014-00](#) Portable heat sealer; 110 VAC, 60 Hz

[GH-03014-02](#) Portable heat sealer; 240 VAC, 50/60 Hz

[GH-06831-15](#) Plastic tie twists, 4"L (10.2 cm). Pack of 1000

Amber Bags, PE

Don't let light get to your samples

- Protect photo-sensitive samples during storage and transport
- Choose resealable zip-top or open-end styles



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
Resealable zip-top bags					
2 x 3	5.1 x 7.6	3	GH-60103-00	100	
3 x 5	7.6 x 12.7		GH-60103-01		
4 x 6	10.2 x 15.2		GH-60103-02		
5 x 8	12.7 x 20.3		GH-60103-03		
8 x 8	20.3 x 20.3		GH-60103-04		
8 x 14	20.3 x 35.6	3	GH-60103-05	100	
9 x 12	22.9 x 30.5		GH-60103-06		
12 x 12	30.5 x 30.5		GH-60103-07		
3 x 9	7.6 x 22.9		GH-60103-08		
Open-end bags					
5 x 7	12.7 x 17.8	2	GH-60103-09	100	
6 x 10	15.2 x 25.4		GH-60103-11		
8 x 14	20.3 x 35.6		GH-60103-13		
22 x 18	55.9 x 45.7		GH-60103-14		