

**Cole-Parmer® Economical Square Carboys, HDPE**

Handle and pour contents easily

- Molded-in handle and indented bottom for easy grabbing
- No leaks with leakproof HDPE threaded cap
- Get easier dispensing by drilling and attaching optional faucet
- Autoclavable; withstand temperature from 68 to 248°F (20 to 120°C)



Buy More...Save More!

- 20% off on 120+
- 15% off on 30 to 119
- 10% off on 6 to 29 packs or cases



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	63 mm	<a href="#">GH-06031-80</a>	1	
20 (5)		<a href="#">GH-06031-82</a>		
10 (2½)	63 mm	<a href="#">GH-06031-84</a>	4	
20 (5)		<a href="#">GH-06031-86</a>		

Accessories

[GH-06031-88](#) HDPE faucet, ¾"

**Thermo Scientific™ Nalgene™ Vacuum Carboys with Handles, PP**

Withstand vacuum for up to eight hours

- Save money compared to stainless steel and fiberglass-reinforced containers
- No leaks with polypropylene leakproof cap
- Autoclavable



Buy More...Save More!

- 20% off on 120+
- 15% off on 30 to 119
- 10% off on 6 to 29 packs or cases

Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	83B	<a href="#">GH-06257-50</a>	1	
20 (5)		<a href="#">GH-06257-60</a>		

Accessories

[GH-06258-30](#) Replacement cap; PP, 83B. Pack of 2

**Baritainer® Jerry Cans, HDPE/Quoral®**

More environmentally friendly design

- Great alternative to glass, metal, or fluorinated containers
- Safely transport chemicals, solvents, adhesives, and cleaning solutions
- Use in temperatures from -148 to 248° (-100 to 120°C)



Buy More...Save More!

- 20% off on 120+
- 15% off on 30 to 119
- 10% off on 6 to 29 packs or cases

Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
5 (1¼)	53 mm	<a href="#">GH-06069-01</a>	1	
10 (2½)	53 mm	<a href="#">GH-06069-03</a>	1	
20 (5¼)	63 mm	<a href="#">GH-06069-05</a>	1	

**Thermo Scientific™ Nalgene™ Safety Dispensing Jugs, LDPE**

Just tip and pour

- Get a better grip with top and bottom handles for easier pouring
- Leakproof design for better transport
- Don't misplace spout plug—attached for easier locating



Buy More...Save More!

- 20% off on 120+
- 15% off on 30 to 119
- 10% off on 6 to 29 packs or cases

Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	83B	<a href="#">GH-06077-22</a>	1	
20 (5)		<a href="#">GH-06077-52</a>		

Accessories

[GH-06258-30](#) Replacement cap; PP, 83B. Pack of 2