


Rubbermaid® Storage Containers with Handles, HDPE and PP

No more difficult lifting

- Easy-lifting grip handles, 22-qt containers have removable handles
- Keep track of contents with “date-content” label area
- Feature graduations in quarts and liters and optional tight fitting lids
- Resist temperatures from –40 to 150°F (–40 to 66°C)
- Dishwasher safe, NSF certified



HDPE
63511-01



PP
63512-13

Lids sold separately



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Capacity L (qt)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
White high-density polyethylene (HDPE) containers					
2 (2)	4 x 8½	10.2 x 21.6	GH-63511-01	12	
4 (4)	7¾ x 8½	19.7 x 21.6	GH-63511-03		
6 (6)	7½ x 10	19.4 x 25.4	GH-63511-05	12	
8 (8)	10½ x 10	27.0 x 25.4	GH-63511-07		
12 (12)	8½ x 13½	20.6 x 33.3	GH-63511-09	6	
18 (18)	11½ x 13½	30.2 x 33.3	GH-63511-11		
22 (22)	14 x 13½	35.6 x 33.3	GH-63511-13		
Clear polypropylene (PP) containers					
2 (2)	4 x 8½	10.2 x 21.6	GH-63512-01	12	
4 (4)	7¾ x 8½	19.7 x 21.6	GH-63512-03		
6 (6)	7½ x 10	19.4 x 25.4	GH-63512-05	12	
8 (8)	10½ x 10	27.0 x 25.4	GH-63512-07		
12 (12)	8½ x 13½	20.6 x 33.3	GH-63512-09	6	
18 (18)	11½ x 13½	30.2 x 33.3	GH-63512-11		
22 (22)	14 x 13½	35.6 x 33.3	GH-63512-13		
Yellow low-density polyethylene (LDPE) lids					
For 2- and 4-L containers			GH-63511-53	12	
For 6- and 8-L containers			GH-63511-54	12	
For 12-, 18-, 22-L containers			GH-63511-55	6	

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.co.uk/ChemChart



Keep track of contents by recording date and material on lid.

Round Dipping Baskets, PP

Easily remove contents from solution

- ½" perforations for liquid release
- Easy-grip handle for ease of handling
- Use with crocks sold at right
- NSF-certified



Round dipping basket
06717-50
(shown with heavy crock)



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Fits crock model no.	Dimensions (H x ID)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
06724-20	8 x 8	20.3 x 20.3	GH-06717-20	1	
06724-40	10 x 9	25.4 x 22.9	GH-06717-30		
06724-60	13 x 9	33.0 x 22.9	GH-06717-50		
06724-80	17 x 11½	43.2 x 29.2	GH-06717-70		

Heavy Crocks, HDPE

Replace heavier, older style crocks

- Sturdy design—one piece molded HDPE
- Use with dipping baskets sold at left for easy item retrieval
- Include cover for easy storage
- NSF certified; USDA approved



Heavy crock
06724-60
(cover included)
shown with
round dipping
basket



Buy More...Save More!

20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Capacity L (gal.)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
12 (3.25)	11¾ x 10½	29.8 x 25.7	GH-06724-20	1	
15 (4)	13¼ x 12	33.7 x 30.5	GH-06724-40		
22 (5.75)	16½ x 11	41.9 x 27.9	GH-06724-60		
40 (10.25)	21 x 13	53.3 x 33.0	GH-06724-80		

Square Utility Pan, 304 SS (medium gauge)

Holds up to tough demands

- Sturdy, stainless steel construction for longer life
- NSF certified

Buy More...Save More!

20% off on 120+
15% off on 30 to 119
10% off on 6 to 29 packs or cases

Capacity L (qt)	Dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
1.25 (1¼)	6½ x 6½ x 2½	15.6 x 15.6 x 7.3	GH-07286-01	1	

