

## Cotton Inspection Gloves

### Stop leaving fingerprints

- Minimize cross-contamination from human hands
- Washable and reversible for reusing
- Ideal for liner between liquid-proof gloves and hands
- Ambidextrous



Size	Qty/pk	Catalog number	Price per pack	
			1 pk	10 pk
<b>Lightweight gloves</b>				
Women's	12 pairs	<a href="#">GH-86369-00</a>		
Men's		<a href="#">GH-86369-02</a>		
<b>Medium-weight gloves</b>				
Men's	12 pairs	<a href="#">GH-86369-05</a>		

## 100% Nylon Knit Glove Liners

### Prevent skin irritation

- Prevent irritation during extended glove wearing
- Comfortable, nylon knit is dust-free for cleanroom use
- Reusable—just wash and sterilize



Size	Qty/pk	Catalog number	Price per pack	
			1 pk	10 pk
Small	12 pairs	<a href="#">GH-86382-00</a>		
Medium		<a href="#">GH-86382-01</a>		
Large		<a href="#">GH-86382-02</a>		
X-large		<a href="#">GH-86382-03</a>		

## Nylon Inspection Gloves

### Won't leave lint behind

- Ideal for lint-sensitive areas
- Washables and reusable
- Ambidextrous

Economy gloves are constructed using a 40 denier tricot nylon in a two-piece pattern and feature rolled cuffs. Standard-weight gloves are constructed from a stretch nylon material in a two-piece pattern; cuffs are hemmed using a zig-zag pattern.



Size	Cuff	Qty/pk	Catalog number	Price per pack	
				1 pk	10 pk
<b>40 denier tricot nylon gloves</b>					
Small	Rolled	12 pairs	<a href="#">GH-86369-09</a>		
Large			<a href="#">GH-86369-13</a>		
<b>Standard-weight stretch nylon gloves</b>					
Small	Hemmed	12 pairs	<a href="#">GH-86369-10</a>		
Large			<a href="#">GH-86369-12</a>		

**Wait!**

**There's More at  
[ColeParmer.com](http://ColeParmer.com)**

- ✓ Search 100,000+ products
  - ✓ Find full product details
- And much more!

## Ansell HyFlex® Lite Precision Work Gloves

### Easily handle small parts

- Stretch nylon liner for better tactile sensitivity



Size	Catalog number	Price per pair	
		1 pr	12 pr
7	<a href="#">GH-09719-00</a>		
8	<a href="#">GH-09719-01</a>		
9	<a href="#">GH-09719-02</a>		
10	<a href="#">GH-09719-03</a>		

## Biohazard Autoclave Gloves

### Don't get burned

- Heat resistant to 450°F (232°C)—ideal for removing hot objects from an oven
- Machine washable
- One size fits all



Length	Catalog number	Price per pair	
		1 pr	12 pr
13" (33 cm)	<a href="#">GH-06496-40</a>		
18½" (47 cm)	<a href="#">GH-06496-42</a>		