

Scotsman® TouchFree® Flaked-Ice Maker/Dispenser

**Eliminate
cross-
contamination**

- Touch-free actuation provides convenient one-handed operation
- Up to 25% less fan noise
- Ice only or ice/water selectable
- Easy access for most containers



Condenser unit	Ice/24 hr	Bin capacity	Power (VAC, Hz)	Catalog number	Price
Air cooled	392 lb (178 kg)	12 lb (5.4 kg)	115, 60	GH-55100-06	

Wait!

There's More at
ColeParmer.com

- ✓ Search **100,000+** products
- ✓ Find full **product details**

And much more!

Frigimat® Junior Dry Ice Maker

No more wasted dry ice

- Produces only what you need in brick or snow form
- Simply attach to compressed LCO₂ tank, fill bag, and squeeze if hard pack desired
- Makes ~0.77 kg of flakes per bag

Note: LCO₂ Cylinder must have a dip tube.

What's included: collection bag, clamp, and plastic cone.



Description	Catalog number	Price
Junior dry ice maker	GH-03716-32	

[GH-03716-90](#) Extra bags. Pack of 3

Scotsman® Cubed- and Flaked-Ice Machines

**Don't run out
of ice**

- High-volume production up to 395 lbs per day
- Convenient size fits under standard lab benches
- Get easy access through slide-away door

Cubed-Ice Makers are constructed of durable rust-free rotocast plastic for greater insulation and reduced melting of ice.

Flaked-Ice Makers produce up to 395 lbs of flaked ice per day and store 140 lbs in the self-contained bin. Flaked ice quickly conforms to your container or sample tubes for greater cooling capabilities.



Cubed ice



Flaked ice



Condenser unit	Ice/24 hr	Bin capacity	Power (VAC, Hz)	Catalog number	Price
Cubed-ice makers					
Air cooled	84 lb (38 kg)	24 lb (11 kg)	115, 60	GH-55100-05	
Air cooled	150 lb (68 kg)	86 lb (39 kg)	115, 60	GH-55100-01	
Water cooled				GH-55100-02	
Flaked-ice makers					
Air cooled	395 lb (179 kg)	80 lb (36 kg)	115, 60	GH-55100-03	
Water cooled				GH-55100-04	

ThermoSafe Dry Ice Machine

**Get more dry ice out
of your CO₂**

- Provides 30% more yield from standard 50-lb CO₂ cylinder
- No need to dismantle to get dry ice

- No tip, stainless steel design

- Makes compact 1-lb blocks that measure 6"L x 4"W x 1 1/2"H

What's included: CO₂ supply regulator valve, hose, and fitted brass safety release valve that automatically vents excess gas pressure.



Description	Catalog number	Price
Dry ice machine	GH-03716-03	
Dry ice machine for European cylinders	GH-03716-04	