




Poly-Net™ Protective Netting, HDPE

Minimize impact on glassware

- Simply place over glassware to protect from bumps or flying debris
- Cut to fit your needs



Close mesh Heavy-gauge
(Shown at actual size)

Approximate diameter	Color	Thickness	Maximum temperature	Roll length	Catalog number	Price/roll
Close-mesh netting						
3/8" to 5/8" (1.0 to 1.6 cm)	Orange	0.050"	175°F (80°C)	820 ft (250 m)	GH-09405-15	
5/8" to 1" (1.6 to 2.5 cm)	Blue	0.050"				
1" to 2" (2.5 to 5.1 cm)	Yellow	0.050"	175°F (80°C)	164 ft (50 m)	GH-09405-25	
2" to 4" (5.1 to 10.2 cm)	Red	0.070"				
4" to 6" (10.2 to 15.2 cm)	Green	0.070"	175°F (80°C)	164 ft (50 m)	GH-09405-30	
6" to 8" (15.2 to 20.3 cm)	Brown	0.070"				
8" to 10" (20.3 to 25.4 cm)	Crystal	0.070"			GH-09405-40	
					GH-09405-45	
Heavy-gauge netting—more rigid for extra protection						
1 1/4" to 1 5/8" (3.2 to 4.1 cm)	Yellow	0.110"	175°F (80°C)	82 ft (25 m)	GH-09405-50	
2 1/2" to 3" (6.4 to 7.6 cm)	Green	0.110"				

Chemware® Lab Matting, PTFE Fluoropolymer

Don't let items go down the drain

- Line sinks and shelves to minimize impact
- Easily cut to size
- 1/4" (6 mm) mesh opening
- Autoclavable




(Shown at actual size)

Material	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
White	0.0625"	500°F (260°C)	12" x 10 ft (31 cm x 3 m)	GH-09406-00	

Cole-Parmer® PTFE Film

Stain and stick resistance

- Protect surfaces from chemicals and high temperatures
- Simply wipe away spills
- Cut to size with scissors for custom fit
- Reusable




Color	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
Opaque	0.005" (1.3 mm)	500°F (260°C)	12" x 25 ft (30.5 cm x 7.6 m)	GH-08277-15	

Ribbed Matting, LDPE

Minimize broken glass

- Prevent glass breakage by reducing impact on surfaces
- Solvent resistant



Color	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
Clear	0.036"	176°F (80°C)	30" x 15 ft (76.2 cm x 4.6 m)	GH-06809-10	
			30" x 50 ft (76.2 cm x 15 m)	GH-06809-20	

Questions?

Trust our **expert** Application Specialists to help you find the right product.

Bob, 7 years of experience

Call: 800-323-4340 E-mail: sales@coleparmer.com Chat live: ColeParmer.com