

Cole-Parmer® Mortar and Pestle Sets, Porcelain

No more messy transfers

- Pouring spout for easier pouring
- Complete set includes mortar and pestle



Buy More...Save More!
 20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Mortar			Pestle	Catalog number	Qty/ pk	Price/ pk
Capacity (mL)	Top OD (mm)	Height (mm)	Length (mm)			
80	60	50	64	GH-63100-54	1	
125	75	50	102	GH-63100-56	1	
300	100	90	152	GH-63100-58	1	
500	125	90	165	GH-63100-60	1	

Cole-Parmer® Mortar and Pestle Sets, Agate

No more cross-contamination

- Get optimal results and easier cleaning with highly polished grinding surface
- Naturally hard agate stone—6.5 hardness on Mohs scale
- Complete set (mortar with pestle) is autoclavable

NEW



Buy More...Save More!
 15% off on 10+ packs or cases 10% off on 5 to 9 packs or cases

Mortar			Pestle	Catalog number	Qty/ pk	Price/ pk
Capacity (mL)	Top OD (mm)	Height (mm)	Length (mm)			
3.8	38	25	38	GH-63101-30	1	
8	51	28	50	GH-63101-32		
40	76	34	76	GH-63101-34		
50	89	38	88	GH-63101-36	1	
90	102	41	101	GH-63101-38		
300	152	65	152	GH-63101-40		

CoorsTek® Mortars and Pestles, Porcelain

Get a better grip with oversized pestle

- Glazed for easier cleaning, except for outside bottom of mortar and grinding surface
- Order mortar and pestle separately
- Autoclavable

Buy More...Save More!
 15% off on 10+ packs or cases 10% off on 5 to 9 packs or cases

Piece	Capacity (mL)	Top OD (mm)	Height (mm)	Length (mm)	Catalog number	Qty/ pk	Price/ pk	Catalog number	Qty/ cs	Price/ cs
Mortar	65	80	53	—	GH-17820-02	1		GH-17820-22	18	
Pestle	—	—	—	130	GH-17821-02	1		GH-17821-22	18	
Mortar	145	90	70	—	GH-17820-04	1		GH-17820-24	12	
Pestle	—	—	—	157	GH-17821-04	1		GH-17821-24	12	
Mortar	275	115	70	—	GH-17820-06	1		GH-17820-26	8	
Pestle	—	—	—	180	GH-17821-06	1		GH-17821-26	8	
Mortar	400	130	80	—	GH-17820-08	1		GH-17820-28	8	
Pestle	—	—	—	194	GH-17821-08	1		GH-17821-28	8	
Mortar	750	163	110	—	GH-17820-10	1		GH-17820-30	4	
Pestle	—	—	—	222	GH-17821-10	1		GH-17821-30	4	
Mortar	1900	210	130	—	GH-17820-12	1		GH-17820-32	2	
Pestle	—	—	—	262	GH-17821-12	1		GH-17821-32	2	
Mortar	4000	283	150	—	—	—		GH-17820-34	1	
Pestle	—	—	—	280	—	—		GH-17821-34	1	



CoorsTek® Mortars and Pestles, Alumina

Nonabsorbent for minimal sample contamination

- Glazed for easier cleaning, except grinding surface and outside bottom of mortar
- Order mortar and pestle separately
- Autoclavable

Buy More...Save More!
 15% off on 10+ packs or cases 10% off on 5 to 9 packs or cases

Piece	Capacity (mL)	Top OD (mm)	Height (mm)	Length (mm)	Catalog number	Qty/ pk	Price/ pk	Catalog number	Qty/ cs	Price/ cs
Mortar	7	50	26	—	GH-17822-02	1		GH-17822-22	3	
Pestle	—	—	—	57	GH-17823-02	1		GH-17823-22	3	
Mortar	80	100	47	—	—	—		GH-17822-28	1	
Pestle	—	—	—	100	—	—		GH-17823-28	1	
Mortar	120	120	57	—	—	—		GH-17822-30	1	
Pestle	—	—	—	117	—	—		GH-17823-30	1	

