

Cole-Parmer® Triple Shelf

Get more storage without giving up bench space

- Minimal footprint takes up vertical space instead
- Don't worry about items falling off—safety lip keeps items on the shelf
- Won't scratch delicate glassware



Description	Dimensions (W x H x D)	Catalog number	Price
Triple shelf	24" x 16" x 6" (61 x 41 x 15 cm)	GH-63525-30	

Cole-Parmer® Shelves

Contain spills while keeping chemicals close at hand

- Maximize working space by raising items 8" off benchtop
- Maintain airflow with raised PVC design even in fume hoods



Description	Dimensions (W x H x D)	Catalog number	Price
3-foot shelf	5" x 8" x 36" (13 x 20 x 91 cm)	GH-63527-41	
4-foot shelf	5" x 8" x 48" (13 x 20 x 122 cm)	GH-63527-43	
5-foot shelf	5" x 8" x 60" (13 x 20 x 152 cm)	GH-63527-45	

Cole-Parmer® Open Supply Bins

Keep items close at hand and accessible

- Easy-to-clean PVC construction resists chemicals



Key	Description	Bin/shelf dimensions	Overall dimensions (W x H x D)	Catalog number	Price
A	Step shelf; place over meters	—	18" x 10" x 6" (46 x 25 x 15 cm)	GH-63529-55	
B	4-bin supply shelf	Bottom: 5¾"W; top: 8"D	24" x 12" x 10½" (61 x 31 x 27 cm)	GH-63529-60	
C	5-bin supply shelf	Bottom: 5½"W; middle: 3¾"W; top: 6"D	12" x 16" x 10½" (31 x 41 x 27 cm)	GH-63529-70	
D	8-bin supply shelf with cover	Bottom: 5¾"W; top: 5¾"W	24" x 12" x 10½" (61 x 31 x 27 cm)	GH-63529-65	
E	10-bin supply shelf	Bottom: 5¾"W; middle: 3¾"W; top: 6"D	24" x 16" x 10¾" (61 x 41 x 27 cm)	GH-63529-80	
F	13-bin supply shelf with cover	Bottom: 3¾"W; middle: 2¾"W; top: 1¾"W	12" x 12" x 10½" (31 x 31 x 27 cm)	GH-63529-85	