

## Cole-Parmer® Serological Pipettes, PS

### Fits just about any pipettor

These serological pipettes are made from transparent medical-grade polystyrene (PS), sterilized by gamma irradiation, and designed for quick and accurate fluid dispensing. The highly legible permanent graduations prevent ambiguity and guarantee precise dispensing to  $\pm 2\%$ . Graduations show volume removed or remaining on all pipettes larger than 1 mL. Generous negative graduations provide additional pipetting volume.

Unique filter plug eliminates lint fibers and ensures plug stays in place to prevent overfilling. The uniform mouth piece fits all major pipettors and has color-coded TD rings for easy identification. Each lot is tested for endotoxins and cytotoxicity, and guaranteed nonhemolytic—ideal for your tissue culture applications.

The individually wrapped pipettes are supplied in a fiber-free paper/plastic combination that is easy to open. The dispenser box allows you to remove pipettes without opening the top. Universal stacking design requires less laboratory space and each box is color coded for quick recognition.



**STERILE**

Capacity (mL)	Graduation interval (mL)	Tip style	Color ring	Qty/bag	Qty/cs	Catalog number	Price/cs
<b>Individually wrapped</b>							
1	0.01	Standard	Yellow	—	500	<a href="#">GH-13002-12</a>	
1	0.1	Standard	Red	—	500	<a href="#">GH-13002-14</a>	
2	0.01	Standard	White	—	500	<a href="#">GH-13002-15</a>	
5	0.1	Standard	Blue	—	200	<a href="#">GH-13002-16</a>	
10	0.1	Standard	Orange	—	200	<a href="#">GH-13002-17</a>	
25	0.2	Standard	Green	—	150	<a href="#">GH-13002-19</a>	
50	0.5	Standard	Purple	—	100	<a href="#">GH-13002-20</a>	
<b>Bulk wrapped</b>							
1	0.01	Standard	Yellow	25	1000	<a href="#">GH-13002-24</a>	
1	0.1	Standard	Red	25	1000	<a href="#">GH-13002-26</a>	
2	0.01	Standard	White	25	1000	<a href="#">GH-13002-28</a>	
5	0.1	Standard	Blue	25	500	<a href="#">GH-13002-29</a>	
10	0.1	Standard	Orange	25	400	<a href="#">GH-13002-30</a>	
10	0.1	Wide	Orange	25	400	<a href="#">GH-13002-32</a>	
25	0.2	Standard	Green	10	150	<a href="#">GH-13002-36</a>	
50	0.5	Standard	Purple	10	100	<a href="#">GH-13002-37</a>	

## Pipette Pump™ Controllers

### Dispense more quickly

- Attach pipette and easily dispense via plunger or wheel
- Made of HDPE



For pipettor volume	Color	Catalog number	Price
<b>Pipette pumps</b>			
0.2 mL	Yellow	<a href="#">GH-06221-01</a>	
2.0 mL	Blue	<a href="#">GH-06221-02</a>	
10.0 mL	Green	<a href="#">GH-06221-03</a>	
25.0 mL	Red	<a href="#">GH-06221-04</a>	
<b>Pipette pumps with rapid-release lever</b>			
2.0 mL	Blue	<a href="#">GH-06221-32</a>	
10.0 mL	Green	<a href="#">GH-06221-33</a>	
25.0 mL	Red	<a href="#">GH-06221-34</a>	

## Cole-Parmer® Safetypette Pipette Filler

### Get fast, comfortable suction

- Ergonomic, lightweight filler fits in either hand
- Draw solution in with bulb, control fluid precisely with button
- Use with any 1- to 100-mL pipette
- Hydrophobic filter prevents liquid from entering filler
- Autoclavable



## Silicone Pipette Filler

### Stop mouth pipetting

- Better chemical resistance than rubber

Description	Catalog number	Price
Silicone pipette filler, 50-mL capacity	<a href="#">GH-24805-20</a>	



**Cole-Parmer®**

Description	Catalog number	Price
Safetypette pipette filler	<a href="#">GH-25200-30</a>	

[GH-25200-32](#) Replacement filters. Pack of 5  
[GH-25200-34](#) Replacement pipette holder