

Scienceware® Sterile Sampling Scoops, PS

No need to presterilize

- Come sterilized in individual polyethylene bags
- Ideal for single use




STERILE

Capacity mL (oz)	Scoop (L x W)		Handle length		Catalog number	Qty/ cs	Price/ cs
	in.	cm	in.	cm			
60 (2)	3 x 1	7.6 x 2.5	4	10.2	GH-06260-10	100	
120 (4)	3½ x 1	8.9 x 2.5			GH-06260-20		

Scienceware® Sterile Sampling Spatula, PS

Get more volume per scoop

- Designed to hold twice the volume of standard spatulas
- Break up caked material easily with sturdy edge
- Individually sterilized and bagged



STERILE

Overall length	Catalog number	Qty/cs	Price/cs
9" (22.9 cm)	GH-06262-00	100	

Cole-Parmer® Double-Ended Spoon Set, 304 SS

Sample from both ends

- Minimize discomfort—solid rounded shafts for more comfortable handling

Set contains three spoons: a 5⁷/₈" (14.9 cm) long spoon with 0.5 mL and 2.0 mL bowls, an 8¹/₄" (21 cm) long spoon with 1.6 mL and 5.0 mL bowls, and an 11⁷/₈" (30.2 cm) long spoon with 2.0 mL and 8.0 mL bowls.

Cole-Parmer®

Buy More...Save More!

\$ 20% off on 120+
15% off on 30 to 119
10% off on 6 to 29 packs or cases

Set
contains
three
spoon
sizes



Qty/set	Catalog number	Price/set*
Three spoons	GH-06261-00	

*One set qualifies as one pack for discount savings.

Scienceware® Autoclavable Scoops, PP

Eliminate contamination risks

- Autoclavable at 250°F (121°C) to maintain sample integrity
- Rigid, strong design stands up to your needs



Buy More...Save More!

\$ 20% off on 120+
15% off on 30 to 119
10% off on 6 to 29 packs or cases

Capacity mL (oz)	Scoop (L x W)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
22 (¾)	2 ³ / ₈ x 1 ¹ / ₈	6.0 x 2.9	GH-06275-15	12	
60 (2)	3 x 1 ¹ / ₂	7.6 x 3.8	GH-06275-16	12	
75 (2½)	3 ¹ / ₂ x 1 ¹ / ₂	8.9 x 3.8	GH-06275-17	12	
100 (3½)	3 ⁵ / ₈ x 1 ⁷ / ₈	9.2 x 4.8	GH-06275-18	12	
200 (6¾)	5 x 2	12.7 x 5.1	GH-06275-19	12	
1100 (37)	8 x 4¾	20.3 x 12.1	GH-06275-20	6	

Scienceware® Sterile Sampling Spoons, PS

Stop washing spoons

- Individually sterilized and bagged in a polyethylene bag
- Ideal for single use



STERILE

Capacity	Handle length	Catalog number	Qty/cs	Price/cs
¼ tsp (1.2 mL)	6 ¹ / ₈ " (15.5 cm)	GH-06263-10	200	
½ tsp (2.5 mL)	5 ⁷ / ₈ " (14.9 cm)	GH-06263-20	200	
1 tsp (4.9 mL)	5 ¹ / ₂ " (14.0 cm)	GH-06263-30	200	
1 tbsp (14.8 mL)	5 ¹ / ₄ " (13.3 cm)	GH-06263-40	200	


Scienceware® Lab Spoon/Spatula, 304 SS

Easier sampling from large containers

- Total length of 12" (31 cm)

Buy More...Save More!

\$ 20% off on 120+
15% off on 30 to 119
10% off on 6 to 29 packs or cases



Spoon capacity	Spoon end	Spatula end	Catalog number	Qty/pk	Price/pk
5 mL (0.18 oz)	1 ¹ / ₂ " (3.8 cm) W	1 ³ / ₁₆ " (3 cm) W	GH-01019-14	1	