

Cole-Parmer® Griffin Low-Form Beakers, SS (medium gauge)

Get durable, economical mixing capabilities

- Sturdy vessels for general stirring or measuring of solutions
- Provide good, general chemical resistance



Buy More...Save More!
20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Capacity, mL (oz)	Catalog number	Qty/pk	Price/pk
125 (4.2)	GH-07205-30	1	
250 (8.5)	GH-07205-40	1	
600 (20.3)	GH-07205-50	1	
1200 (40.6)	GH-07205-60	1	
2000 (67.6)	GH-07205-70	1	
3000 (101.4)	GH-07205-80	1	
4000 (135.3)	GH-07205-90	1	

Cole-Parmer® Scoops, SS

Won't corrode like other materials—ideal for food-grade sampling

- Choose either 201 or 304 stainless steel—the 304 SS provides higher corrosion resistance



Buy More...Save More!
20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Capacity mL (oz)	Scoop (L x dia)	Catalog number	Qty/pk	Price/pk
201 Stainless steel scoops				
175 (5.5)	5½" x 3" (14 x 7.6 cm)	GH-07205-01	1	
1550 (52)	10" x 5¼" (25.4 x 14.6 cm)	GH-07205-03†	1	
304 Stainless steel scoop				
1550 (52)	7⅞" x 5¼" (18.1 x 13.3 cm)	GH-07210-00	1	

†Model 07205-03 has a hook on the end for convenient storage.

Cole-Parmer® Deep Trays and Lids, SS (medium gauge)

Won't fail under the heat

- Handles high temperatures up to 650°F (182°C)
- Protect contents with optional flat cover


Buy More...Save More!
20% off on 120+ 15% off on 30 to 119 10% off on 6 to 29 packs or cases

Capacity L (gal.)	Dimensions (L x W x H)		Tray		Optional Cover	
	in.	cm	Catalog number	Price/pk of 1	Catalog number	Price/pk of 1
1 (¼)	6⅞ x 6⅞ x 2½	17.5 x 16.2 x 6.4	GH-07242-10		GH-07242-60	
2 (½)	10⅞ x 6⅞ x 2½	26.4 x 16.2 x 6.4	GH-07243-10		GH-07234-60	
2.5 (⅝)	12¾ x 6⅞ x 2½	32.4 x 17.5 x 6.4	GH-07244-10		GH-07236-60	
3 (¾)	10⅞ x 6⅞ x 4	26.4 x 16.2 x 10.2	GH-07245-10		GH-07234-60	
4 (1)	12¾ x 6⅞ x 4	32.4 x 17.5 x 10.2	GH-07246-10		GH-07236-60	
4.25 (1⅙)	10⅞ x 6⅞ x 6	26.4 x 16.2 x 15.2	GH-07234-10		GH-07234-60	
4.5 (1¼)	16½ x 10 x 2½	41.9 x 25.4 x 6.4	GH-07235-05		GH-07235-60	
6 (1½)	12¾ x 6⅞ x 6	32.4 x 17.5 x 15.2	GH-07236-10		GH-07236-60	
7 (1¾)	12¾ x 10⅞ x 4	32.4 x 26.4 x 10.2	GH-07235-10		GH-07238-60	
9 (2¼)	20¾ x 12¾ x 2½	52.7 x 32.4 x 6.4	GH-07237-10		GH-07240-60	
10 (2½)	12¾ x 10⅞ x 6	32.4 x 26.4 x 15.2	GH-07238-10		GH-07238-60	
10 (2½)	12¾ x 10⅞ x 8	32.4 x 27.3 x 20.3	GH-07239-05†		GH-07239-60	
15 (3¾)	20¾ x 12¾ x 4	52.7 x 32.4 x 10.2	GH-07239-10		GH-07240-60	
22 (5½)	20¾ x 12¾ x 6	52.7 x 32.4 x 15.2	GH-07240-10		GH-07240-60	
29 (7¾)	20¾ x 12¾ x 8	52.7 x 32.4 x 20.3	GH-07241-00		GH-07240-60	
30 (8)	21¾ x 13¾ x 8	55.2 x 34.9 x 20.3	GH-07240-95†		GH-07241-60	



Find even more online! Go to [ColeParmer.com](#) to search more than 100,000 products.

Wait! There's More...

-  Beakers
Pages 64–69.
-  Bottles
Pages 72–90.
-  Bowls & Containers
Pages 178–183.
-  Pails & Buckets
Page 484.
-  Scoops & Spoons
Pages 733–734.
-  Trays
Pages 861–863.