

## Beaker and Bottle Brushes

No more dirty glassware—clean it easily

- Choose the right brush for your glassware from shape to bristle type
- Radial tip brushes reduce scratches and provide effective cleaning
- Contoured brushes allow you to reach angled areas with ease
- Tufted end brushes fold to get in narrow-mouth vessels and then open once inside
- Choose angle and swivel head brushes to easily clean curves
- Volumetric models have a bendable handle to fit container's shape

**Wait!**

There's More at  
[ColeParmer.com](http://ColeParmer.com)

Find more than  
**100,000** products  
to choose from!



Key	Description	Brush bristles		Handle length		Catalog number	Qty/pk	Price
		Size (L x W)	Material	in.	cm			
<b>Beaker/bottle/flask brushes</b>								
A	Radial tip, small	3" x 1 3/4"	7.6 x 3.2 cm	Polyester	9	22.9	<a href="#">GY-84550-80</a>	12
	Radial tip, large	6" x 3 1/8"	15.2 x 7.9 cm		10	25.4	<a href="#">GY-84550-81</a>	12
	Radial tip, XL	6" x 3 1/8"	15.2 x 7.9 cm		13	33	<a href="#">GY-84550-82</a>	6
B	Contoured, large	10" x 3 1/8"	25.4 x 7.9 cm		7	17.8	<a href="#">GY-84550-83</a>	6
C	Fan tip, small	6" x 1 1/2"	15.2 x 3.8 cm	Natural	6	15.2	<a href="#">GY-84550-86</a>	12
D	Tufted end, small	4 1/4" x 1 1/2"	11.4 x 3.8 cm		7 1/2	19.1	<a href="#">GY-84550-87</a>	12
	Tufted end, XL	4 1/4" x 1 3/4"	11.4 x 4.5 cm		19 1/2	49.5	<a href="#">GY-84550-88</a>	12
E	Angled head, small	3" x 1"	7.6 x 2.5 cm	Polyester	13	33	<a href="#">GY-84550-76</a>	1
	Angled head, large	4" x 1 1/2"	0.2 x 3.8 cm	Natural	12	30.5	<a href="#">GY-84550-89</a>	12
	Volumetric, large	4" x 1 1/2"	10.2 x 3.8 cm	Nylon	10	25.4	<a href="#">GY-84550-90</a>	12
<b>Bottle/cylinder brushes</b>								
F	Radial tip, white bristles	6 1/2" x 3 3/8" dia	16.5 x 7.9 cm dia	Nylon	9	22.9	<a href="#">GY-84550-00</a>	1
	Radial tip, X-shaped yellow bristles						<a href="#">GY-84550-03</a>	1
	Radial tip, X-shaped green bristles						<a href="#">GY-84550-04</a>	1
G	Radial tip, white bristles			Polyester	9	22.9	<a href="#">GY-84550-70</a>	1
	Radial tip, red bristles						<a href="#">GY-84550-72</a>	1
	Radial tip, blue bristles						<a href="#">GY-84550-74</a>	1
H	Radial tip, small	4 1/2" x 1 3/4"	11.4 x 4.5 cm	Natural	7	17.8	<a href="#">GY-84550-15</a>	12
	Radial tip, medium	5" x 2"	12.7 x 5.1 cm		11	27.9	<a href="#">GY-84550-16</a>	12
	Radial tip, large	6" x 2 3/8"	15.2 x 6.1 cm	Nylon	14	35.6	<a href="#">GY-84550-17</a>	12
	Radial tip, XL	6" x 2 3/8"	15.2 x 6.1 cm		18	45.7	<a href="#">GY-84550-18</a>	12
	Radial tip, XXL	6" x 2 3/4"	15.2 x 7 cm		24	60.9	<a href="#">GY-84550-19</a>	12
<b>Glassware brushes</b>								
I	Foam-tipped	7" x 2"	17.8 x 5.1 cm	Polypropylene/foam	13	33	<a href="#">GY-84552-14</a>	1
		5" x 2"	12.7 x 5.1 cm		7	17.8	<a href="#">GY-84552-16</a>	1
J	Cylinder/Funnel brush	16" x 3 1/8"	40.6 x 7.9 cm	Polyester	10 1/2	26.7	<a href="#">GY-84550-71</a>	12
		24" x 3 1/8"	60.9 x 7.9 cm		17 1/2	44.5	<a href="#">GY-84550-73</a>	12
K	Swivel head brush	1 3/4" dia	3.2 cm	Natural	13	33	<a href="#">GY-84550-55</a>	1