


Disposable 50-mm Sterile Petri Dishes, PS

No need to presterilize

- Sterilized dishes are ready to use
- Use only what you need—come in convenient packs of 20
- Polystyrene dishes come with or without 47-mm absorbent cellulose pad that hold 1.8 to 2.0 mL of liquid media
- Ridges for easy stacking



Description	Dia x H (mm)	Qty/cs	Catalog number	Price/cs
Petri dish with pad	50 x 11	100	GY-14005-20	
		500	GY-14005-28	
Petri dish without pad	50 x 11	100	GY-14005-24	
		500	GY-14005-26	

[GY-14005-22](#) Presterilized absorbent pads. Box of 100 individual packages

Petri Dishes, PFA

Maintain purity of your culture

- Ideal for contamination-free biological applications
- Inert, nonstick, noncontaminating surface
- Lids not included
- Autoclavable



Capacity (mL)	Dia x H (mm)	Qty/pk	Catalog number	Price/pk
60	50 x 15	1	GY-06148-00	
100	100 x 15	1	GY-06148-02	

Disposable Gridded Sterile Contact Plates, PS

Identify and monitor cultures easily

- For determining surface microbial contamination
- Dishes have a raised agar bed for surface contact
- Stackable lids are flared for easy handling and automatic filling



Dia x H (mm)	Catalog number	Qty/cs	Price/cs
65 x 15	GY-14006-50	500	

The grid consists of sixteen 10 x 10 mm squares. Disposable polystyrene (PS). Radiation sterilized.

Petri Dishes, Glass

Save money—autoclavable for repeated use

- Dishes withstand repeated wet or dry sterilization
- Easy sorting—tops marked in blue, bottoms marked in white



Dia x H (mm)	Catalog number	Qty/pk	Price/pk	Catalog number	Qty/cs	Price/cs
60 x 15	GY-34551-01	12		GY-34551-21	72	
100 x 10	GY-34551-02	12		GY-34551-22	72	
100 x 15	GY-34551-03	12		GY-34551-23	72	
100 x 20	GY-34551-04	12		GY-34551-24	72	
150 x 15	GY-34551-05	12		GY-34551-25	24	
150 x 20	GY-34551-06	12		GY-34551-26	24	

Sterile Petri Dishes, PS

Stacks won't tip over—convenient stacking ridges for better stability

- Use what you need—packaged in bags of 20 dishes
- Save time—already sterilized by gamma irradiation



Dia x H (mm)	Qty/pk	Catalog number	Price/pk
60 x 15	500	GY-06139-01	
100 x 15	500	GY-06139-00	
150 x 20	100	GY-06139-02	

Disposable Cell Spreaders

Easy spreading with L-shaped design

- Won't cut into agar—upturned end to protect from puncture

Made of durable polystyrene, these sterile, disposable spreaders provide safe, convenient sample plating without flaming.



Description	Qty/pk	Qty/cs	Catalog number	Price/cs
Disposable cell spreaders	10	500	GY-14201-00	
	1	500	GY-14201-10	