

Single barb fittings seal better than multi-barb fittings because they allow the tubing to relax around the barb.




FREE
certificate of compliance
See page 313.


Tubing ID	Material	Key	Catalog number	Qty/ pk	Price/ pk
Y connectors—continued					
14 mm	Clear PP	AM	GY-06295-60	20	
14–16 mm	White PP	AN	GY-06296-60	20	
5/8"	HDPE	AJ	GY-30726-08	10	
5/8"	Nylon	AK	GY-30726-28	10	
5/8"	PP	AL	GY-30726-48	10	
5/8"	PVDF	AO	GY-30726-68	10	
16 mm	Clear PP	AM	GY-06295-70	20	
Cross connectors					
1/16"	316 SS	AR	GY-31208-51	1	
3/32"	316 SS	AR	GY-31208-52	1	
1/8"	316 SS	AR	GY-31208-53	1	
4 mm	Clear PP	AQ	GY-06297-10	20	
3/16"	316 SS	AR	GY-31208-54	1	
1/4"	Black PP	AS	GY-31407-11	5	
1/4"	Nylon	AT	GY-31408-11	5	
6 mm	Clear PP	AQ	GY-06297-20	20	
8 mm	Clear PP	AQ	GY-06297-30	20	
3/8"	316 SS	AR	GY-31208-55	1	
10 mm	Clear PP	AQ	GY-06297-40	20	
12 mm	Clear PP	AQ	GY-06297-50	20	
14 mm	Clear PP	AQ	GY-06297-60	20	
NPT male threaded adapters, NPT(M) x tubing ID					
1/16" x 1/16"	PVDF	AZ	GY-30704-00	10	
1/16" x 1/16"	316 SS	BC	GY-31208-59	1	
1/16" x 3/32"	316 SS	BC	GY-31208-60	1	
1/16" x 1/8"	PVDF	AZ	GY-30704-02	10	
1/16" x 1/8"	316 SS	BC	GY-31208-61	1	
1/16" x 3/16"	PVDF	AZ	GY-30704-04	10	
1/16" x 3/16"	316 SS	BC	GY-31208-62	1	
1/16" x 1/4"	316 SS	BC	GY-31208-63	1	
1/8" x 1/16"	Clear PP	AY	GY-06365-41	25	
1/8" x 1/16"	PVDF	AZ	GY-30704-05	10	
1/8" x 1/16"	316 SS	BC	GY-31208-65	1	
1/8" x 3/32"	Clear PP	AY	GY-06365-42	25	
1/8" x 3/32"	PVDF	AZ	GY-30704-06	10	
1/8" x 3/32"	316 SS	BC	GY-31208-66	1	
1/8" x 1/8"	HDPE	AU	GY-30621-26	10	
1/8" x 1/8"	Nylon	AV	GY-30621-27	10	
1/8" x 1/8"	PP	AX	GY-40621-28	10	
1/8" x 1/8"	Clear PP	AY	GY-06365-43	25	
1/8" x 1/8"	PVDF	AZ	GY-30704-07	10	
1/8" x 1/8"	316 SS	BC	GY-30621-30	1	
1/8" x 5/32"	HDPE	AU	GY-31517-21	10	
1/8" x 5/32"	Nylon	AV	GY-30621-24	10	
1/8" x 5/32"	PP	AX	GY-40621-22	10	
1/8" x 5/32"	PVDF	AZ	GY-30704-08	10	
1/8" x 5/32"	PFA	BA	GY-06383-90	1	
1/8" x 3/16"	HDPE	AU	GY-31702-04	10	
1/8" x 3/16"	Nylon	AV	GY-31702-02	10	
1/8" x 3/16"	Black PP	AW	GY-06444-16	10	
1/8" x 3/16"	316 SS	BC	GY-31208-68	1	
1/8" x 1/4"	HDPE	AU	GY-30621-31	10	
1/8" x 1/4"	Nylon	AV	GY-30621-32	10	
1/8" x 1/4"	PP	AX	GY-40610-36	10	
1/8" x 1/4"	PVDF	AZ	GY-30704-09	10	
1/8" x 1/4"	PFA	BA	GY-06383-91	1	
1/8" x 1/4"	Brass	BB	GY-30904-01	5	
1/8" x 1/4"	316 SS	BC	GY-30621-33	1	

Tubing ID	Material	Key	Catalog number	Qty/ pk	Price/ pk
NPT male threaded adapters, NPT(M) x tubing ID—cont.					
1/8" x 5/16"	HDPE	AU	GY-30401-72	10	
1/8" x 5/16"	Nylon	AV	GY-30622-31	10	
1/8" x 5/16"	PP	AX	GY-41517-27	10	
1/8" x 5/16"	PVDF	AZ	GY-31517-26	10	
1/8" x 3/8"	HDPE	AU	GY-30621-35	10	
1/8" x 3/8"	Nylon	AV	GY-30621-36	10	
1/8" x 3/8"	PP	AX	GY-40610-37	10	
1/8" x 3/8"	PVDF	AZ	GY-30704-10	10	
1/8" x 3/8"	Brass	BB	GY-30904-00	5	
1/8" x 3/8"	316 SS	BC	GY-31208-70	1	
1/4" x 1/16"	316 SS	BC	GY-31208-74	1	
1/4" x 3/32"	316 SS	BC	GY-31208-75	1	
1/4" x 1/8"	HDPE	AU	GY-30621-42	10	
1/4" x 1/8"	Nylon	AV	GY-30621-44	10	
1/4" x 1/8"	PP	AX	GY-40621-45	10	
1/4" x 1/8"	PVDF	AZ	GY-30621-43	10	
1/4" x 1/8"	316 SS	BC	GY-31208-76	1	
1/4" x 5/32"	HDPE	AU	GY-30621-38	10	
1/4" x 5/32"	Nylon	AV	GY-30621-40	10	
1/4" x 5/32"	PP	AX	GY-40621-41	10	
1/4" x 5/32"	PVDF	AZ	GY-30621-39	10	
1/4" x 5/32"	316 SS	BC	GY-31208-77	1	
1/4" x 3/16"	HDPE	AU	GY-31702-08	10	
1/4" x 3/16"	Nylon	AV	GY-31702-06	10	
1/4" x 3/16"	PP	AX	GY-41517-47	10	
1/4" x 3/16"	Black PP	AW	GY-06444-18	10	
1/4" x 3/16"	PVDF	AZ	GY-31517-46	10	
1/4" x 3/16"	316 SS	BC	GY-31208-78	1	
1/4" x 1/4"	HDPE	AU	GY-30621-46	10	
1/4" x 1/4"	Nylon	AV	GY-30621-47	10	
1/4" x 1/4"	PP	AX	GY-40610-38	10	
1/4" x 1/4"	PVDF	AZ	GY-30704-11	10	
1/4" x 1/4"	PFA	BA	GY-06383-92	1	
1/4" x 1/4"	Brass	BB	GY-30904-10	5	
1/4" x 1/4"	316 SS	BC	GY-31208-79	1	
1/4" x 5/16"	HDPE	AU	GY-30621-48	10	
1/4" x 5/16"	Nylon	AV	GY-30621-50	10	
1/4" x 5/16"	PP	AX	GY-40621-51	10	
1/4" x 5/16"	PVDF	AZ	GY-30621-49	10	
1/4" x 3/8"	HDPE	AU	GY-30621-52	10	
1/4" x 3/8"	Nylon	AV	GY-30621-53	10	
1/4" x 3/8"	PP	AX	GY-40610-39	10	
1/4" x 3/8"	PVDF	AZ	GY-30704-12	10	
1/4" x 3/8"	PFA	BA	GY-06383-93	1	
1/4" x 3/8"	Brass	BB	GY-30904-11	5	
1/4" x 3/8"	316 SS	BC	GY-31208-80	1	
1/4" x 1/2"	HDPE	AU	GY-30621-54	10	
1/4" x 1/2"	Nylon	AV	GY-30621-55	10	
1/4" x 1/2"	PP	AX	GY-40610-40	10	
1/4" x 1/2"	PVDF	AZ	GY-30704-13	10	
3/8" x 3/32"	316 SS	BC	GY-31208-83	1	
3/8" x 1/8"	316 SS	BC	GY-31208-84	1	
3/8" x 1/4"	HDPE	AU	GY-30621-56	10	
3/8" x 1/4"	Nylon	AV	GY-30621-57	10	
3/8" x 1/4"	Clear PP	AY	GY-40610-41	10	
3/8" x 1/4"	PVDF	AZ	GY-30704-14	10	
3/8" x 1/4"	Brass	BB	GY-30904-20	5	
3/8" x 1/4"	316 SS	BC	GY-31208-85	1	
3/8" x 5/16"	HDPE	AU	GY-30621-58	10	
3/8" x 5/16"	Nylon	AV	GY-30621-60	10	
3/8" x 5/16"	PP	AX	GY-40621-61	10	
3/8" x 5/16"	PVDF	AZ	GY-30621-59	10	


Y connectors (continued)



Cross connectors



NPT male threaded adapters



See next page for more...