

Cole-Parmer® Serological Pipettes, PS

Fit just about any pipettor

These serological pipettes are made from transparent polystyrene (PS), sterilized by gamma irradiation, and designed for quick and accurate fluid dispensing. The highly legible permanent graduations prevent ambiguity and guarantee precise dispensing to $\pm 2\%$. Graduations show volume removed or remaining on all pipettes larger than 1 mL. Generous negative graduations provide additional pipetting volume.

Unique filter plug eliminates lint fibers and ensures plug stays in place to prevent overfilling. The uniform mouth piece fits all major pipettors and has color-coded TD rings for easy identification. Each lot is tested for endotoxins and cytotoxicity, and guaranteed nonhemolytic—ideal for your tissue culture applications.

The individually wrapped pipettes are supplied in a fiber-free paper/plastic combination that is easy to open. The dispenser box allows you to remove pipettes without opening the top. Universal stacking design requires less laboratory space and each box is color coded for quick recognition.



STERILE

Capacity (mL)	Graduation interval (mL)	Tip style	Color ring	Qty/bag	Qty/cs	Catalog number	Price/cs
Individually wrapped							
1	0.01	Standard	Yellow	—	500	GY-13002-12	
1	0.1	Standard	Red	—	500	GY-13002-14	
2	0.01	Standard	White	—	500	GY-13002-15	
5	0.1	Standard	Blue	—	200	GY-13002-16	
10	0.1	Standard	Orange	—	200	GY-13002-17	
25	0.2	Standard	Green	—	150	GY-13002-19	
50	0.5	Standard	Purple	—	100	GY-13002-20	
Bulk wrapped							
1	0.01	Standard	Yellow	25	1000	GY-13002-24	
1	0.1	Standard	Red	25	1000	GY-13002-26	
2	0.01	Standard	White	25	1000	GY-13002-28	
5	0.1	Standard	Blue	25	500	GY-13002-29	
10	0.1	Standard	Orange	25	400	GY-13002-30	
10	0.1	Wide	Orange	25	400	GY-13002-32	
25	0.2	Standard	Green	10	150	GY-13002-36	
50	0.5	Standard	Purple	10	100	GY-13002-37	

Pipette Helper™ Controller

Chemical resistant

- Single-lever operation for even flow
- Fits pipettes 0.1 to 100 mL
- PP, PTFE, and silicone wetted parts

ISO9001:2008
CERTIFIED SUPPLIER

Description	Catalog number	Price
Pipette Helper™ controller	GY-24549-00	

[GY-24549-50](#) Replacement membrane filter, 3 µm. Pack of 10

Pipette Controllers

Protect against contamination

- Color-coded for quick identification
- Single-lever operation
- Fits pipettes 0.1 to 100 mL
- Autoclavable

ISO9001:2008
CERTIFIED SUPPLIER

Colors	Catalog number	Price
Neutral	GY-24550-49	
Blue	GY-24550-51	
Magenta	GY-24550-53	
Green	GY-24550-57	



Cole-Parmer® Safetypette Pipette Filler

Get fast, comfortable suction

- Ergonomic, lightweight filler fits in either hand
- Draw solution in with bulb, control fluid precisely with button
- Use with any 1- to 100-mL pipette
- Hydrophobic filter prevents liquid from entering filler
- Autoclavable

Cole-Parmer®

Description	Catalog number	Price
Safetypette pipette filler	GY-25200-30	

[GY-25200-32](#) Replacement filters. Pack of 5

[GY-25200-34](#) Replacement pipette holder



Silicone Pipette Filler

Stop mouth pipetting

- Better chemical resistance than rubber

Description	Catalog number	Price
Silicone pipette filler, 50-mL capacity	GY-24805-20	

