



Rack 67035-02 shown with drain basket and glove dispenser

Rack 67035-10 shown with pipette rack and screen insert

## Modular Drying Racks, 304 SS

### Don't settle—customize to your needs in seconds

- Easy-to-clean, corrosion-resistant stainless steel racks and polypropylene (PP) pegs
- No mess—catch drain directs drainage to drip trough
- Autoclavable
- Include 6" L pegs, SS drain tube, and wall-mounting bracket



Color of pegs (included)	Catalog number	Price/ea
<b>32-peg racks</b> measure 24"W x 30"H x 4"D (61 x 76.2 x 10.2 cm)		
White	<a href="#">GY-67035-00</a>	
Assorted†	<a href="#">GY-67035-02</a>	
<b>50-peg racks</b> measure 30"W x 30"H x 4"D (76.2 x 76.2 x 10.2 cm)		
White	<a href="#">GY-67035-10</a>	
Assorted†	<a href="#">GY-67035-12</a>	
Black	<a href="#">GY-67035-14</a>	
<b>60-peg racks</b> measure 36"W x 30"H x 4"D (91.4 x 76.2 x 10.2 cm)		
White	<a href="#">GY-67035-20</a>	
Assorted†	<a href="#">GY-67035-22</a>	
Black	<a href="#">GY-67035-24</a>	

†Assorted peg colors include yellow, orange, pink, red, green, light blue, dark blue, lilac, and purple.

### Add-On Accessories, 304 SS

These accessories include the necessary "snap-in" support pegs to quickly and easily attach to the rack. **Drain basket** holds small items such as stir bars, stoppers, and pipette tips. **Drain shelf** is ideal for drying and storing various glassware items. **Pipette rack** hangs at an angle to allow for complete drainage of pipettes. **Funnel rack** holds four funnels, volumetric flasks, or small cylinders. **Cylinder yoke** holds two cylinders. **Flask holder** accepts one large flask, cylinder, or separatory funnel. **Glove dispenser** provides a convenient spot for one box of gloves. **Screen inserts** fit into the drip trough and increase the drying capacity of the rack.

Description	Catalog number	Price/ea
Drain basket, 6"W x 2½"H (15.2 x 6.4 cm)	<a href="#">GY-67035-70</a>	
Drain shelf, 14"W (35.6 cm)	<a href="#">GY-67035-71</a>	
Pipette rack, 12½"L (31.8 cm)	<a href="#">GY-67035-72</a>	
Funnel rack, 14"W (35.6 cm)	<a href="#">GY-67035-73</a>	
Cylinder yoke, 14"W (35.6 cm)	<a href="#">GY-67035-74</a>	
Flask holder, 8"W (20.3 cm)	<a href="#">GY-67035-75</a>	
Glove dispenser, 10"W x 5¼"H (25.4 x 14.6 cm)	<a href="#">GY-67035-76</a>	
Screen insert, 24"W (61 cm)	<a href="#">GY-67035-77</a>	
Screen insert, 30"W (76.2 cm)	<a href="#">GY-67035-78</a>	
Screen insert, 36"W (91.4 cm)	<a href="#">GY-67035-80</a>	

### Plugs and Pegs, Molded Polypropylene

Removable "snap-in" plugs and pegs come in packs of six.

Description	Catalog number	Price/pk
White hole plugs	<a href="#">GY-67035-79</a>	
White tube pegs, 4"L x 5/16" dia (10.2 x 0.8 cm)	<a href="#">GY-67035-90</a>	
White long pegs, 9"L x 5/8" dia (22.9 x 1.6 cm)	<a href="#">GY-67035-91</a>	
White long pegs, 12"L x 5/8" dia (30.5 x 1.6 cm)	<a href="#">GY-67035-92</a>	

