

Sterile Sampling Scoops, PS

No need to presterilize

- Come sterilized in individual polyethylene bags
- Ideal for single use



STERILE

Capacity mL (oz)	Scoop (L x W)		Handle length		Catalog number	Qty/ cs	Price/ cs
	in.	cm	in.	cm			
60 (2)	3 x 1	7.6 x 2.5	4	10.2	GY-06260-10	100	
120 (4)	3½ x 1	8.9 x 2.5			GY-06260-20		

Sterile Sampling Spatula, PS

Get more volume per scoop

- Designed to hold twice the volume of standard spatulas
- Break up caked material easily with sturdy edge
- Individually sterilized and bagged



STERILE

Overall length	Catalog number	Qty/cs	Price/cs
9" (22.9 cm)	GY-06262-00	100	

Cole-Parmer® Double-Ended Spoon Set, 304 SS

Sample from both ends

- Minimize discomfort—solid rounded shafts for more comfortable handling

Set contains three spoons: a 5 7/8" (14.9 cm) long spoon with 0.5 mL and 2.0 mL bowls, an 8 1/4" (21 cm) long spoon with 1.6 mL and 5.0 mL bowls, and an 11 7/8" (30.2 cm) long spoon with 2.0 mL and 8.0 mL bowls.

 **Cole-Parmer®**

Set contains three spoon sizes



Qty/set	Catalog number	Price/set
Three spoons	GY-06261-00	

Autoclavable Scoops, PP

Eliminate contamination risks

- Autoclavable at 250°F (121°C) to maintain sample integrity
- Rigid, strong design stands up to your needs



Capacity mL (oz)	Scoop (L x W)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
22 (¾)	2 3/8 x 1 1/8	6.0 x 2.9	GY-06275-15	12	
60 (2)	3 x 1 1/2	7.6 x 3.8	GY-06275-16	12	
75 (2½)	3 1/2 x 1 1/2	8.9 x 3.8	GY-06275-17	12	
100 (3½)	3 5/8 x 1 7/8	9.2 x 4.8	GY-06275-18	12	
200 (6¾)	5 x 2	12.7 x 5.1	GY-06275-19	12	
1100 (37)	8 x 4¾	20.3 x 12.1	GY-06275-20	6	

Sterile Sampling Spoons, PS

Stop washing spoons

- Individually sterilized and bagged in a polyethylene bag
- Ideal for single use



STERILE

Capacity	Handle length	Catalog number	Qty/cs	Price/cs
¼ tsp (1.2 mL)	6 1/8" (15.5 cm)	GY-06263-10	200	
½ tsp (2.5 mL)	5 7/8" (14.9 cm)	GY-06263-20	200	
1 tsp (4.9 mL)	5 1/2" (14.0 cm)	GY-06263-30	200	
1 tbs (14.8 mL)	5 1/4" (13.3 cm)	GY-06263-40	200	

Lab Spoon/Spatula, 304 SS

Easier sampling from large containers

- Total length of 12" (31 cm)



Spoon capacity	Spoon end	Spatula end	Catalog number	Qty/pk	Price/pk
5 mL (0.18 oz)	1 1/2" (3.8 cm) W	1 3/16" (3 cm) W	GY-01019-14	1	