

### Microcentrifuge Tubes

#### A. Cole-Parmer® Tubes

- Boil-proof tubes have attached lid and are guaranteed to stay closed
- Frosted writing area and assorted colors provide easy sample identification
- Certified RNase-, DNase-, and endotoxin-free

#### B. Cole-Parmer® Maximum Recovery Tube

- Same features as "A" plus free from cavities or occlusions causing this effect
- Ideal for long term storage of long chain molecules

#### C. Eppendorf® Safe-Lock™ Tubes

- Secure closure with included cap that opens easily with one hand
- Feature marking area and graduations



#### D. Eppendorf® Biopur® Tubes

- Same features as "A" and "C", plus certified sterile and wrapped individually to maintain sterility

#### E. Cole-Parmer® Screw-Cap Tubes

- Preserve sample purity with external threads
- Ideal for freezing, boiling, autoclaving, or irradiation
- Order caps with leakproof O-rings separately

#### F. BackSpring™ Tube

- Eliminate stress and sore thumbs with back-open flip cap—flips open with thumb pressure
- Boil-proof and freeze-proof to -90°C
- RNase-, DNase-, and pyrogen-free
- Autoclavable



Key	Capacity (mL)	Description	Tube material	Max RCF (x g)	Color	Cap included	Sterile	Catalog number	Qty/pk	Price/pk
A	0.6	Cole-Parmer	PP	20,000	Natural	Yes	No	<a href="#">GY-06333-40</a>	1000	
A	0.6	Cole-Parmer	PP	20,000	Natural	Yes	Yes	<a href="#">GY-06333-42</a>	500	
A	0.6	Cole-Parmer	PP	20,000	Amber	Yes	No	<a href="#">GY-06333-43</a>	1000	
A	1.5	Cole-Parmer	PP	20,000	Natural	Yes	No	<a href="#">GY-06333-50</a>	500	
A	1.5	Cole-Parmer	PP	20,000	Assorted	Yes	No	<a href="#">GY-06333-51</a>	500	
A	1.5	Cole-Parmer	PP	20,000	Amber	Yes	No	<a href="#">GY-06333-56</a>	500	
A	1.5	Cole-Parmer	PP	20,000	Assorted	Yes	Yes	<a href="#">GY-06333-55</a>	250	
B	1.5	Cole-Parmer max recovery	PP	20,000	Assorted	Yes	No	<a href="#">GY-06333-53</a>	250	
C	1.5	Eppendorf Safe-Lock	PP	25,000	Natural	Yes	No	<a href="#">GY-02560-33</a>	500	
D	1.5	Eppendorf Biopur	PP	25,000	Natural	Yes	Yes	<a href="#">GY-02560-35</a>	100	
E	1.5	Cole-Parmer, screw-cap†	PP	30,000	Natural	No	Yes	<a href="#">GY-06245-50</a>	1000	
A	1.7	Cole-Parmer	PP	20,000	Natural	Yes	No	<a href="#">GY-06333-60</a>	500	
A	2.0	Cole-Parmer	PP	20,000	Natural	Yes	No	<a href="#">GY-06333-70</a>	500	
A	2.0	Cole-Parmer	PP	20,000	Assorted	Yes	No	<a href="#">GY-06333-71</a>	500	
A	2.0	Cole-Parmer	PP	20,000	Natural	Yes	Yes	<a href="#">GY-06333-72</a>	250	
E	2.0	Cole-Parmer, screw-cap†	PP	30,000	Natural	No	No	<a href="#">GY-06245-51</a>	1000	
E	2.0	Cole-Parmer, screw-cap, skirted†‡	PP	30,000	Natural	No	No	<a href="#">GY-06245-52</a>	1000	
C	2.0	Eppendorf Safe-Lock	PP	25,000	Natural	Yes	No	<a href="#">GY-02560-36</a>	500	
D	2.0	Eppendorf Biopur	PP	25,000	Natural	Yes	Yes	<a href="#">GY-02560-38</a>	100	
F	1.5	BackSpring	Homopolymer Plastic	30,000	Clear	Yes	No	<a href="#">GY-22999-00</a>	500	

†Models 06245-50, -51, and -52 require screw caps (sold separately below).

‡Model 06245-52 requires microcentrifuge adapter (sold separately below).

[GY-06245-55](#) Caps for screw-cap tubes, blue. Pack of 1000

[GY-06245-56](#) Caps for screw-cap tubes, yellow. Pack of 1000

[GY-06245-57](#) Caps for screw-cap tubes, red. Pack of 1000

[GY-06245-58](#) Caps for screw-cap tubes, green. Pack of 1000

[GY-06245-60](#) Caps for screw-cap tubes, violet. Pack of 1000

[GY-06245-53](#) Microcentrifuge adapter for skirted tube 06245-52

### Cole-Parmer® Conical-Bottom Disposable Centrifuge Tubes

Fit your needs—choose sterile or nonsterile

- Blue graduations and writing patch won't rub off
- Self-standing versions (06344-25 and -26) have permanent, molded-in graduations



Capacity (mL)	Tube material	Max RCF (x g)	Cap style (cap included)	Sterile	Catalog number	Packaging	Qty/pk	Price/pk
15	PS	1800	Screw cap/rim seal	No	<a href="#">GY-06344-20</a>	1 bag of 500	500	
15	PS	1800	Screw cap/rim seal	Yes	<a href="#">GY-06344-21</a>	10 trays of 50	500	
15	PPCO	6400	Screw cap/rim seal	No	<a href="#">GY-06344-22</a>	1 bag of 500	500	
15	PPCO	6400	Screw cap/rim seal	Yes	<a href="#">GY-06344-24</a>	20 bags of 25	500	
50	PPCO	10,000	Flat top/rim seal	No	<a href="#">GY-06344-27</a>	1 bag of 500	500	
50	PPCO	10,000	Flat top/rim seal	Yes	<a href="#">GY-06344-29</a>	20 bags of 25	500	
50†	PPCO	N/A	Flat top/rim seal	No	<a href="#">GY-06344-26</a>	20 bags of 25	500	
50†	PPCO	N/A	Flat top/rim seal	Yes	<a href="#">GY-06344-25</a>	20 bags of 25	500	

†Self-standing.

