



Certificate of Compliance



Product Description: Serological Pipettes

Cat. No's. : PB01,PB1,PMB11, PB2,PB22, PB5, PB10,PB10W,
PB25,PB50,PS01,PS1, PMS11, PS2, PS22, PS5,
PS10, PS10W, PS25, PS50, PS100, PM5, PM10,
PM25, BPS5, BPS10, BPS25

Hereby it is confirmed that all the above mentioned products have been manufactured on behalf of Argos Technologies and meets the specifications as listed below. The product is developed, produced and distributed according to a Quality Management System that is certified for compliance with ISO 13485:2003& ISO 9001:2008.

This document certifies that the above mentioned products do not contain, make use of, or involve, at any point of their respective manufacturing process, raw materials of animal origin (including animal proteins). The product is not stored with products of animal origin (including animal proteins) in its manufacturing facilities.

The items have been given appropriate inspection and thus meets Argos Technologies quality control standards. The following parameters are checked on a regular basis.

- Material: Polystyrene - Appearance
- Leak Tightness
- Graduation Accuracy +/- 2%
- Pyrogen Free
- Sterilization Certificate
- Cytotoxicity

York Zhang
Quality Manager

Argos Technologies, Inc.
1551 S. Scottsdale Ct. S200
Elgin, IL 60123

Ph: 800-886-8675 / 847-622-0456
Fx: 847-622-0454
www.argos-tech.com