

1. **Part Number:** 06043-11
2. **Date Issued:** Sept 6, 2022



3. Product Information:

- Description:** Viscous Sample Collector, HDPE, 500mm length, 100mL; 20/PK
- Use:** Single use, disposable liquid sampler (for high viscosity liquids)
- Material of Construction:** High Density Polyethylene (HDPE) - virgin
- Conforms to FDA 21 CFR 177.1520
 - Conforms to EU Regulations 10/2011
 - Conforms to EC Regulations 1935/2004
 - Conforms to REACH regulation 1907/2006
 - BSE/TSE free
- O-ring Seal – virgin Silicone Rubber
- Conforms to FDA CFR 177.2600
- Method of Construction:** Ultrasonic welding
- Nominal Length**
- Fully closed: 523mm
 - Fully extended: 881mm
- Diameter of Outer Tube:** 21mm
- Diameter of hole in tip:** 10mm
- Width of Handle:** 103mm
- Nominal Weight of Sampler:** 83g
- Assembly & Primary Packing Environment:** Class 100,000 Cleanroom
- Individually Bagged?** Yes (heat sealed PE/PET bag)
- Method of Sterilisation:** N/A
- Number of Samplers per Box:** 20
- Recommended Storage Conditions:** Dry and ambient temperature
- Shelf Life:** 5 years from the date of manufacture
4. **Country of Origin:** Made and packed in the UK.
5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.