

**1. Part Number:** 06043-12

**2. Date Issued:** Sept 6, 2022



**3. Product Information:**

**Description:** Viscous Sample Collector, Sterile, HDPE, 500mm length, 100mL; 20/PK

**Use:** Single use, disposable liquid sampler (for high viscosity liquids)

**Material of Construction:** High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- Conforms to REACH regulation 1907/2006
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

**Method of Construction:** Ultrasonic welding

**Nominal Length**

Fully closed:	523mm
Fully extended:	881mm

**Diameter of Outer Tube:** 21mm

**Diameter of hole in tip:** 10mm

**Width of Handle:** 103mm

**Nominal Weight of Sampler:** 83g

**Assembly & Primary Packing Environment:** Class 100,000 Cleanroom

**Individually Bagged?** Yes (heat sealed PE/PET bag)

**Method of Sterilisation:** Gamma Irradiation (25 to 45 kGy)

**Number of Samplers per Box:** 20

**Recommended Storage Conditions:** Dry and ambient temperature

**Shelf Life:** 5 years from the date of manufacture

**4. Country of Origin:** Made and packed in the UK.

**5. Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.