


1. Part Number: 06040-14

2. Date Issued: 22-Sep-22



3. Product Information:

Description: Viscous Sample Collector, Sterile, HDPE, 1000mm length, 190mL;
20/PK

Use: Single use, disposable liquid sampler (for high viscosity liquids)

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- Conforms to REACH regulation 1907/2006
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

Method of Construction: Ultrasonic welding

Nominal Length

Fully closed: 991mm

Fully extended: 1818mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: 10mm

Width of Handle: 103mm

Nominal Weight of Sampler: 158g

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25kGy minimum)

Number of Samplers per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. Country of Origin: Made and packed in the UK.

5. Repeated Use: This is a single-use, disposable sampler. Repeated use will lead to failure.