


1. Part Number: 06240-01

2. Date Issued: 22-Sep-22



3. Product Information:

Description: Volumetric Spoon, Sterile, HDPE, 2.5mL; 100/PK

Use: Single use, disposable measuring spoon

Material of Construction: High Density Polyethylene (HDPE) – rigid, natural

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant
- BSE/TSE free

Method of Construction: Single piece, injection moulded

Nominal Bowl Volume: 2.46ml

Length: 158mm

Bowl Diameter: 23mm

Nominal Weight: 4g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Quantity per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested

4. Country of Origin: Made and packed in the UK.

5. Repeated Use: This is a single-use, disposable device. Repeated use will lead to failure.