


**1. Part Number:** 06264-02

**2. Date Issued:** 22-Sep-22



**3. Product Information:**

Description:	Volumetric Spoon, Sterile, HDPE, 5mL; 100/PK
Use:	Single use, disposable measuring spoon
Material of Construction:	High Density Polyethylene (HDPE) – rigid, natural <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Bowl Volume:	5ml
Length:	165mm
Bowl Diameter:	29mm
Nominal Weight:	4g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

**4. Country of Origin:** Made and packed in the UK.

**5. Repeated Use:** This is a single-use, disposable device. Repeated use will lead to failure.