


1. Part Number: 06264-04

2. Date Issued: 22-Sep-22



3. Product Information:

Description:	Volumetric Spoon, Sterile, HDPE, 20mL; 100/PK
Use:	Single use, disposable measuring spoon
Material of Construction:	High Density Polyethylene (HDPE) – rigid, natural <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Bowl Volume:	20ml
Length:	180mm
Bowl Diameter:	43mm
Nominal Weight:	6g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Country of Origin: Made and packed in the UK.

5. Repeated Use: This is a single-use, disposable device. Repeated use will lead to failure.