

1. **Part Number:** 06264-05  
2. **Date Issued:** 22-Sep-22



**3. Product Information:**

- Description:** Volumetric Spoon, Sterile, HDPE, 30mL; 100/PK
- Use:** Single use, disposable measuring spoon
- Material of Construction:** High Density Polyethylene (HDPE) – rigid, natural
  - FDA acceptable CFR 177.1520
  - EU Regulations 10/2011 compliant
  - EC regulation 1935/2004 compliant
  - BSE/TSE free
- Method of Construction:** Single piece, injection moulded
- Nominal Bowl Volume:** 30ml
- Length:** 187mm
- Bowl Diameter:** 50mm
- Nominal Weight:** 7g
- Moulding & Primary Packing Environment:** Class 100,000 Cleanroom
- Individually Bagged?** Yes (heat sealed PE/PET bag)
- Method of Sterilisation:** Gamma Irradiation (25-45kGy)
- Quantity per Box:** 100
- Recommended Storage Conditions:** Dry and ambient temperature
- Shelf Life:** 5 years from the date of manufacture
- ATEX:** This product has not been ATEX tested
- 4. Country of Origin:** Made and packed in the UK.
- 5. Repeated Use:** This is a single-use, disposable device. Repeated use will lead to failure.