

1. **Part Number:** 06264-10
2. **Date Issued:** 22-Sep-22



3. Product Information:

Description:	Long Handled Scoop, HDPE, 50mL; 50/PK
Use:	Single use, powder scoop
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Width:	40mm
Overall length:	371mm
Handle length:	300mm
Nominal Scoop Volume:	50ml
Nominal Weight:	30g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested
4. Country of Origin:	Made and packed in the UK.
5. Repeated Use:	This is a single-use, disposable sampler. Repeated use will lead to failure.