

**1. Part Number:** 06264-19

**2. Date Issued:** 06/09/22



**3. Product Information:**

Description: Powder Spatula, Sterile, PS, 230mm length; 100/PK

Use: Single use, disposable powder spatula

Material of Construction: High Impact Polystyrene (HIPS) – rigid, white

- FDA acceptable CFR 177.1640
- EU Regulations 10/2011 compliant
- EC Regulations 1935/2004 compliant
- BSE/TSE free

Method of Construction: Single piece, injection moulded

Blade Shape: U

Blade Style: Open Tip

Nominal Length: 230mm

Blade Length: 145mm

Blade Width: 18mm

Nominal Weight: 20g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

**4. Country of Origin:** Made and packed in the UK.

**5. Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.