


**1. Part Number:** 06264-22

**2. Date Issued:** 06/09/22



**3. Product Information:**

Description:	Micro Spatula VN, Sterile, PS, .3mL; 100/PK
Use:	Single use, powder spatula
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	0.3ml
Overall Length	170mm
Blade Length:	40mm
Blade Width:	4mm
Nominal Weight:	3g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (Minimum 25 kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested
<b>4. Country of Origin:</b>	Made and packed in the UK.
<b>5. Repeated Use:</b>	This is a single-use, disposable sampler. Repeated use will lead to failure.