


1. Part Number: 06264-23

2. Date Issued: 06/09/22



3. Product Information:

Description: Micro Spatula N, Sterile, PS, .5mL; 100/PK

Use: Single use, disposable spatula

Material of Construction: High Impact Polystyrene (HIPS)
 - FDA acceptable CFR 177.1640
 - EU Regulations 10/2011 compliant
 - EC Regulations 1935/2004 compliant
 - BSE/TSE free

Method of Construction: Single piece, injection moulded

Nominal Volume: 0.5ml

Overall Length: 170mm

Spatula Blade Length: 40mm

Blade Width: 6mm

Nominal Weight: 4g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Scoops per Box: 100

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. Country of Origin: Made and packed in the UK.

5. Repeated Use: This is a single-use, disposable sampler. Repeated use will lead to failure.