

1. Part Number: 06286-19

2. Date Issued: 06/09/22



3. Product Information:

Description: Slot Sampler, 3 openings, Sterile, HDPE, 620mm length, 100mL;
20/PK

Use: Single use, disposable powder sampler

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding


Overall Length: 620mm

Diameter of Outer Tube: 26mm

Diameter of Inner Blade: 22mm

Nominal Weight of Sampler: 148g

Slot Positions:



Number of Sampling Slots: 3

Sampling Slot Width: 8mm

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of Samplers per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. Country of Origin: Made and packed in the UK.

5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.