

1. **Part Number:** 06286-20
2. **Date Issued:** 09/07/22



3. Product Information:

Description: Slot Sampler, 3 openings, Sterile, HDPE, 970mm length, 170mL; 20/PK

Use: Single use, disposable powder sampler

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding


Overall Length: 968mm

Diameter of Outer Tube: 26mm

Diameter of Inner Blade: 22mm

Nominal Weight of Sampler: 234g

Slot Positions:



- Number of Sampling Slots:** 3
- Sampling Slot Width:** 8mm
- Assembly & Primary Packing Environment:** Class 100,000 Cleanroom
- Individually Bagged?** Yes (heat sealed PE/PET bag)
- Method of Sterilisation:** Gamma Irradiation (25-50 kGy)
- Number of Samplers per Box:** 20
- Recommended Storage Conditions:** Dry and ambient temperature
- Shelf Life:** 5 years from the date of manufacture
4. **Country of Origin:** Made and packed in the UK.
5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.