

1. **Part Number:** 06286-21

2. **Date Issued:** 09/07/22



3. Product Information:

Description: Slot Sampler, 1 opening, Sterile, HDPE, 620mm length, 50mL;
20/PK

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 620mm

Diameter of Blade: 22mm

Diameter of Outer Tube: 26mm

Number of Sampling Slots: 1

Slot Positions:



Sampling Slot Width: 8mm

Sampling Slot Length: 60mm

Nominal Weight of Sampler: 167g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25 – 50 kGy)

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested

- 4. **Country of Origin:** Made and packed in the UK.
- 5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.