

1. **Part Number:** 06286-30

2. **Date Issued:** 09/07/22



3. Product Information:

Description: Powder Sample Collector, Sterile, HDPE, 500mm length, 10mL;
20/PK

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 511mm

Diameter of Outer Tube: 21mm

Nominal Tip Size: 10ml

Width of Handle: 103mm

Length 'L': 63mm

Nominal Weight of Sampler: 93g

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (minimum dose 25 kGy)

Number of Samplers per Box: 20

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. **Country of Origin:** Made and packed in the UK.

5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.