

1. **Part Number:** 06286-44

2. **Date Issued:** Sept 6, 2022



3. Product Information:

Description: Powder Lance, Sterile, HDPE, 400mm length, 75mL; 20/PK

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 401mm

Blade Length: 300mm

Diameter of Blade: 22mm

Nominal Weight of Sampler: 63g

Moulding & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. **Country of Origin:** Made and packed in the UK.

5. **Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.